

HEARTSTART MRX

HEARTSTART MRX ALS MONITOR



PRODUCT DESCRIPTION

The first thing you'll notice about the HeartStart MRx is its large, color display. Look further and you'll see that it has much more. This combination multi-parameter monitor, defibrillator, AED, and 12-Lead unites Philips' industry-leading monitoring and display technologies with superior diagnostic measurements and our patented resuscitation therapies.

Monitoring starts once a patient cable is connected to the device. Equipped for 3- and 5-Lead ECG monitoring with arrhythmia detection, and optional 12-Lead ECG, pulse oximetry, noninvasive blood pressure, and end-tidal CO2, HeartStart MRx is prepared for today's needs and upgradeable to meet tomorrow's.

Its therapies - manual and semi-automatic defibrillation and synchronized cardioversion - feature Philips' patented low-energy SMART Biphasic waveform, which is proven effective in emergency resuscitation and for minimizing post-resuscitation heart dysfunction. No other external defibrillation waveform is supported by more peer-reviewed clinical data. Transcutaneous pacing can be added and the MRx will pace in either demand or fixed mode.

HeartStart MRx logically displays measurements and patient care data on its easy-to-read, backlit, 8.4-inch screen --- the largest in its class!

Numerics and waveforms can be reconfigured, and the screen reorganized, enabling you to quickly locate the information you need most. With wide viewing angles, it displays an event timer, event markers, numeric vital signs, and up to four waves, as well as text prompts, alarms, and battery status indicators. On-screen menus simplify navigation for configuring data, setting and responding to alarms, and accessing additional functionality. And, automated self-tests, straightforward operational checks, data collection, and two long-life batteries make the device easy to operate.

All of these features, measurements, and therapies, plus its compact size, low weight (13.2 pounds), and balanced shape mean that HeartStart MRx has the capabilities you need and the performance you demand for rapid intervention, thorough care, and positive patient outcomes --- that's the big picture.



Features

Standard Features

- ▶ ST/AR Basic algorithm for arrhythmia detection
- ▶ ECG monitoring through monitoring electrodes and defibrillation pads
- Synchronized cardioversion
- ▶ Adjustable ECG size and autogain
- Manual and AED operation
- ▶ SMART Biphasic waveform for defibrillation therapy
- ▶ Large 4-wave color display
- Strip chart printer
- Individual, adjustable volume of QRS beeper, voice prompts, and alerts
- Event summary
- ▶ Configuration mode
- Service mode
- Operational checks
- ▶ Automated self-tests with "ready-for-use" indicator
- ▶ Lithium ion battery with capacity gauge
- ▶ ECG out capability

Optional Features

- ▶ SpO₂ with Fourier Artifact Suppression Technology (FAST)
- ▶ Noninvasive Blood Pressure
- Microstream™ EtCO₂
- Noninvasive Pacing
- ▶ 12-Lead ECG with Philips interpretation algorithm

Standard Accessories

- ▶ Lithium ion battery with capacity gauge
- ▶ Hands-free multifunction electrode cable
- 3-Lead ECG cable
- Disposable monitoring electrodes
- Strip chart recorder paper
- ▶ Carrying case

- Defibrillator test load
- ▶ Instructions for use (User's guide)
- Quick reference cards
- ▶ User training workbook
- Application notes

Package

M3536A HeartStart MRx ALS Monitor

Ordering Option Information

A01	SpO ₂

A02 SpO₂ and NBP

A03 SpO₂, NBP, and EtCO₂

A04 EtCO₂

B01 External Pacing

B02 12-Lead ECG acquisition

C02 Water Resistant External Paddles

C03 Data Card

C05 Additional Battery

C06 AC Power Module

C07 Barrel style Pad Cable (Replacement for Standard Pad

Cable)

C15 5-Lead ECG Cable

LP1 Instructions for Use

LP2 User Training Video (English only)

LPK Alternate Label Set, AED/Energy Select

SM1 Service Manual (English only)

SM3 Service Training Video (English only)



Upgrades/Supplies/Accessories

Upgrades

 M3530A
 SpO2

 M3531A
 NBP

 M3532A
 EtCO2

 M3533A
 Pacing

M3534A 12-Lead ECG

Option B02 - Acquisition
Option B03 - Transmission
Option B04 - Wider Printer

M5527A External Paddles with Paddle Tray

Option C02 - Water Resistant Paddles

External Paddles

M3543A Water Resistant External Paddles

M4745A Sterilizable External Paddles

Internal Paddles

7.5 cm Switchless M1741A M1742A 6.0 cm Switchless M1743A 4.5 cm Switchless M1744A 2.8 cm Switchless M4741A 7.5 cm Switched M4742A 6.0 cm Switched 4.5 cm Switched M4743A M4744A 2.8 cm Switched 989803127121 Large Disposable Switched 989803127131 Medium Disposable Switched 989803127141 Small Disposable Switched 989803127151 Large Disposable Switchless 989803127161 Medium Disposable Switchless 989803127171 Small Disposable Switchless

Internal Paddles Adapter Cable

Multifunction Electrode Pads

M4740A

M3713AAdult PlusM3716AAdult RadiolucentM3717APediatric Plus

M3718A Adult Radiotransparent/Reduced Skin

Multifunction Electrode Pads

M3719A Pediatric Radiotransparent/Reduced Skin

Pads Cable

M3508A Defib Hands-free, plug style

ECG Monitoring Electrodes

M2202A High-Tack Foam, 5 electrodes/pack (60 packs/case)

ECG Cables

12-Lead Cable Set

M3525A 10-lead ECG Trunk Cable, 12-pin Connector (for 3-Lead

and 12-Lead use)

M3526A 3-lead ECG Set and Plug with Snap (AAMI)
M3527A Add 7-lead ECG Set for 12-Lead use (AAMI)
M3528A 3-lead ECG Set and Plug with Snap (IEC)
M3529A Add 7-lead ECG Set for 12-Lead use (IEC)

M1949A 10-lead ECG Patient Trunk Cable, 12-pin ECG Input

Connector (for 5-Lead and 12-Lead use)

M1968A 10-electrode Cable Set, Extremities, Grabber (use with

M1976A) (AAMI)

M1976A 10-electrode Cable Set, Chest, Grabber (use with

M1968A) (AAMI)

M1971A 10-electrode Cable Set, Extremities, Grabber (use with

M1978A) (IEC)

M1978A 10-electrode Cable Set, Chest, Grabber (use with

M1971A) (IEC)

3-Lead Cable Set

M1500A 3-lead ECG Trunk Cable (AAMI)

M1605A 3-lead ECG Lead Set with Snaps (AAMI)

M1510A 3-lead ECG Trunk Cable (IEC)

M1615A 3-lead ECG Lead Set with Snaps (IEC)

5-Lead Cable Set

M1520A 5-lead ECG Trunk Cable (AAMI)

M1625A 5-lead ECG Lead Set with Snaps (AAMI)

M1530A 5-lead ECG Trunk Cable (IEC)

M1635A 5-lead ECG Lead Set with Snaps (IEC)

Sync Cables

M1783A Sync Cable (8 ft.)

M5526A Sync Cable (25 ft.)



SpO ₂ Sens	sors/Cables	NBP				
M1191A	Reusable SpO ₂ Sensor - Adult Finger		Disposable Blood Pressure Cuffs			
M1192A	Reusable SpO ₂ Sensor - Pediatric/Small Adult Finger	M4572A	Soft Single-Patient Disposable Cuff - Infant			
M1193A	Reusable SpO ₂ Sensor - Adult Thumb	M4573A	Soft Single-Patient Disposable Cuff - Pediatric			
M1194A	Reusable SpO ₂ Sensor - Adult Ear Clip	M4574A	Soft Single-Patient Disposable Cuff - Small Adult			
M1195A	Reusable SpO ₂ Sensor - Infant	M4575A	Soft Single-Patient Disposable Cuff - Adult			
	, 2	M4576A	Soft Single-Patient Disposable Cuff - Adult X-Long			
M1941A	SpO ₂ Extension Cable, 2 m	M4577A	Soft Single-Patient Disposable Cuff - Large Adult			
M1943A	Nellcor SpO ₂ Sensor Adapter Cable, 1 m (use with	M4578A	Soft Single-Patient Disposable Cuff - Large Adult X-			
WIISTSA	M1902/3/4B)		Long			
		M4579A	Soft Single-Patient Disposable Cuff - Thigh			
NBP		EtCO ₂				
	Interconnect Cable	_	Intubated Circuits			
M1598B	Adult Pressure (5 ft.)	M1920A	Filter Line Set - Adult/Pediatric (25 sets/case)			
M1599B	Adult Pressure (10 ft.)	M1921A	Filter H Set - Adult/Pediatric (25 sets/case)			
		M1923A	Filter H Set - Infant/Neonatal (yellow, 25 sets/case)			
	Reusable Blood Pressure Cuffs					
40400 A	Reusable NBP Cuff Kit, 3 sizes (pediatric, adult, large adult)		Non-Intubated Dual Purpose Circuits (CO ₂ + O ₂)			
40400B	Reusable NBP Cuff Kit, 5 sizes	M2520A	Smart CapnoLine - Pediatric			
	(infant, pediatric, adult, large adult, thigh)	M2521A	Smart CapnoLine - Intermediate			
40401A	Traditional Reusable NBP Cuff - Infant	M2522A	Smart CapnoLine - Adult			
40401B	Traditional Reusable NBP Cuff - Pediatric					
40401C	Traditional Reusable NBP Cuff - Adult		Non-Intubated Single Purpose Circuits (CO ₂)			
40401D	Traditional Reusable NBP Cuff - Large Adult	M2524A	Smart CapnoLine - Pediatric			
40401E	Traditional Reusable NBP Cuff - Thigh	M2525A	Smart CapnoLine - Intermediate			
		M2526A	Smart CapnoLine - Adult			
M4552A	Antimicrobial Reusable NBP Cuff - Infant					
M4553A	Antimicrobial Reusable NBP Cuff - Pediatric	Power				
M4554A	Antimicrobial Reusable NBP Cuff - Small Adult	M3538A	Lithium Ion Battery with Capacity Gauge			
M4555A	Antimicrobial Reusable NBP Cuff - Adult	M3539A	AC Power Module			
M4557A	Antimicrobial Reusable NBP Cuff - Large Adult	M5529A	DC Power Module			
M4559A	Antimicrobial Reusable NBP Cuff - Thigh	Paper				
		40457C	50 mm Chemical Thermal, Gray Grid (10 rolls)			
M1572A	Multi-Patient Comfort Cuffs - Pediatric	40457D	50 mm Chemical Thermal, Gray Grid (80 rolls)			
M1573A	Multi-Patient Comfort Cuffs - Small Adult					
M1574A	Multi-Patient Comfort Cuffs - Adult					
B84575A	M III D II 10 (10 (1 A I II					

M1575A

Multi-Patient Comfort Cuffs - Large Adult



Color Handle

M5521AGreen and color label setM5522ABlue and color label setM5523AYellow and color label setM5524ARose and color label setM5525AGray and color label set

Miscellaneous

M3725A Test Load for use with M3508A Pad Cable

M3541A Carrying Case (includes 3 accessory pouches and

shoulder strap)

M3545A Data Card and Tray

M5528A Vehicle Wall Mount



HEARTSTART MRX

Specifications

General

Dimensions with pads: 12.4 in. (W) x 8.3 in. (D) x 11.7 in. (H)

(31.5 cm x 21.0 cm x 29.5 cm)

Dimensions with paddles: $13.4 \text{ in. (W)} \times 8.3 \text{ in. (D)} \times 13.6 \text{ in. (H)}$

(34.0 cm x 21.0 cm x 34.5 cm)

Weight: 13.2 lbs. including pads, pads cable, full roll of paper, and battery. Incremental weight of external standard paddles and paddle tray is 2.5 lbs. Additional battery weighs 1.6 lbs.

Defibrillator

Waveform: Biphasic Truncated Exponential. Waveform parameters adjusted as a function of patient impedance

Shock Delivery: Via multifunction defib electrode pads, or paddles

Delivered Energy Accuracy:

Selected Energy	Nomi	Accuracy						
	25	50	75	100	125	150	175	
1J	1.2	1.3	1.2	1.1	1.0	0.9	0.8	±2J
2J	1.8	2.0	2.0	1.9	1.7	1.6	1.5	±2J
3J	2.8	3.0	3.0	3.1	3.0	2.9	2.7	±2J
4J	3.7	4.0	4.0	4.1	4.2	4.2	4.0	±2J
5J	4.6	5.0	5.1	5.1	5.2	5.2	5.0	±2J
6J	5.5	6.0	6.1	6.2	6.3	6.3	6.1	±2J
7J	6.4	7.0	7.1	7.2	7.3	7.3	7.1	±2J
8J	7.4	8.0	8.1	8.2	8.4	8.3	8.1	±2J
9J	8.3	9.0	9.1	9.3	9.4	9.4	9.1	±2J
10J	9.2	10	10	10	10	10	10	±2J
15J	14	15	15	15	16	16	15	±15%
20J	18	20	20	21	21	21	20	±15%
30J	28	30	30	31	31	31	30	±15%
50J	46	50	51	51	52	52	50	±15%
70J	64	70	71	72	73	73	71	±15%
100J	92	100	101	103	104	104	101	±15%
120J	110	120	121	123	125	125	121	±15%
150J	138	150	152	154	157	156	151	±15%
170J	156	170	172	175	177	177	172	±15%
200J	184	200	202	206	209	209	202	±15%

Charge Time: Less than 5 seconds to 200 Joules with a new, fully charged Lithium Ion battery pack at 25°C. Less than 15 seconds when operating without a battery, using the M3539A AC power module alone at 90% rated mains voltage

Patient Impedance Range:

Minimum: 15 Ohm (internal defibrillation); 25 Ohm (external defibrillation)

Maximum: 180 Ohm

Note: Actual functional range may exceed the above values.

Manual Mode

Manual Output Energy (Selected): 1-10, 15, 20, 30, 50, 70, 100, 120, 150, 170, 200 Joules

Controls: On/Off Therapy Knob, Charge, Shock, Sync, Print, Mark Event, ECG Lead Select, Alarm Pause, Event Review

Energy Selection Control: Front panel therapy knob

Charge Control: Front panel button; button on external paddles

Shock Control: Front panel button; buttons on external or switched internal paddles

Synchronized Control: Front panel SYNC button

Indicators: Text prompts, audio alerts, QRS beeper, battery status, Ready For Use, external power/battery charging, Sync mode

Armed Indicators: Charging tone, charged tone, flashing shock button, and energy level indicated on display

AED Mode

AED Energy Profile: 150 Joules nominal into a 50 Ohm test load

Text and Voice Prompts: Extensive text/audible messages guide user through configured protocol

AED Controls: On/Off, Shock

Indicators: Monitor display messages and prompts, voice prompts, battery status, Ready For Use, external power/battery charging

Armed Indicators: Flashing shock button, charged tone, and voice prompts

ECG Analysis: Evaluates patient ECG and signal quality to determine if a shock is appropriate and evaluates connection impedance for proper defibrillation pad contact

Shockable Rhythms: Ventricular fibrillation and certain ventricular tachycardias, including ventricular flutter and polymorphic ventricular tachycardia

Shock Advisory Algorithm Sensitivity and Specificity: Meets AAMI DF-



ECG and Arrhythmia Monitoring

Inputs: Up to four (4) ECG waves may be viewed on display and up to 2 waves printed simultaneously. Lead I, II, or III is obtained through the 3-lead ECG cable. With a 5-lead cable, leads aVR, aVL, aVF, and V can also be obtained. Pads ECG is obtained through 2 multifunction defibrillation electrode pads.

Lead Fault: Leads Off message and dashed line appear on the display if an electrode or lead becomes disconnected. Lead Off indicator in wave sector

Pad Fault: Dashed line appears on the display if a pad becomes disconnected.

Heart Rate Display: Digital readout on display from 15 to 300 bpm, with an accuracy of $\pm 10\%$

Heart Rate/Arrhythmia Alarms: HR, Asystole, VFIB/VTACH, VTACH, Extreme Tachy, Extreme Brady, PVC rate, Pacer not capture, Pacer not pacing

Hands Free Defibrillation Cable Length:

M3508A = 5 ft. (1.6 m)M3507A = 7 ft. (2.2 m)

ECG Cable Length: 9 ft. (2.7 m)

Common Mode Rejection: Greater than 90 dB measured per AAMI

standard for cardiac monitors (EC 13)

ECG Size: 2.5, 5, 10, 20, 40 mm/mV, autogain

Frequency Response:

AC Line Filter: 60 Hz or 50 Hz

Pads ECG for Display: Monitor (0.15-40 Hz) or EMS (1-30 Hz) Pads ECG for Printer: Monitor (0.15-40 Hz) or EMS (1-30 Hz) Leads ECG for Display: Monitor (0.15-40 Hz) or EMS (1-30 Hz)

Leads ECG for Printer: Diagnostic (0.05-150 Hz) or Monitor (0.15-40 Hz) or

EMS (1-30 Hz)

Patient Isolation (defibrillation proof):

ECG: Type CF SpO₂: Type CF EtCO₂: Type CF NBP: Type CF

External Defib: Type BF Internal Defib: Type CF

Display

Size: 8.4 in. diagonal (128 mm x 171 mm)

Type: TFT Color LCD

Resolution: 480 x 640 pixels (VGA)

Sweep Speed: 25 mm/s nominal (stationary trace; sweeping

erase bar) for ECG and SpO₂; 6.25 mm/s for CO₂

Wave Viewing Time: 5 seconds (ECG)

Battery

Type: 6.3 Ah, 14.8 V, rechargeable, Lithium Ion

Dimensions: 6.5 in. (H) x 3.8 in. (W) x 1.6 in. (D) (165 mm x 95 mm x 42

mm)

Weight: 1.6 lb. (0.73 kg)

Charge Time: Approximately 3 hours to 100%. Approximately 2 hours to 80%, indicated by battery gauge indicator. Charging the battery at temperatures above 45°C may degrade battery life.

Capacity: At least 5 hours of monitoring with ECG, SpO_2 , and CO_2 monitored continuously and NBP measured every 15 minutes, or at least 50 200-joule shocks (on a new, fully charged battery at room temperature, 25° C). At least 3.5 hours while pacing at 180 ppm at 160 mA and monitoring as described above

Battery Indicators: Battery gauge on battery, capacity indicator on display; flashing RFU indicator, chirp, and LOW BATTERY message appears on display for low battery condition*

*Low battery condition triggered when 10 minutes of monitoring time and 6 maximum energy discharges remain (with a new battery at room temperature, 25° C)

Storage: Storing the battery for extended periods at temperatures above 40°C will reduce battery capacity and degrade battery life.



HEARTSTART MRX

Thermal Array Printer

Continuous ECG Strip: The Print key starts and stops the strip. The printer can be configured to run real time or with a 10-second delay. The strip prints the primary ECG lead with event annotations and measurements.

Auto Printing: The printer can be configured to automatically print on Marked Events, Charge, Shock, and Alarm. When an alarm condition occurs, the unit prints the primary ECG wave and the alarming wave, if configured.

Reports: The following can be printed: Event Summary, 12-Lead, Operational Check, Configuration, Status Log, and Device Information

Speed: 25 or 50 mm/s with an accuracy of \pm 5%

Amplitude Accuracy: ± 5% or ± 40 uV, whichever is greater

Paper Size: 50 mm (W) by 30 m (100 ft.) (L)

Noninvasive Pacing

Waveform: Monophasic Truncated Exponential

Current Pulse Amplitude: 10 mA to 160 mA (5 mA resolution);

accuracy 10% or 5 mA, whichever is greater

Pulse Width: 40 ms with \pm 10% accuracy

Rate: 30 ppm to 180 ppm (10 ppm increments); accuracy \pm 1.5%

Modes: Demand or Fixed Rate

Refractory Period: 340 msec (30 to 80 ppm); 240 msec (90 to 180 ppm)

SpO₂ Pulse Oximetry

Range: SpO2: 0-100%; pulse rate: 30 to 300 bpm

Resolution: 1%

Display Update Period: 1 sec. typical numeric update rate

Accuracy with:

M1191A sensor - 1 standard deviation 70% to 100%, \pm 2.5% M1192A sensor - 1 standard deviation 70% to 100%, \pm 2.5% M1193A sensor - 1 standard deviation 70% to 100%, \pm 2.5% M1194A sensor - 1 standard deviation 70% to 100%, \pm 4.0% M1195A sensor - 1 standard deviation 70% to 100%, \pm 4.0% NELLCOR sensors - 1 standard deviation 80% to 100%, \pm 3.0%

Pulse Rate Accuracy: 2% or 1 bpm (whichever is greater)

Alarm Range:

Low Limit: 50 to 99% (Adult/Pediatric) High Limit: 51 to 100% (Adult/Pediatric)

Alarm Delay: 10 seconds

Noninvasive Blood Pressure

Pressure Range:

Systolic: 40-260 mmHg Diastolic: 20-200 mmHg

Initial Pressure: 160 mmHg (Adult); 120 mmHg (Pediatric)

Maximum Pressure: 280 mmHg

Overpressure Safety Limits: Maximum of 300 mmHg
Cuff Inflation Time: 50 second maximum (pediatric or adult)

Accuracy: ±3 mmHg

Alarm Range:

Systolic high limit: 35 - 270 (Adult), 35 - 180 (Pediatric) Systolic low limit: 30 - 265 (Adult), 30 - 175 (Pediatric) Diastolic high limit: 15 - 245 (Adult), 15 - 150 (Pediatric) Diastolic low limit: 10 - 240 (Adult), 10 - 145 (Pediatric) Mean high limit: 25 - 255 (Adult), 25 - 160 (Pediatric) Mean low limit: 20 - 250 (Adult), 20 - 155 (Pediatric)

Rated Life: 50,000 measurement cycles (36/day for 2.3 years)

Auto Mode Repetition Time: 1, 2.5, 5, 10, 15, 30, 60, or 120 minutes

Measurement Time: Auto/manual mode: 30 seconds (average) @ HR > 60

bpm, 170 seconds (maximum)
Interconnect Tube Length:
M1598B Connect tubing (1.5 m)
M1599B Connect tubing (3 m)

End-Tidal Carbon Dioxide

Range: 0 to 99 mmHg

Resolution: 1mmHg (0.1 kPa)

Accuracy: For values between 0 and 38 mmHg: ± 2 mmHg. For values between 39 and 99 mmHg: $\pm 5\%$ of reading + 0.08% for every 1 mmHg (above 40 mmHg). Values read at sea level after ≥ 20 min. warm up

Alarm Range:

Low Limit: 10 to 94 mmHg (Adult/Pediatric) High Limit: 20 to 95 mmHg (Adult/Pediatric)

Sample Size: 50 ml per min.



Airway Respiration Rate

Range: 0 to 150 rpm Resolution: 1 rpm

Accuracy:

0 to 40 rpm ±1 rpm 41 to 70 rpm ±2 rpm 71 to 100 rpm ± 3% 101 to 150 rpm ±5%

Alarm Range:

Low Limit: 0 to 99 rpm (Adult/Pediatric)
High Limit: 10 to 100 rpm (Adult/Pediatric)

Apnea Alarm Time: 10-40 seconds, in increments of 5

Calibration Gas for CO₂ Measurement

Ingredients: 5% Carbon Dioxide, 21% Oxygen, 74% Nitrogen

Cylinder Size: BD

Method of Preparation: Gravimetric

Blend Tolerance: 0.03%
Accuracy: 0.03% absolute
Moisture: 10 PPM Maximum
Expiration Period: 2 years
Pressure: 144 PSIG, Volume: 1L

12-Lead ECG

Inputs: With a 10-Lead cable, leads I, II, III, aVR, aVL, aVF, V/C1-V/C6 can be obtained. All 12-Lead ECG waves can be viewed on the display simultaneously. All 12 leads can be printed on the strip chart printer in 3x4 format.

12-Lead ECG Transmission: Data card to a PC

Receiver Type: HeartStart Event Review Pro

Patient Data Storage

events AND 50 12-Lead ECG reports

Internal Event Summary: 8 hours of continuous ECG waveforms and events AND 50 12-Lead ECG reports

Data Card Event Summary: 8 hours of continuous ECG waveforms and

Environmental

Temperature: 0° C to 45°C operating, -20° to 70°C storage

Humidity: Up to 95% relative humidity

Altitude:

Operating: 0 to 15,000 ft. (0 to 4,500 m) Storage: 0 to 15,000 ft. (0 to 4,500 m)

Mechanical Shock:

Operating Impact: Half-sine waveform, duration < 3 ms, velocity change > 111 in/s \pm 5% 3 times on all six faces (i.e., axes)

Non-operating Shipping: Trapezoidal waveform, acceleration \geq 30 g, velocity change=742 cm/s \pm 10% on all six faces (i.e., axes)

Bump: EN60068-2-29 Bump (Half-sine, 40 g peak, 6 msec duration, 1,000 bumps \times 3 axes)

Free fall: IEC 68-2-32 Free fall. Drops on all faces onto a steel surface (excluding bed rail hook)

- ? in. (? cm) with carrying case
- ? in. (? cm) without carrying case
- ? in. (? cm) onto handle

Vibration:

Operating: MIL STD 810E 514.4 Category 6 Helicopter, General Storage, IJH60

Non-Operating:

- IEC 68-2-6 Vibration (sinusoidal) (10-57 Hz, +/- 0.15mm; 58-150 Hz, 2g; 20 sweeps x 3 axes)
- IEC 68-2-64 Vibration, broad-band random (10-20 Hz, 0.05 g^2 /Hz; 20-150 Hz, -3 dB/octave; 150Hz, 0.0065 g^2 /Hz; 1.5 hours x 3 axes)

Solids/Water Resistance: IP24. Water testing performed with patient cables connected to the device

EMC: Meets EN 60601-1-2:2001

Safety: Meets EN 60601-1, UL 2601-1, CAN/CSA C22.2 No. 601-1

Other Considerations: Device not suitable for use in the presence of concentrated oxygen or a flammable anesthetic mixture with air, oxygen, or nitrous oxide

Mode of Operation: Continuous

AC Line Powered:

Input: 100-240 VAC, 50-60 Hz, 1-0.46 A

Output: 18 V, 5 A, 90 W

Battery Powered: 14.8 V Rechargeable, Lithium Ion



Philips Medical Systems is part of Royal Philips Electronics

UNITED STATES:

Philips Medical Systems Cardiac and Monitoring Systems 3000 Minuteman Road Andover, MA 01810 Tel: (800) 934-7372

CANADA:

Philips Medical Systems 2660 Matheson Blvd. E. Mississauga, Ontario L4W 5M2 Tel: (800) 291-6743

. . . (000, 20. 07. .

Philips Medizinsysteme Böblingen GmbH Cardiac and Monitoring Systems Hewlett-Packard Str. 2 71034 Böblingen Germany

EUROPE, MIDDLE EAST AND AFRICA:

Fax: (+49) 7031 463 1552

LATIN AMERICA HEADQUARTERS:

Philips Medical Systems 5200 Blue Lagoon Drive 9th Floor Miami, FL 33126

Tel: (305) 267-4220

ASIA PACIFIC HEADQUARTERS:

Philips Medical Systems 26/F Hopewell Centre 17 Kennedy Road Wanchai Hong Kong Tel: (+852) 2821 5888

Fax: (+852) 2821 5888 Fax: (+852) 2527 7136

www.medical.philips.com

© 2003 Koninklijke Philips Electronics North America. All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder. Koninklijke Philips Electronics N.V. reserves the right to make changes in specifications or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Published December 2003 4522 981 85091