Laerdal BaXstrap Spineboard - Directions for Use

REF 982500, 982599, 982600 and 982699

Brand Name: Laerdal BaXstrap Spineboard
Common Name: Backboard/Spineboard/Longboard
Description: Full-body length spinal immobilization device

Intended Use: A base structure to be used with other adjunct cervical spine and head immobilization devices to facilitate in-line, neutral immobilization and transport of adult and pediatric patients.

Indications for Use: For use with patients having known or suspected spine or spinal cord injuries when patient immobilization and transportation is indicated.

Contraindications: None known

Notes

- Always follow local protocol for patient movement and immobilization techniques.
- Not made with natural rubber latex.
- The product is X-ray, CT, and MRI compatible.

Warnings and Cautions: Read and thoroughly understand these directions. These directions are not intended as, and do not in any way constitute medical treatment advice, instruction in or teaching of medical practice, or imply or convey any rights, medical direction or authority. This and other spinal immobilization adjunct devices should only be used by personnel who have received adequate training. Inspect the BaXstrap prior to initial use and before each patient use. Pay particular attention to the underside of the head end. If the product develops punctures, broken pins or cracks then discontinue use immediately.

Made in USA for:
LAERDAL MEDICAL AS
P.O. Box 377
Tanke Svilandsgate 30,
4002 Stavanger,
Norway
Tel. +47 51 51 17 00
Fax. +47 51 52 35 57


Distributed in USA by:
LAERDAL MEDICAL CORPORATION
167 Myers Corners Road,
Wappingers Falls New York 12590
Tel. +1 (800) 431-1055, +1 (845) 297-7770
Fax +1 (800) 227-1143
E-mail: customerservice@laerdal.com

Distributed in Canada by:
LAERDAL MEDICAL CANADA LTD.
151 Nashdene Rd., Unit #45,
Toronto, ON, Canada, M1V 4C3
Tel. +1 (416) 298-9600
Fax +1 (416) 298-8016
E-mail: savelives@laerdal.ca

US Patent 5,950,627

BaXstrap is a trademark of Laerdal Medical AS.
All rights reserved.

EN 1865 compliant

Directions for Use:

- Immobilize the cervical spine using a proper Cervical Immobilization Device
- Position the patient on the BaXstrap using appropriate number of personnel.
- Properly fasten/strap the patient to the BaXstrap to protect against patient movement in any direction (up, down, left or right).
- Secure the head to the BaXstrap using a proper Head Immobilization Device.
- Assure in-line neutral immobilization of the head, neck, torso and pelvis.

Cleaning: Non-sterile, non-critical infection risk. Withstands cleaning in soap and water (1:16) and 1:100 household bleach: water.

Environmental Conditions:
Temperature (Operating and Storage):
-34 °C to +52 °C (-29 °F to +126 °F)
Humidity (Operating and Storage):
5-95% RH non-condensing.

The product has a 5 year limited warranty. See the enclosed "Laerdal Global Warranty" for terms and conditions.
The warranty is also available at www.laerdal.com.

This product is in compliance with the essential requirements of EU Council directive 93/42/EEC as amended by EU Council directive 2007/47/EC.

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© LAERDAL 2019. All rights reserved. 5885 Rev E Printed in USA