# SimMom

Important Product Information

**Limbs & Things**  
BRINGING SKILLS TRAINING TO LIFE

**Laerdal**  
helping save lives

## CAUTIONS AND WARNINGS

### Disclaimer

Use of the SimMom patient simulator system to train personnel should be undertaken under supervision of suitably trained technical or medical personnel with an understanding of educational principles as well as recognized medical protocols.

As with all patient simulators or other such training devices there may be approximations, variations and inaccuracies in anatomical features and the physiological modeling. This being the case, Laerdal Medical does not guarantee that all features are completely accurate.

### General Simulator Handling

Use only SimMom & PROMPT Birthing Lubricant (Laerdal Catalogue No. 377-18850, Limbs & Things No. 10193) to lubricate baby, cervix, birth canal, and modules.

Do not use any other lubricant. Do not use PROMPT Birthing Lubricant (Laerdal Catalogue No. 376-02950, Limbs & Things No. 50181) or *former* SimMom Birthing Lubricant (Laerdal Catalogue No. 377-14450, Limbs & Things No. 10191). Non-approved lubricants can damage the system.

### Warnings

- SimMom & PROMPT Birthing Lubricant is not for personal use.
- Avoid pinch hazards - Do not use the patient simulator without the external skins.
- Do not use automated chest compression machines on the patient simulator.

### Caution Latex

This product contains Natural Rubber latex, which may cause allergic reactions when in contact with humans.

## Defibrillation Hazards

A conventional defibrillator may be used on SimMom. During defibrillation, the defibrillator and patient simulator may present a shock hazard. When using a defibrillator on the simulator, take the following safety precautions:

### Warnings

- Observe all standard safety precautions for the use of defibrillators.
- Follow defibrillation protocol by avoiding contact between the external paddles and any of the electrode sites while defibrillating.
- Failure to follow safety measures could result in injury or death.
- When performing defibrillation, use the defibrillator connectors or the zap plate mounted on the patient simulator's chest. Do not use the ECG connectors, which are designed exclusively for ECG monitoring. Defibrillation on the ECG connectors will damage the internal electronics of the patient simulator and may cause personal injury.
- Do not defibrillate the patient simulator when it is turned OFF or if it is not functioning normally.
- Do not perform defibrillation when the patient simulator is resting on a wet surface.

### Warnings

- The patient simulator torso must always be kept dry. Sudden changes in temperature may result in condensation collecting on electronic components, which could pose a shock hazard. Allow the patient simulator to acclimate before defibrillating.
- The patient simulator must not come into contact with electrically conductive surfaces or objects during defibrillation.
- Avoid use in all flammable environments. For example, high levels of pure oxygen should be avoided during defibrillation. Ensure good ventilation if concentrated oxygen is used near the patient simulator:
- To prevent torso skin electrode pitting, do not apply conductive gel or conductive defibrillation pads intended for patient use.
- Pressing down too hard on the defibrillation connectors during defibrillation may also cause arcing and pitting.
- Do not defibrillate the patient simulator if the torso skin is not in place.

## REGULATORY INFORMATION

### Federal Communications Commission Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1 This device may not cause harmful interference, and
- 2 This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and the receiver

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected
- Consult the dealer or an experienced radio/TV technician for help.

This Class B digital apparatus complies with Canadian ICES-003.

## EU

**CE** The product is in compliance with the essential requirements of Council Directive 2004/108/EC on electromagnetic compatibility (EMC). The product is in compliance with Council Directive 2011/65/EU on restriction of the use of certain hazardous substances (RoHS).

## WASTE HANDLING

Dispose of in accordance with local requirements and regulations.

---

## SPECIFICATIONS

### Size and Weight

Length/width (patient simulator only)	69 in x 19 in (175.3 cm x 48.3 cm)
Width at shoulders (patient simulator only)	19 in (48.3 cm)
Width at hips (patient simulator only)	19 in (48.3 cm)
Weight (patient simulator only)	61.5 lbs (28 kg)

### Patient Simulator Power

Link Box battery	5 V
------------------	-----

### Air & CO<sub>2</sub> Pressure

External air connection	Max 16 psi
-------------------------	------------

### Temperature Limits

Operating temperatures	+4 °C to 40 °C (39 °F to 104 °F)
Storage temperatures	-15 °C to 50 °C (5 °F to 122 °F)

### Environment - Patient Simulator Only

Relative humidity	20% -90% (non-condensing)
DO NOT use outdoors in wet conditions	
Not tested with salt spray	

### Material Chart for Patient Simulator

Clothes	Cotton, Nylon
Skins and airways	Silicone, PVC
External hard plastics	PVC, PC-ABS
Inner plastics	Silicone, PVC, ABS, POM, Nylon, Polyurethane, Latex (pneumothorax bladders)
Metal components	Aluminum, Brass, Steel

### Recommended Hardware Specifications

Intel i-core generation 3 or newer
Score over 3000 PassMark - CPU Mark
4 GB RAM
120GB hard disk space
1366x768 or better
Headset w/microphone
USB ports – 2

### Software Requirements

Windows 7 or Windows 8
100% DPI

### Cleaning Fluids

To clean the patient simulator use only a mild solution of liquid soap and water
--

### Simulated IV Fluids

Only use distilled or deionized water to simulate IV and IO fluids
--

### Reservoir Capacities

Urine reservoir	400 ml
Blood reservoir	800 ml

### Wig size

Use a large wig
-----------------

---

## Airway Adjustment Sizes

7.5 endotracheal tube
#4 LMA
Large Adult or Trainer Combitube
KING LT – #4

## Urine Catheter Sizes

12 French
-----------

## Blood Pressure

Pressure range:	0-300 mmHg
Accuracy	+/-4 mmHg
Calibration	Pressure sensor must be calibrated to sphygmomanometer
Available pulses	Carotid, brachial, and radial. Only two pulses can be palpated at once.
Palpated BP	Radial pulse turns off at systolic BP. Brachial pulse turns off at 20 mmHg, to prevent noise generation in auscultation area.

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© 2014 Laerdal Medical AS. All rights reserved.
SimMom™ is a trademark and SimMom® is a
registered trademark of Laerdal Medical AS.