EN

FR

DE

ES

IT

BR

NL

PL

JA

ZH

KO

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User Guide


Table of Contents

Introduction 1
 Cautions and Warnings (Intended Use) 1
 Items Included 2
 Skills Taught 2
 Simulation Software 2
 Unpack 3
 Attaching the Legs 3
 SimPad and LLEAP Connections 4
 Getting Started 4
 Cricothyroid Membrane/Neck Skin Collar 4
 Pupil Replacement 4
 In Use 4
 Airway Management 4
 Obstructed Airway 4
 Tension Pneumothorax Decompression 4
 IV Cannulation 5
 Intramuscular Injections 5
 Urinary Catheterization 5
 Simulator Functionalities 6
 Cardiac Related Skills 6
 Blood Pressure (BP) Arm 6
 Heart, Breath and Bowel Sounds and Speech 7
 Spontaneous Breathing and Airway Closures 7
 Maintenance 7
 Cricothyroid Membrane/Neck Skin Collar 7
 IV Cannulation 7
 Intramuscular Injections 7
 Urinary Catheterization 7
 Pneumothorax 7
 A) Bilateral Mid Clavicular Sites 7
 B) Mid Axillary Site (RIGHT) 8
 Chest Drain 8
 Attaching/Removing Arms 8
 IV Arm Skin and Vein Replacement 8
 After Use 9
 Cleaning 9
 Extend Cleanliness 9
 Technical Data 9
 General 9
 ALS Simulator 9
 Replacement Parts 9

Introduction

The Advanced Life Support (ALS) Simulator is a realistic interactive training simulator for simulating a wide range of advanced life saving skills in pre-hospital emergencies. The simulator responds to clinical intervention, instructor control, and comprehensive pre-programmed scenarios for effective practice of diagnosis and treatment of a patient. With spontaneous breathing, airway control, voice, sounds, ECG, and many other clinical features, the ALS Simulator is the fully functional pre-hospital simulator.

The manufacturing quality of this simulator should provide many sessions of training when reasonable care and maintenance are practiced.

Laerdal Recommends
 Endotracheal Tube – Size 7.5
 Laryngeal Mask Airway – # 4
 Combitube- Large Adult or Trainer
 KING LT – # 4
 Needle or Catheter – Size 21-gauge or smaller

 **Cautions and Warnings**
 This product contains Natural Rubber latex which may cause allergic reactions when in contact with humans.

Cautions and Warnings (Intended use)

- There are electronic components mounted inside the airway head. The following techniques should not be performed on this simulator due to the inability to properly sanitize the airway if they are performed:
 - Mouth-to-mouth/Mouth-to-mask ventilation
 - Insertion of simulated vomit for suctioning
- Lubricate the oral and nasal airways liberally with the lubricant provided prior to inserting any instrument, tube or airway device into the airway. Additionally, instruments and tubes should also be lubricated prior to use.
- Do not allow the simulator’s skin to come in direct contact with ink or photocopied paper, as this can permanently stain the skin. Avoid using colored plastic gloves when handling the simulator, as they may also cause discoloration.
- Care should be taken when palpating the pulses to not use excessive force as this will result in no pulse being felt. Only two individual pulses can be felt at the same time.
- If a training session involves the administration of fluids and/or drugs into the IV arm, empty the arm immediately following the training session.
- **Defibrillation**
 - Only apply the defibrillator to a defibrillation post or zap plate which is properly mounted on the simulator’s chest.
 - Do not provide more than 2 x 360J defibrillator discharges per minute as an average over a period of time to prevent overheating.
 - The simulator chest must be kept dry. Special attention should be taken when using IV Arm.
 - Do not apply conductive gel or conductive defibrillation pads intended for patient use to prevent chest skin pitting.
 - Do not use cables or connectors with visible damage.
 - Observe all normal safety precautions for use of defibrillators.
 - Defibrillation must be performed over the two defibrillation connectors only.

Air Tank Precautions: The right thigh of the ALS Simulator contains a compressed air tank. Contents under pressure:

- Do not puncture
- Do not store in excessive heat
- Release pressure in air tank prior to shipment

Do not pressurize above 120 psi or remove any safety devices to increase pressure. The tank is protected with safety devices not allowing the pressure to exceed 150psi. Do not disassemble or attempt to repair the internal parts of the right thigh. In the event of failure, contact the Laerdal Service Center immediately for instructions.

This device generates, uses and possibly radiates radio-frequency energy. If it is not installed and used in accordance with the instructions, it may cause harmful interference to radio communications. In that case, the user is encouraged to attempt correction of the interference by:

- Reorienting or relocating the receiving antenna
- Increasing the distance between the device and receiver
- Connecting the device to an outlet on a circuit different from that to which the receiver is connected
- Consulting the dealer or an experienced radio/TV technician for help

Electromagnetic radiation from other radio transmitters or other electronic equipment may cause noise in the head speaker. To eliminate this noise, move simulator away from the radiation source or turn the head speaker volume to zero.

Items Included:

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- (1) Adult, Full Body Simulator
- (1) Articulating IV Training Arm
- (1) Left Lower IO Leg
 - (1) IO Leg Skin
 - (1) Tibial IO Pad
 - (1) Tibial IO Bag Set
 - (1) Red Simulated Blood
- (1) Blood Pressure Training Arm
- (1) Simulator Lubricant
- (1) Interchangeable Pupil Set
- (6) Neck Skin Collars
- (1) Cricothyroid Membrane Tape Roll
- (1) Set of Replacement Pneumothorax Bladders
 - (2) Mid-Clavicular Bladders
 - (2) Mid Axillary Bladders
- (6) Chest Drain Modules
- (1) Blood Pressure Cuff
- (1) Air Pump
- (1) Set of Defibrillation Posts
- (1) Carry Case
- (1) Jacket
- (1) Pants
- (1) Directions for Use
- (1) Laerdal Global Warranty Booklet

Skills Taught:

Airway Management Skills

- Assessment of respirations
- Obstructed airway
 - Tongue Edema
 - Right Lung Blockage
 - Left Lung Blockage
 - Both Lung Blockage
- Endotracheal intubation
- Nasotracheal intubation
- Digital intubation
- Oropharyngeal airway insertion
- Nasopharyngeal airway insertion
- Bag Valve Mask
- Retrograde intubation
- Lightwand intubation
- Laryngeal Mask Airway insertion
- Laryngeal Tube insertion
- Combitube insertion
- Trans-Tracheal Jet Ventilation
- Surgical cricothyrotomy
- Needle cricothyrotomy
- Suctioning techniques
- Stomach auscultation to verify proper airway positioning
- Head Tilt/Chin Lift
- Modified Jaw thrust
- Ventilations

Drug Administration

- IV insertion, infusion and bolus into peripheral veins of forearm, antecubital fossa and dorsum of the hand
- Subcutaneous and intramuscular injections

Tibial I.O Capabilities

- Proximal tibial tuberosity (left)

Tension Pneumothorax Decompression

- Bilateral mid-clavicular sites
- Mid-axillary site (Right)

Chest Tube Insertion

- Mid-axillary site (Left)

Cardiac Related Skills

- 3 – 4 lead ECG with an extensive ECG library available for interpretation.
- Pacing with variable threshold, with or without capture
- Defibrillation capabilities (25 -360 joules)
- Programmable scenario base algorithms for instructor control
- Closed chest compressions

Circulatory Skills

- Bilateral carotid pulse
- Measurement of non-invasive blood pressure
 - Brachial and radial pulses
 - Auscultation or palpation of non-invasive blood pressure measurement
 - Values to be set for systolic, diastolic pressures, heart rate, auscultatory gap and volume

Sounds and Speech

Refer to SimPad Directions for Use or LLEAP Help Files

Logging and Scenario Functions

Refer to SimPad Directions for Use or LLEAP Help Files

Simulation Software

If you have purchased relevant licenses, you will have access to a number of PC applications that facilitate your simulation. To run a simulation, LLEAP (Laerdal Learning Application) must be started from Laerdal Simulation Home on the Instructor PC.

Laerdal Simulation Home

Laerdal Simulation Home is an application from where LLEAP and other Laerdal programs related to patient simulation can be found and started. Also, the help files can be opened from here. Laerdal Simulation Home is located in the Laerdal Medical folder under the Windows start menu (Windows 7).

Software used in a simulation session can be divided in the following main applications:

- LLEAP (Laerdal Learning Application)
- Voice Conference Application
- Patient Monitor
- SimView Server or Session Viewer

In addition SimDesigner and other applications are used for designing or preparing a simulation.

LLEAP

LLEAP is the instructor's application from where the simulation session is run, controlled, and monitored. LLEAP can be operated in Automatic or Manual mode. Automatic mode is used for pre-programmed scenarios while Manual mode allows the instructor full manual control over the simulation session. Running simulations in Manual mode generally requires some medical expertise to create clinically sound simulations.

Voice Conference Application (VCA)

The VCA software is needed to transmit all vocal sounds used during simulation. It allows the instructor to communicate through the simulator during the session. VCA can also be used to communicate with other instructors on a network, and create separate channels that only members can talk and listen to.

Patient Monitor

The Patient Monitor application emulates a typical hospital patient monitor. It is the learner's console and can be set up and controlled by the instructor, as well as by the learner, through on- screen touch menus.

Session Viewer and SimView Server

Session Viewer and SimView Server are applications that record video and patient monitor screen captures during simulation, in addition to providing an interface to debrief your session. After a session is ended, log files generated in LLEAP are transferred and merged with the video files in Session Viewer or SimView Server for the debrief.

Session Viewer typically runs locally on the same computer used for LLEAP, SimView Server runs on a dedicated server in the local network. During the first start-up of LLEAP you are prompted to select a debriefing system available on your computer or local network. This can be changed later.

Other Applications

There are also other programs that are used in conjunction with the simulation sessions, for example License Manager for handling program licenses and Simulator Firmware & Network Fixer for updating the firmware of the simulators or troubleshooting network problems.

SimDesigner

SimDesigner is an application for configuring your own pre-programmed scenarios. It can also be used to analyze and print out a graphical representation of a scenario.

SimDesigner must be installed to allow conversion of legacy instructor application files to LLEAP compatible file formats.

For a full overview of all applications and their help files, start LLEAP Home.

Web Downloads

Visit www.laerdal.com/downloads to download the latest User Guide and Software.

Unpack

The ALS Simulator is packaged without the legs attached. Remove the upper body and legs from the packaging and attach the legs to the torso.

Attaching the Legs:

To attach the Right Leg (with enclosed air tank):

1. Remove bladder reservoir from pelvis.
2. Remove blank genitalia module from pelvis.
3. Carefully slide the tubing and wiring exiting the fixed bolt from the upper thigh of the right leg through the connector hole and into internal cavity of pelvic area. (Photo 1)


Photo 1

4. Continue to advance the bolt and wiring through the connector hole until the thigh is flush with the pelvis of the simulator.
5. From inside pelvic area, thread the longer air connection tube first through the provided washer, spring, and thumb nut in that order. Then, thread the shorter electronic wire through the washer, spring, and thumb nut.
6. Connect the wire and cable inside the pelvis.
7. Tighten the thumb nut until desired articulation is achieved.
8. Replace bladder reservoir in torso.
9. Replace blank genitalia module in pelvic region of torso.

To attach the Left Leg:

1. Remove left thigh injection pad (A), genitalia module (B) and bladder reservoir assembly (C)

2. Select and prepare hardware (Photo 2), pass spring and small washer over bolt (Photo 3)
3. To install pass bolt assembly through thigh, then pelvic opening from outside inward (Photo 4), select large fender washer and pass over bolt from inside simulator abdominal cavity. Screw wing nut onto bolt (Photo 5)
4. Replace bladder reservoir and genitalia module in pelvis and injection pad in thigh. (See Step 1)

NOTE: Tighten wing nut as desired to simulate leg range of motion.

SimPad and LLEAP Connections

1. Connect ALS Simulator to Link Box (Photo 6), via cable located on lower right side of simulator.
2. Connect the clear tubing exiting the simulators lower right side to the Link Box
3. Attach blood pressure cuff to simulator's left arm.
4. Attach clear pneumatic tubing from blood pressure cuff (Photo 7) to corresponding clear pneumatic tubing exiting the simulator's left shoulder.

Photo 6

Photo 7

Getting Started

Cricothyroid Membrane/Neck Skin Collar:

1. Cut a 5 cm (2 inch) strip of Cricothyroid Membrane Tape.
2. Adhere tape to edges of cricoid opening, covering opening with tape, to create actual membrane.
3. Lay a Neck Skin Collar into molded track around neck area of simulator. (Figure 1)
4. Attach collar using Velcro® strips.


Fig. 1

Laerdal Recommends

A tight seal will enhance actual feel and sound of penetrating cricothyroid membrane, and facilitate full rise and fall of the chest during bag valve mask ventilation.

Pupil Replacement

ALS Simulator is delivered with a set of normal pupils mounted in the head. A separate case contains 3 sets of plastic pupil inserts (normal, constricted and dilated pupils) for use in simulating other conditions.

To change the pupils in the airway head:

1. Open the eyelids wide, take care not to rip the face skin.
2. Using the suction cup provided or the edge of your fingernail, remove the pupil from the eye.
3. Replace with the pupil of choice, using the small suction cup tool provided or press into place with finger.

Filling the Thigh Air Tank

1. Attach the provided air pump to the Schrader valve located laterally on the posterior of the right thigh. (Photo 8)
2. Fill the air container using the air pump to 120 psi. **NOTE:** An over pressure

Photo 8

- release valve is included inside the thigh and will be activated at pressure above 10 Bar (150 psi). This valve has a built in reset mechanism.
3. Detach air pump from the Schrader valve.

NOTE: Most products designed for automobile/bike tire inflation will be sufficient to fill the thigh tank. i.e. electric pumps, air compressors, air tanks, etc. Pressure limitations vary with each unit and should be considered prior to purchase.

In Use

Laerdal Recommends

- Endotracheal tube - size 7,5
- LMA # 4
- Combitube Large Adult or Trainer
- KING LT – # 4

Airway Management:

1. Spray inside of pharynx, nostrils and all intubation tubes with a liberal amount of provided simulator lubricant or liquid soap. Sit the simulator upright to allow lubricant to coat passages.
2. Perform further intubation procedures as per your local training protocol.

NOTE: Should incorrect ventilation be performed, air will pass through esophagus, causing distention of the abdomen.

Obstructed Airway

Obstructed Airway (Manual Inflation of Tongue):

1. Locate inflation bulb at lower left side of torso labeled "tongue."
2. Tighten air release knob.
3. Pump bulb several times to inflate tongue to desired size.

To deflate tongue, loosen air release knob.

Obstructed Airway (Right/Left/Both Lung Blockages):

See SimPad Directions for Use or LLEAP Help Files for complete airway blockage functionality and operating information.

Tension Pneumothorax Decompression (Manual Inflation of Pneumothorax):

1. Locate inflation bulb which exits lower left side of simulator; labeled "pneumo."
2. Tighten air release knob so that as you pump the bulb, no air escapes.
3. Pump bulb several times to prepare three (3) sites for decompression.
4. Decompress the site, bilateral mid-clavicular or right side mid-axillary or per your local training protocol.

When the site is "decompressed," air can be heard/observed passing through catheter.

IV Cannulation:**Laerdal Recommends**

Use a 21G or smaller catheter for IV cannulation to extend the life of your simulator's IV Arm. When using a catheter or other venipuncture device, spray with lubricant spray for easier insertion.

NOTE: Do not use iodine as this will permanently stain the simulator.

1. Attach an IV bag with tubing. **NOTE:** Make sure the clamp is closed.
2. Use a syringe to inject blood concentrate and 500cc of water into an IV bag with tubing. **NOTE:** Make sure the clamp is closed. This will serve as the "blood source" bag. (Photo 9)

Photo 9

3. Attach the "blood source" bag with IV tubing to one of the latex tubes exiting the simulator's right IV Arm.
4. Attach an empty bag with IV tubing to the second latex tube exiting the simulator's IV Arm. This bag will serve as the collection reservoir. (Photo 10)
5. Control flow of blood from arm, via clamp on collection reservoir tubing. Put collection bag on floor to allow the now closed system to gravity feed.
6. Hang "blood source" bag from IV pole and open clamp to allow concentrate to flow through arm.

NOTE: A blood pressure cuff can be used as a pressure infuser. This will allow the blood bag to be placed on the floor or between the mattress and frame. (Photo 10)


Photo 10

7. When the "blood source" bag is empty, switch the "blood source" bag and reservoir: Hang the now full reservoir and place empty reservoir bag on floor to reverse flow of concentrate.

NOTE: Be sure to adjust the flow regulator. The collection reservoir clamp should always be used to control flow rate.

The arm is now ready for cannulation. Train as per your local protocol.

Intramuscular Injections:

Medications can be administered via intramuscular injections bilaterally in the deltoids and thighs and in the gluteal and ventro-gluteal locations.

Laerdal Recommends

Medications may be administered via local protocol. All precautions and safety measures must be followed during training.

Urinary Catheterization (with genitalia upgrade kit):

1. Remove blank genitalia from the simulator.
 2. Lift and rotate the pelvic pin to one side.
- NOTE:** Failure to replace the pelvic pin will result in pelvic spread.
3. Lift the urinary reservoir.
 4. Using the urinary valves provided, attach threaded end of the connector valve to the corresponding threaded urinary connector on the genitalia (Photo 11-13).
 5. Place the C Clamp over the valve, leaving loose.
 6. Install the valve onto the black connector of the reservoir. Ensure that it sits flush to the reservoir.

Photo 11 (Clamp not yet in place.)

7. Move the C Clamp over the connector and the valve.

Photo 12 (Area where clamp is placed.)

8. Squeeze to tighten the C Clamp. This will lock the C Clamp in place.

Photo 13 (Apply clamp and squeeze to tighten.)

9. Replace the genitalia assembly into the simulator torso.
10. Replace Pelvic Pin.

NOTE: Failure to replace the pelvic pin will result in pelvic spread.

11. Remove plug in urinary reservoir.
12. Using a 100cc syringe, fill reservoir with colored water.
13. Perform urinary catheterization as per your local training protocol.

Simulator Functionalities

Cardiac Related Skills:

ALS Simulator, when used with SimPad, features an extensive library of ECG variations. External pacing with or without capture, defibrillation and synchronized cardioversion can also be performed. **See SimPad Directions for Use or LLEAP Help Files**

NOTE: Prior to training, ensure the simulator cable exiting the simulator's right side is attached to the Link Box.

1. Compressions may be performed with one or two hands.
2. For monitoring purposes, lead placement snaps are located on simulator as follows:
 - i. RA (Right Arm)
 - ii. LA (Left Arm)
 - iii. LL (Left Leg)
 - iv. RL (Right Leg)

NOTE: If the ECG snaps/defibrillation connectors do not attach to the clinical monitor/defibrillator/external pacer used by your institution, additional ECG snaps and or defibrillation/pacing cables are available for purchase. **Please contact your representative for correct adaptor.**

3. Two Zap Plates (discs) and two hands free adapters have been designed for use with defibrillators. Screw discs or hands free adaptors into post sockets located on apex and sternum of simulator. Rhythm may be read from simulator using a monitor by applying paddles to discs with firm pressure. **(Do not use gel or hands free self adhesive pads.)**
4. For manual defibrillation, place paddles on zap plates and press down firmly for good contact. **NOTE:** Monophasic and Biphasic defibrillators may be used.

Laerdal Recommends

Caution: Follow defibrillation protocol by avoiding contact between the paddle and any of the electrode sites while defibrillating. A load of up to **360** joules can be delivered.

See SimPad Directions for Use or LLEAP Help Files for complete cardiac functionality and operating information.

Safety Precautions used during Defibrillation of a Patient/Trainer

1. Read and follow all safety and operating instructions provided with your defibrillator and associated equipment.
2. This trainer can be shocked with actual voltages and current used on a live patient. All precautions and safety measures must be used during the defibrillation and pacing phases of training. Failure to follow safety measures could result in injury or death to operators, students or onlookers not heeding these warnings.

Blood Pressure (BP) Arm:

Left Arm Installation

The simulator is packaged with the left BP Arm attached. The BP Arm is designed to rotate approximately 220°.

Caution: To prevent damage, do not over rotate the left BP Arm.


Photo 14

Photo 15


Photo 16

Photo 17 – Inside View
Torso Left Side

Photo 18 – Outside View Left
Shoulder Standard or Wound Arm

To detach the BP Arm:

1. Detach the chest skin, lift the chest plate and push it aside to remove the stopper screw (photo 14) and release the BP Arm.
2. Take out the arm stopper, disconnect the BP cables and remove the BP Arm. **NOTE:** Properly store the BP Arm and the attachment accessories for future use.

To install the Left Arm (standard or wound model):

1. The simulator includes the necessary arm hardware (photo 15) Select the replacement arm, remove the deltoid pad
2. Insert the arm adapter (photo 16) located inside the torso on the left side
3. Pass the small washer over the bolt, thread the bolt and washer assembly through the openings of the arm shoulder, torso and arm adapter (photo 17 & 18) **NOTE:** Tighten the wing nut as desired to simulate range of motion
4. Replace the chest plate, the chest skin and the shoulder deltoid pad.

To re-attach the BP Arm, reverse the processes. Properly store the arm, the adapter and the hardware set for future use.

Auscultating and palpating blood pressure:

1. Place the supplied Blood Pressure cuff on the BP Arm (left).
2. Attach clear tubing on cuff to matching clear pneumatic tubing exiting underneath the simulator's left arm. (Photo 19)

Photo 19

3. Ensure that the simulator cable and the clear pneumatic tubing, both exiting the simulator's right side are attached to the Link box.

Korotkoff sound volume can be controlled, an auscultatory gap heard and pressure set in steps of 2 mmHg.

See SimPad Directions for Use or LLEAP Help Files for complete Blood Pressure Arm functionality and operating information.

Laerdal Recommends

Do not insert needles into the BP Arm.

Heart, Breath and Bowel Sounds and Speech

The ALS Simulator has heart, breath and bowel sounds, both normal and abnormal, with adjustable instructor controlled volume.

Ensure that the ALS Simulator is connected to the Link Box via the simulator cable located on lower left side of simulator.

See SimPad Directions for Use or LLEAP Help Files for complete Heart, Breath and Bowel Sound Auscultation and Speech functionality and operating information.

Spontaneous Breathing and Airway Closures

The ALS Simulator has spontaneous breathing (visible chest rise and fall) with variable breathing rate. The breathing is generated via an enclosed air tank in the right thigh that is pressurized using a standard air pump. (See "Filling the High Air Tank" located in the "Getting Started" section of the DFU for more information).

The left lung and right lung can be closed independently or together to create a partial or complete airway obstruction.

See SimPad Directions for Use or LLEAP Help Files for complete breathing and airway blockage functionality and operating information.

Maintenance:**Cricothyroid Membrane/Neck Skin Collar:**

1. Remove neck skin collar, by detaching Velcro®.
2. Remove used Cricothyroid Membrane Tape from cricoid opening.
3. Cut a 2 inch/ 5 cm strip of Cricothyroid Membrane Tape.

Fig. 2

4. Adhere tape to edges of cricoid opening, covering the opening with tape, to create the actual membrane.
5. Lay a Neck Skin Collar into molded track around neck area of the simulator. (Figure 2) **NOTE:** If the used collar is still in good repair simply slide the collar in either direction for a fresh site. If the collar is showing wear, discard and replace with a new one.
6. Attach collar using Velcro® strips.

Laerdal Recommends

The neck skin collar is designed to provide multiple sites for needle and surgical techniques. When a fresh site is needed, slide the collar in either direction. If the neck skin collar is showing wear, discard and replace with a new collar.

IV Cannulation**After IV training is complete:**

1. Close both clamps on IV tubing and detach bag tubing from arm tubing. Fill a 100cc syringe with water. Flush IV Arm with syringe. **NOTE:** This should be done over a basin to catch any fluid.

2. Using the empty syringe, push air through venous system to remove excess water.

Intramuscular Injections

After intramuscular injection training is complete:

1. Remove sponge from inside of injection pad skin. **NOTE:** Do not store wet foam pads in skin. This will cause shrinkage of injection pad skin.
2. Rinse sponge with water, allowing it to dry overnight.
3. Insert sponge back into injection pad skin. **NOTE:** A small amount of powder applied to sponge will ease insertion.
4. Replace the pad in the simulator.

Urinary Catheterization (with genitalia upgrade kit):

After urinary catheterization is complete:

1. Remove complete genitalia assembly
2. Disassemble the genitalia assembly
3. Working over a basin empty the reservoir
4. Rinse parts and reservoirs thoroughly and allow to dry overnight

Pneumothorax:**(A) Bilateral Mid Clavicular Sites:**

1. Detach the chest skin from torso from tabs at the shoulder and back and remove.
2. Remove the hard chest plate from the torso.
3. Disconnect bladder hose from Y connector, located on the underside of the chest plate. (Photo 20)

Photo 20

4. Remove the bladder by pulling it through the opening between the 2nd and 3rd intercostal spaces on the top side of the hard chest plate. (Photo 21)

Photo 21

5. Trim bladder tubing to match original tubing length to fit this simulator style.
6. Insert the new pneumothorax bladder into the top side of the chest plate through the 2nd and 3rd intercostal space opening with the bladder tubing exiting through the back side of the chest plate. **NOTE:** The narrow edge of the bladder should be inserted closest to the sternum. (Photo 22)

Photo 22

7. Reconnect hose to Y hose connector.
8. Return chest plate to proper position on torso.
9. Replace Chest Skin over torso, being sure to secure skin at shoulders and back.
10. Fill puncture marks on the exterior of the Chest Skin with the wax included in the Bladder Replacement Kit.

(B) Mid Axillary Site (RIGHT):

1. Remove Chest Skin from tabs at the shoulder and back.
2. Remove pneumo bladder box from right side of torso. (Figure 3)


Fig. 3

3. Disconnect bladder hose from inline hose connector, being sure hose does not fall through hole and into torso.
4. Remove bladder from pneumo pad and discard.
5. Trim bladder tubing to match original tubing length to fit this simulator style. Connect new bladder hose to inline hose connector.
6. Fold and insert new pneumo bladder into pneumo box.
7. Replace Chest Skin over torso, being sure to secure skin at shoulder and back tabs.
8. Fill puncture marks on exterior of Chest Skin with the modeling wax supplied in the Bladder Replacement Kit.

NOTE: Bladders can be repaired using vulcanized or contact cement glue (not provided). Cover punctured surface of pneumothorax bladders with glue. Allow bladders to air dry completely overnight before replacing in the chest plate.

Chest Drain:

Mid Axillary Site (LEFT)

1. Remove Chest Skin from tabs at shoulders and back.
2. Remove the Chest Tube Insertion Module from the left mid-axillary site of the simulator. (Photo 23)


Photo 23

3. Insert a new Chest Tube Insertion Module.
4. Replace Chest Skin over the torso being sure to secure the skin at the shoulders and both sides.

Attaching/Removing Arms:

1. Remove Deltoid Injection Pad from upper arm.
2. Detach Chest Skin from tabs at shoulders and back.
3. Remove skin to reveal internal portion of upper chest.
4. Lift hard chest plate to reveal inside sockets for arm connections.
5. Place one washer on threaded 4 1/2" bolt.
6. Insert threaded bolt through upper arm and through torso so that it is visible in chest cavity. (Figure 4)


Fig. 4

7. From inside of chest area, slip a washer, a spring and another washer on bolt.
8. Screw a wing nut on bolt and tighten until desired articulation is achieved.
9. Replace hard chest plate.
10. Pull Chest Skin back onto shoulder area.
11. Connect Chest Skin to torso at shoulder and back tabs.
12. Replace Deltoid Injection Pad on upper arm.

To remove arms(s), reverse procedure.

IV Arm Skin and Vein Replacement

Laerdal Recommends

When excessive leaking occurs at the puncture sites, a new vein system and skin should be installed to reduce loss of fluid.

We recommend working at a basin when replacing skin and veins.

Replacing Skin and Veins:

1. Cut off skin. This can be done with a sharp knife or scalpel. (Figure 5)

Fig. 5

2. Discard skin.
3. Remove tubing from track in mandrel. Glue may need to be scraped away.
4. Rinse and dry vein grooves well and swab with alcohol. Be sure to remove any excess glue.
5. Place new veins along grooves, (Figure 6) spot gluing as needed. (We recommend a fast-drying glue applied sparingly every 2 to 3 inches.)


Fig. 6

6. Generously lather arm mandrel with liquid soap.
7. Slide hand into skin. (Photo 24)


Photo 24

8. Work skin over fingers, as with a glove. (Photo 25)

Photo 25

9. Work arm skin up, over mandrel. (Photos 26 & 27.)

Photo 26

Photo 27

After Use:

Cleaning

1. Clean with mild soap and water; do not submerge the simulator or parts in cleaning fluids or water.
2. Use only on clean surface. Avoid felt tipped markers, ink pens, acetone, iodine or other staining products and avoid placing the simulator on newsprint or inked lines of any kind.
3. To ensure longevity, each simulator should be cleaned after each training session and a general inspection should be conducted regularly.
4. Modules and all other parts should be drained and air-dried thoroughly before storage and disinfected when needed. After use of injection pads (use water only), accumulated water should be squeezed out. **NOTE:** Do not store wet foam pads in the skin. To prevent mildew or mold, pads can be soaked in a mild solution of disinfectant and water or bleach and water. Squeeze excess solution from pads, allow them to dry, then store or reinsert in simulator.
5. Articulating parts will benefit from a light application of talcum powder prior to training sessions.

Extend Cleanliness

To extend the cleanliness of your simulator:

1. Use a small towel and lightly dust the simulator head and chest with a small amount of powder (provided). Wipe away any excess.
2. Use gloves when handling the simulator.

Laerdal Recommends

Use only on clean surface. Avoid felt tipped markers, ink pens, acetone, iodine or other staining products and avoid placing the simulator on newsprint or inked lines of any kind. These will stain the simulator.

Technical Data

General

Operating Temperature:

0°C – 40°C at 90% relative humidity, non-condensing

Storage Temperature:

-15°C – 50°C at 90% relative humidity, non-condensing

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Defibrillation: Average of 720 J /minute max

Emergency Cardiac Capabilities

- Synchronized variable heart rate, rhythm, abnormalities and duration
- Defibrillation (25 -360 J)

IV Arm

Accessible veins include median, basilic and cephalic

Standards/Approvals

The product is in compliance with the essential requirements of Council Directive 2004/108/EC on electromagnetic compatibility (EMC).

The product is in compliance with Council Directive 2011/65/EU on restriction of the use of certain hazardous substances (RoHS).

Minimum Computer Requirements

- Intel i-core generation 3 or newer
- Score over 3000 PassMark - CPU Mark
- 4 GB RAM
- 120GB hard disk space
- I366x768 or better
- Headset w/microphone
- USB ports – 2

Minimum Software Requirements

- Windows 7 or Windows 8
- 100% DPI

Replacement Parts:

For latest version of Spare Parts and Accessories, visit www.laerdal.com

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- (1) Adult, Full Body Simulator
- (1) Articulating IV Training Arm
- (1) Left Lower IO Leg
 - (1) IO Leg Skin
 - (1) Tibial IO Pad
 - (1) Tibial IO Bag Set
 - (1) Red Simulated Blood
- (1) Blood Pressure Training Arm
- (1) Simulator Lubricant
- (1) Interchangeable Pupil Set
- (6) Neck Skin Collars
- (1) Cricothyroid Membrane Tape Roll
- (1) Set of Replacement Pneumothorax Bladders
 - (2) Mid-Clavicular Bladders
 - (2) Mid Axillary Bladders
- (6) Chest Drain Modules
- (1) Blood Pressure Cuff
- (1) Air Pump
- (1) Set of Defibrillation Posts
- (1) Carry Case
- (1) Jacket
- (1) Pants
- (1) 사용설명서
- (1) Laerdal 글로벌 제한보증 책자

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© 2014 Laerdal Medical AS. All rights reserved.
Manufacturer: Laerdal Medical Corporation