
- EN
- FR
- DE
- ES
- IT
- BR
- NL
- PL
- JA
- ZH
- KO

SimBaby

Directions for Use

INTRODUCTION	4
SimBaby Patient Simulator	4
Items Included	4
Link Box	4
CAUTIONS AND WARNINGS	5
General Simulator Handling	5
Antivirus and Firewalls	6
File Security and Backup of Data	6
Regulatory Information	6
SPECIFICATIONS	7
FEATURES	8
SimBaby Overview	8
General	8
Laerdal Simulation Software	9
Blood System	9
Airway	9
Breathing	10
Air and CO ₂ Source	10
Circulation	10
Sounds	11
Body Movement	11
Drugs and IV	11
SETUP	11
Circulation	11
Connecting IV Arm	11
Connecting IV/IO Legs	11
Filling the IV Bag	12
Connecting the SpO ₂ Probe	12
Connecting to the Link Box	12
MAINTENANCE	13
Daily Maintenance	13
Replacing Tension Pneumothorax Bladder	13
Replacing Chest Tube Module	14
Replacing Esophagus Filter	14
Replacing Pupils	14
Replacing Skin and Veins in IV Arm	15
Replacing IV/IO Leg	15
SPARE PARTS AND ACCESSORIES	17

SimBaby Patient Simulator

For Clinical Simulation

SimBaby is a patient simulation system that facilitates training of Basic and Advanced Life Support. This system allows the instructor to effectively assess the learner's individual and team skills based on realistic clinical situations.

SimBaby allows observation and recognition of most vital signs. This is achieved through direct interaction with the Patient Simulator and observation of the Patient Simulator's status which can also be viewed on the Patient Monitor PC (optional).

The SimBaby Simulation System Features Include:

- An advanced configurable airway which allows difficult airway management skills to be simulated.
- Eyes with interchangeable pupils.
- The instructor can modify both the lung compliance and the airway resistance.
- A range of pre-programmed vocal sounds like crying and coughing to support the simulation scenario.
- Vascular access routes include an IV Arm and on both legs. IO access is also available on both legs.
- A wide range of spontaneous breathing patterns including variable rate and depth, and complications like subcostal retractions and see-saw breathing.
- Automatic Simulation Control based on pre-programmed and validated Patient Cases.

The Main Components of the SimBaby

SimBaby is a 6 months old, life-size Patient Simulator. Interventions by learners are registered by the system in the session log, and used for later debriefing.

The Instructor PC controls the simulation. The headset allows the instructor to simulate interactive voice communication between the patient and the learner.

Software available for use with the simulator includes LLEAP for controlling scenarios, SimDesigner for creating and editing scenarios, SimView Server or SessionViewer for debriefing simulation sessions with video capture and a Patient Monitor application.

Items Included

The Laerdal SimBaby includes the following main components:

- LLEAP (Laerdal Learning Application) Software License
- SpO₂ Pulse Oximetry Probe
- Set of plates for Manual Defibrillation
- Airway Lubricant
- Esophagus filter assy (50)
- Pneumothorax Kit, consisting of:
 - Extra Pneumothorax Bladder
 - Extra Pleura Sleeves (10)
 - Extra Chest Skin

- IV/IO Kit, consisting of:
 - Tubing, arm
 - Extra arm skin
 - Lower IV/IO leg right
 - Lower IV/IO leg left
 - Simulated Blood Concentrate
 - Blood system
- Talcum Powder
- Blood Pressure Measuring Kit

Additional items required for use:

- (configurations vary)
- Laptop computer
 - Simulated Patient Monitor
 - Link Box

Connection Cables and Tubing, consisting of:

- Power supply, Simulated Patient Monitor
- AC Line Cord, Link Box
- Cable, Link Box to Simulator (15-pin)
- Cables, PC to Simulated Patient Monitor (signal and USB)
- Cable, Link Box to PC (9-pin)
- Cable Audio, PC to Link Box and Patient Monitor
- Tube, Simulator Air and CO₂ Supply

Pressurized Air and CO₂ source (to be purchased separately if not included):

- Compressor Unit (included in some configurations) or Regulator Unit.

Link Box

SimBaby operates with Link Box (SimMan version). Link Box connects the simulator to the computer. See *Setup* Section for more information.

Disclaimer

Use of the Patient Simulator system to train personnel should be undertaken under supervision of suitably trained medical personnel with an understanding of educational principles as well as recognized medical protocols.

As with all Patient Simulators or other such training devices there may be approximations, variations and inaccuracies in anatomical features and the physiological modeling. This being the case, Laerdal Medical does not guarantee that all features are completely accurate.

⚠ Cautions and Warnings

A **Caution** identifies conditions, hazards, or unsafe practices that can result in minor personal injury or damage to the Simulator.

A **Warning** identifies conditions, hazards, or unsafe practices that can result in serious personal injury or death.

General Simulator Handling

⚠ Warnings

- Do not ventilate the Patient Simulator with oxygen enriched air or flammable gases.
- Only use fluids as directed in the *Directions for Use*, otherwise this may damage the Patient Simulator and its components.
- Do not introduce humidified air into the system during ventilation.
- Never perform mouth-to-mouth or mouth-to-nose rescue breathing on the Patient Simulator. The simulator's airways are not designed for cleaning or disinfection.
- Do not use the Patient Simulator if the internal tubing and cabling is disconnected.
- Never use the Patient Simulator in temperatures exceeding 40 °C (104 °F), as this may cause overheating and shut down.
- Never use the Patient Simulator in temperatures below 10 °C (50 °F).
- The Patient Simulator should never be stored in temperatures below -15 °C (5 °F) or over 50 °C (122 °F).
- Using a defibrillator in temperatures over 35 °C (95 °F) may cause overheating and shut down.
- Avoid pinch hazards - Do not remove protective bushings from the Patient Simulator's joints or use it without the external skins.
- Avoid all sharp edges on the Patient Simulator to prevent personal injury.

Transportation and Storage

Ensure that the Patient Simulator is properly secured during transportation to prevent personal injury or damage to the product.

Defibrillation Hazards

A conventional defibrillator may be used on SimBaby. During live defibrillation, the defibrillator and Patient Simulator may present a shock hazard. All standard safety precautions must be taken when using a defibrillator on the Patient Simulator. For more information, consult your defibrillator's User Guide.

⚠ Warnings

- The ECG connectors are designed exclusively for ECG monitoring and must not be used for defibrillation. Defibrillation on the ECG connectors will damage the internal electronics of the Patient Simulator and may cause personal injury.
- Defibrillation must be performed on the defibrillator connectors only.
- The Patient Simulator must not come into contact with electrically conductive surfaces or objects during defibrillation.
- Do not defibrillate the Patient Simulator in a flammable or oxygen enriched atmosphere.
- The Patient Simulator torso must always be kept dry. Allow the Patient Simulator to acclimate before defibrillating. Sudden changes in temperature (moving the Patient Simulator from a cold environment to a warm environment and vice versa) may result in condensation collecting on the base board and pose a shock hazard.
- If defibrillation is performed over one or more of the ECG connectors, high voltages may be present on the remaining connectors during the shock. Defibrillation attempts via the ECG connectors will also damage the internal electronics requiring that they have to be replaced.
- To prevent overheating, do not provide more than three (3) defibrillator discharges (max 360J) in a sequence. Do not exceed an average of two (2) defibrillator discharges per minute during the training session.
- To prevent chest skin electrode pitting, do not apply conductive gel or conductive defibrillation pads intended for patient use.
- Do not use cables or connectors which have visible damage.

General Use

⚠ Cautions

- Do not introduce any fluids (except airway lubricant in small amounts to lubricate the airway adjuncts) into the esophagus or trachea of the simulator.
- After each session involving airway lubricant the esophagus filter must be replaced. See "Maintenance" for more information.
- If a training session involves the administration of fluids and/or drugs into the IV Arm and IO/IV Legs, empty the arm and legs immediately following the training session.
- Electromagnetic radiation from other radio transmitters or other electronic equipment may cause noise in the head speaker. To eliminate this noise move simulator away from the radiation source or turn the head speaker volume to zero.

⚠ Warning: The veins in the IV Arm, lower legs and Pneumothorax Bladder contain latex. Users who suffer from latex allergy should take precautions while using or handling the latex parts by wearing non-latex protective gloves.

⚠ Do not use the Patient Simulator if:

- Limbs are not attached to the torso
- Skins are torn or not properly fastened
- Internal or external cables, tubes or connectors are damaged
- There is fluid leakage inside the Patient Simulator torso
- There are unusual sounds indicating air leakage or mechanical damage
- There are signs of electrical malfunction, such as an unresponsive Patient Simulator, an unusual smell or smoke

Hygiene

- To maintain Patient Simulator skins, wash hands before use and place the Patient Simulator on a clean surface.
- Wear gloves as required during simulation scenarios.
- Use only *Laerdal Airway Lubricant* in the Patient Simulator's airway. Do not spray airway lubricant into the Patient Simulator. Lubricate airway adjuncts only.

Prevent Stains on Patient Simulator Skins

Avoid using colored plastic gloves, as they may cause discoloration of the Patient Simulator skin.

Do not use felt-tipped markers, ink pens, acetone, iodine or other staining medications near the Patient Simulator. Take care not to place the Patient Simulator on newsprint or colored paper. Staining may be permanent.

Latex

The veins in the IV Arm and in the lower legs and Pneumothorax Bladder contain latex. Users who suffer from latex allergy should take precautions while using or handling the latex parts by wearing non-latex protective gloves.

Antivirus and Firewalls

The Patient Simulator and PCs are not supplied with Antivirus programs. Windows firewall is activated by default. It is the customer's responsibility to protect the simulation system components from unauthorized access.

The Patient Simulator will revert to factory settings each time the power is switched off.

The customer should install all recommended Windows updates from Microsoft. General security measures should be taken before browsing the internet.

It is recommended that the Instructor PCs are used only as Patient Simulator controllers. Downloading other Software programs onto these machines may introduce unexpected errors.

File Security and Backup of Data

The customer is responsible for file security and backup routines for all simulation session data. All use and storage of simulation session data should be in accordance with local rules, regulations or laws and is the sole responsibility of the customer.

Regulatory Information

CE The product is in compliance with the essential requirements of Council Directive 2014/30/EU on electromagnetic compatibility (EMC); and with the essential requirements of Council Directive 2014/35/EU relating to electrical equipment designed for use within certain voltage limits (LVD); and with Council Directive 2011/65/EU on restriction of the use of certain hazardous substances (RoHS).

Waste Handling

 Dispose of in accordance with local requirements and regulations.

Link Box

 Link Box is in compliance with applicable U.S and Canadian safety standards and has been certified by Canadian Standards Association.

Specifications

Size and Weight

Dimensions	650 mm x 23 mm x 160 mm (25.59 in x 0.91 in x 6.30 in)
Weight	4 kg (8.8 lbs)

Power Consumption

Link Box and simulator Input	110/230V AC 1.4A
Link Box and Simulator Output	12V DC 5.0A
Compressor Unit Input	110V AC 1.9A
Compressor Unit Output	230V AC 1.0A

Air & CO₂ Pressure

External air connection	Max 1.4 bar
External CO ₂ to Patient Simulator	Max 1.4 bar

Temperature Limits

Operating temperatures	10 °C to 40 °C (50 °F to 104 °F)
Storage temperatures	-15 °C to 50 °C (5 °F to 122 °F)

Environment - Patient Simulator Only

Relative humidity	15% to 90% (non-condensing)
For indoor use only	

Material Chart for Patient Simulator

Clothes	Cotton, Nylon
Skins and airways	PVC (DEHP free)
External hard plastics	PP, PA, PC, PC/PET
Inner plastics	Silicone, TPU, TPE, PVC, Nitrile PA, PA+GF, PC, ABS, POM, HDPE, PET
Metal components	Aluminum, Brass, Steel

Recommended Hardware Specifications

Intel i-core generation 3 or newer
Score over 3000 PassMark - CPU Mark
4 GB RAM
120 GB hard disk space
1366x768 or better
Headset w/microphone
USB ports – 2

Minimum Software Requirements

Windows 7 or Windows 8
100% DPI

Acceptable Fluids for Patient Simulator

To clean the simulator's fluid and blood system use one of the following:	
Distilled water or De-ionized water	
Cleaning Fluids	60% - 70% Isopropanol

Simulated IV Fluids


Only use purified water to simulate IV fluids:	
Simulated IV fluids	- Distilled Water - De-ionized Water

Link Box

Dimensions	2.4 kg (5.3 lbs) 210 mm x 100 mm x 220 mm (8.27 in x 3.94 in x 8.66 in)
Operating temperature	10 °C to 40 °C (50 °F to 104 °F)
Storage temperature	-15 °C to 50 °C (5 °F to 122 °F)
Humidity	15% to 90% RH (non-condensing)

SimBaby Overview

- Bulging Fontanel
- ECG Connectors
- ⚡ Defib Connectors
- Pulses

General

Main anatomical features

- Bulging fontanel can be simulated.
- Interchangeable pupils for simulation of three different pupil sizes: contracted, normal and dilated.
- Anatomically modeled from live babies, the torso demonstrates normal anatomical surface landmarks.
- SimBaby has correct anatomical landmarks for external chest compressions. Chest compressions will produce compression artifacts on the Simulated Patient Monitor and palpable pulse.
- Simulator is equipped with separate defibrillation and ECG monitoring connectors. The system provides real 3-lead ECG readings (4 connectors).
- Head tilt and jaw thrust maneuvers can be performed.

Mobility of Joints

Neck:	Movement can be on a 3-axis movement of head. Range of movement can be restricted. See <i>Airway Features</i> section
Shoulders:	3-axis rotation
Lumbar:	1-axis
Elbows:	Fixed, no mobility
Wrists:	Fixed, no mobility
Thumbs:	Fixed, no mobility
Hip joints:	3-axis rotation
Knees:	1-axis rotation
Ankles:	Fixed, no mobility

Note: Do not remove protective bushings at shoulder or lower back. These are present to protect users from pinch points.

Laerdal Simulation Software

To run a simulation, LLEAP (Laerdal Learning Application) must be started from Laerdal Simulation Home on the Instructor PC.

Laerdal Simulation Home

Laerdal Simulation Home is an application from where LLEAP and other Laerdal programs related to patient simulation can be found and started. The help files can also be opened from here. Laerdal Simulation Home is located in the Laerdal Medical folder under the Windows start menu (Windows 7).

Software used in a simulation session comprises the following main applications:

- LLEAP (Laerdal Learning Application)
- Voice Conference Application
- Patient Monitor
- SimView Server or Session Viewer

In addition SimDesigner and other applications are also used for designing or preparing a simulation.

LLEAP

LLEAP is the instructor's application from where the simulation session is run, controlled and monitored. LLEAP can be operated in Automatic or Manual mode. Automatic mode is used for pre-programmed scenarios while Manual mode allows the instructor full manual control of the simulation session. Running simulations in Manual mode requires some medical expertise to create clinically sound simulations.

Voice Conference Application (VCA)

The VCA Software allows the instructor to communicate through the simulator during the session. VCA can also be used to communicate with other instructors on a network, and create separate channels where only members can communicate.

Patient Monitor

The Patient Monitor application emulates a typical hospital patient monitor. It is the learner's console and can be set up and controlled by the instructor, as well as by the learner through on-screen touch menus.

Session Viewer And SimView Server

Session Viewer and SimView Server are applications that record video and patient monitor screen captures during simulation, in addition to providing an interface to debrief your session. After a session is ended, log files generated in LLEAP are transferred and merged with the video files in Session Viewer or SimView Server for the debrief.

Session Viewer typically runs locally on the same computer as used for LLEAP. SimView Server runs on a dedicated server in the local network. During the first start-up of LLEAP, you are prompted to select a debriefing system available on your computer or on a local network. This can be changed later.

Other Applications

There are also other programs that are used in conjunction with the simulation sessions, for example License Manager for handling program licenses and Simulator Firmware & Network Fixer for updating the firmware of the simulators or troubleshooting network problems.

SimDesigner

The SimDesigner application allows you to configure your own pre-programmed scenarios. It can also be used to analyze and print out a graphical representation of a scenario.

SimDesigner must be installed to allow conversion of legacy instructor application files to LLEAP compatible file formats.

For a full overview of all applications and their help files, start LLEAP Home.

Web Downloads

Visit www.laerdal.com/download to download the latest Direction for Use and Software.

Simulated Patient Monitor

The Simulated Patient Monitor allows any or all of the patients' status data (ECG, Heart Rate, Arterial BP waveform, NIBP, Temperature, SpO₂, CO₂, CVP etc.) to be displayed.

Blood System

The blood system includes a blood bag and three tubes with valves that can be connected to the veins in the arm and both legs. The "blood pressure" in the system can be altered by moving the blood bag up or down.

Warning: Never leave the simulator with pressurised air present in the hose as this may ruin some of the simulators functions.

Airway

Airway and Lungs/Stomach

Airway anatomy:

The infant simulator provides an anatomical realistic upper airway including:

- Oropharynx
- Nasopharynx
- Laryngopharynx
- Carina
- Trachea

The simulator simulates spontaneous respiration with the following features:

- Chest rise and fall
- The variable respiratory rate is synchronized to the Simulated Patient Monitor display and lung sounds

- The tidal volumes dynamically alter in line with the selected respiratory rate

SimBaby accepts a wide range of airway management devices and techniques. Some examples are:

- Bag-Valve-Mask (BVM)
- Oral/nasal pharyngeal airways
- Endotracheal tubes - nasal and oral [recommended ET tube size is 3.5. We also recommend the use of a malleable stylette. Care should be taken, that the stylette does not extend beyond the end of the tube, as with any direct intubation.]
- Laryngeal Mask Airways (LMA) [The Laerdal SimBaby will allow use of the LMA Classic and LMA Unique. Recommended size is 1.5.]
- Fiberoptic procedures
- Nasogastric tube insertion
- Correct form and technique are required to perform direct laryngoscopy and Endotracheal intubation
- Correct use of a variety of airway adjuncts will successfully ventilate the Patient Simulator.

The simulator contains two lungs

Too deep intubation will result in unilateral lung filling. This usually occurs on the right side, due to the accurate anatomical modeling of the tracheobronchial junction.

The airway contains a number of instructor-controlled airway complications. Using the computer user interface, the following airway functions can be activated and deactivated:

- Pharyngeal Obstruction
- Tongue Edema
- Laryngospasm
- Decreased Lung Compliance
- Increased right and/or left lung resistance
- Pneumothorax
- Stomach Decompression
- Exhale CO₂
- Variable breathing rate
- Variable breathing pattern
 - See-saw breathing
 - Subcostal retractions
 - Unilateral breathing
- Apnea
- Variable pulse oximetry display
- Breath sounds

 **Warning:** Prior to using airway adjuncts, spray all airway management devices to be inserted with a small amount of the provided airway lubricant. Use only a minimal amount of airway lubricant in the simulator's airway. Make sure to replace the Esophagus filter after every session or course if airway lubricant has been used.

Abdominal distension

Abdominal distension occurs with too high ventilation pressure while using Bag-Valve-Mask. NG Tube can be inserted. Operator needs to activate gastric distension through the software to release the air from the stomach.

Breathing

 **Warning:** To avoid damaging the spontaneous breathing bladder, do not perform chest compressions while the spontaneous breathing function is activated.

Pneumothorax

Pneumothorax on left side can be simulated through the software. The students will then see unilateral chest rise on the simulator.

Chest Drain/Chest Tube

Chest Tube insertion can be performed at the left mid-axillary site. A cut can be made at left mid-axillary line at the 4th and 5th intercostal space.

Needle Decompression

Needle decompression can be performed at the left mid-clavicular line, 2nd intercostal space. We recommend a 22-gauge needle for decompression of the chest. Using a smaller gauge needle increases the longevity of the chest skin and bladder.

Air and CO₂ Source

Compressed air is provided by a compressor or other type of pressurized air source via a regulator unit, allowing many functions to take place:

- Airway complications
- Spontaneous breathing
- Tension pneumothorax inflation
- Carotid pulse

The compressor unit operates at 110 or 230 – 240V AC. It can be connected to a CO₂ source for exhaled CO₂ to take place.

Circulation

Defibrillation

 **Warning**

- Observe all normal safety precautions for use of defibrillators.
- Connectors for external pacing are connected to the simulator's defibrillation connectors.
- Patient pads should not be used, as they do not guarantee sufficient contact.

The system has a variable pacing threshold and the ability to “ignore” pacing. Pacing capture results in a pulse synchronized with the heart rate and the display of a paced rhythm on the Simulated Patient Monitor.

The Patient Simulator is equipped with two defibrillation connectors. ECG signal can also be monitored across these connectors. Instruc-

tor can select, via appropriate keyboard command, the “Ignore Defib” function. This determines if the defibrillation shock results in conversion to a selected waiting rhythm. Manual paddle adapters (Manual Defib. Plates) are supplied for use with manual defibrillators.

 **Note:** The ECG connectors are designed for ECG monitoring only. If defibrillation is attempted over any of the ECG connectors, high voltages may be present on one or more of the uncovered connectors during the shock. (See “Cautions and Warnings” section.) Defibrillation attempts via the ECG connectors will also damage the internal electronics requiring that they have to be replaced.

Pulses

- SimBaby has palpable pulses:
 - Bilateral femoral pulse
 - Left radial and brachial pulse
- The pulses are synchronized to the simulated ECG and, when activated, the external pacemaker upon capture.
- Pulses, once activated, will remain on for approximately five seconds before reactivation is required.

 **Note:** Care should be taken when palpating pulses. Use of excessive force results in the inability to feel the pulse.

Sounds

The torso contains a number of hidden speakers that allow the realistic auscultation of sounds:

- Lung, left and right
- Heart, upper and lower

Body Movement

Body movement can be simulated and requires that the simulator lie onto a rigid base. Place a piece of cardboard or similar under the simulator if it is placed on a soft surface (bed, mattress etc.).

Drugs and IV

Right Arm is dedicated to IV skills and allows:

- Cannulation
- Phlebotomy
- Drug administration
 - Infusion
 - Legs: both legs have IO and IV access with replaceable lower legs that include the IO and IV access points.

Circulation

Blood Pressure Arm

Left Arm is a Blood Pressure Arm with radial and brachial pulses and Korotkoff sounds. Blood pressure can be measured using the supplied blood pressure cuff attached to the Link Box and installed onto the Left Arm (BP Arm) brachial position.

Blood pressure settings are controlled using the computer, and linked to the ECG functionality, so if you are changing the rhythm from a perfusing rhythm to a non-perfusing rhythm, this will also be reflected on the blood pressure settings, which will be changed according to the new type of rhythm. A non-perfusing rhythm will change respiratory rate (RR) to zero. When changing to a perfusing rhythm, the blood pressure will remain at 0/0 until changed.

Auscultation gap can also be simulated.

Connecting IV Arm

The veins are self-sealing, which allows for multiple uses; however, repetitive insertions in the same area will result in leakage sooner than if the cannulations had been spread over a wider area. Using smaller size needles will also extend the “life-time” of the veins. The venous system and the skin sleeve are both replaceable.

The IV Arm can be used with simulated blood by attaching one of the supplied blood system tubes to one of the two latex vein openings near the top of the arm. Using the supplied simulated blood concentrate, mix the desired volume of simulated blood and add to the IV bag (see *Filling the IV Bag*). Release simulated blood until it runs via the tubing into the arm and out the other latex vein. Once fluid is running freely out of the second vein, seal it using a clamp. Using a 22 gauge (or smaller) needle for IV training increases the life of the IV Arm skin. If you want the student to infuse medicines, attach a second IV bag for free flow.

 **Warning:** If a training session involves the administration of fluids and/or drugs into the IV Arm, empty the arm immediately following the training session.

Connecting IV/IO Legs

The IV/IO Legs can be used with simulated blood. Using the supplied simulated blood concentrate, mix the desired volume of simulated blood with water and add to the IV bag.

Before IV/IO practice can start, fill the lower legs with simulated blood. Use a syringe and fill through the drain plugs at the back of the lower leg (drain plug to be the highest point to avoid air remaining in the leg). Attach two of the supplied blood system tubes to the latex veins openings near the top of the lower legs and release simulated blood until it runs via the tubing into the legs. Make sure the drain plugs are not closed before this procedure is done. Close drain plug when fluid is running freely. Using a 22 gauge (or smaller) needle for IV training increases the life of the IV Leg skin and vein.

For IO training a 14 gauge (or smaller) needle is recommended. Replace the lower leg after each course if used for IO cannulations.

Filling the IV Bag

Use a syringe to inject fluid into IV bags, and fill to desired level. Control flow of blood into arm and legs, via valves.

Connecting the SpO₂ Probe

The SpO₂ probe is made up of a light diode and light sensor. When the beam between the diode and sensor is broken, the Patient Monitor Application registers that the SpO₂ probe is connected.

When the probe cable is connected to the Link Box, SpO₂ will not be displayed until the probe is placed onto the simulator. If the probe is not attached to the Link Box, SpO₂ will only be displayed automatically when selected by the user through the instructor panel on the PC.

Connecting to the Link Box

Attach simulator cable to right lower side of simulator's torso and to the connector marked "Simulator" on the back of the Link Box.

1. Connect serial cable to connector marked "PC" on the back of the Link Box and to serial port on back of your computer.
2. Connect clear tubing from blood pressure cuff to the inlet marked "BP cuff" on the back of the Link Box.
3. Connect the SpO₂ cable to the SpO₂ connector on the back of the Link Box.
4. Connect one end of the audio cable to the connector marked "Audio input" on the back of the Link Box and one to the Patient Monitor, and the other end with the mini-jack plug into the headphone outlet of your computer.
5. Plug the Link Box AC power cable into a power supply (110-240 V AC). If you are using the Portability Kit, attach to the 12 V DC input according to the Portability Kit instructions.
6. Connect external speakers, if used, to connector marked "Ext. speaker" on back of Link Box (external speakers not included).

Daily Maintenance

 **Warning:** Ensure that compressed air is not present in the simulator when training session is over.

To Clean the Simulator

Clean simulator skin with mild soap and water. Make sure that no water gets inside the simulator. Do not submerge simulator.

SimBaby Clothing

 **Note:** For washing instructions see care labels.

Modules and all other parts

Modules and all other parts should be drained, flushed through with water and air-dried thoroughly prior to storage. Disinfectants should be used where appropriate.

Esophagus Filter Assy

Replace Esophagus Filter Assy after each session (see Maintenance section for details).

Lower Legs

Replace lower legs if IV/IO skills have been practiced.

Storage

Store properly between teaching sessions. If stored in hard cases disconnect cables and tubings.

Air Filter


The air filter in the compressor should be replaced every two years. Directions for replacement are part of the Technical/Service manual.

Replacing Tension Pneumothorax Bladder

To replace bladder

Remove skin at torso sides and loosen inside tap at the umbilical area from its seat.

Lift the chest plate Assy exposing its underside.

Disconnect bladder hose from the nipple.

Remove bladder from site cavity by sliding it sideways out of the chest.


Insert new Pneumothorax Bladder into site cavity and reconnect to the nipple.
Replace chest plate Assy.
Replace chest skin over torso, fasten the tap inside the skin to its seat and secure skin at the shoulders and both sides.

Replacing Chest Tube Module

To replace pleura sleeve

Loose chest skin at left side or remove chest skin. (See *To Replace Bladder*.)

Remove chest tube insertion module from mid-axillary side of simulator.

Replace the pleura (foil) sleeve.
Put chest tube insertion module back into the torso.
Secure chest skin or replace with a new skin (see *To Replace Bladder* section).

Replacing Esophagus Filter

Note: After each session/course the esophagus filter assy shall be replaced by a new one.

To replace esophagus filter assy

Remove chest skin (see *To Replace Bladder* section).


Lift the chest plate assy exposing the internal parts and find the filter holder.

Disconnect the filter assy and replace with a new one.
Replace chest plate assy and chest skin (see *To Replace Bladder* section).

Replacing Pupils

There are three extra sets of eyes included with the simulator; normal (default with the simulator), contracted and dilated.


If you want to add some of the attached pupils to the simulator; open up the eyelids with your fingers, remove the installed pupil from its seat using the eye lens replacement device included and install the desired pupil into its seat using the same tool.


Replacing Skin and Veins in IV Arm

When excessive leaking occurs at the puncture sites, a new vein and skin should be installed to reduce loss of fluid.


To replace the vein

Pull the skin and vein off the arm. While doing this make sure you hold the simulator's shoulder (not the torso) with your other hand.

Route a new vein along the groove in the arm, leaving two ends of approximately equal lengths. Support the vein so it does not fall out of the groove while you pull on a new arm skin.

Slide a new arm skin over the arm and vein. Make sure that the two holes line up with the ends of the groove close to the shoulder. Use of talcum powder may ease this procedure.

Pull the open ends of the vein through the holes in the arm skin.

Replacing IV/IO Leg

When excessive leaking occurs at the IV/IO puncture sites, the lower leg should be replaced by a new one.

To replace lower leg

Press the bolt in the knee part out of its seat.

Replace the lower leg with a new one.

Secure the leg by pressing the bolt into its seat.

Spare Parts and Accessories

For latest version of Spare Parts and Accessories, visit www.laerdal.com

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© 2016 Laerdal Medical. All rights reserved.