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👷 Resusci Anne Wireless SkillReporter	👷 Resusci Anne Wireless SkillReporter	波 Resusci Anne Wireless SkillReporter	🗶 Resusci Anne Wireless SkillReporter				
Welcome to the Resusci Anne Wireless SkillReporter Setup Wizard	License Agreement	Select Installation Folder	Confirm Installation				
The installer will guide you through the steps required to install Resusci Anne Wireless SkillReporter on your computer.	Please take a moment to read the licence agreement now. If you accept the terms below, click '1 Agree'; then "Next". Otherwise click. "Cancet".	The installer will install Resusci Anne Wireless SkillReporter to the following folder. To install in this folder, click "Next". To install to a different folder, enter it below or click "Brows	The installer is ready to install Resusci Anne Wireless ShillRepoter on your computer. Click "Next" to start the installation.				
	LICENSE AGREEMENT FOR SOFTWARE	Eolder:					
	IMPORTANT - READ CAREFULLY: This License Agreement ("License") is a legal and non-exclusive	Disk.Cost.					
WARNING: This computer program is protected by copyright law and international treaties.	agreement between you (either an individual or a legal entity) and Laerdal Medical AS ("Laerdal") for the software product Resusci Anne Wireless	Install Resusci Anne Wreless SkilReporter for yourself, or for anyone who uses this compute					
Unauthorized duplication or distribution of this program, or any portion of it, may result in severe civil or criminal penalties, and will be prosecuted to the maximum extent possible under the law.	SkillReporter, comprised of computer software and associated media and printed materials, which may also include "online" or electronic	Eveyone					
	I Do Not Agree	O Just me					
Cancel Specie Next >	Cancel (Back Next)	Cancel (gack Next)	Cancel Gack Net>				
Resusci Anne Wireless SkillReporter	Resusci Anne Wireless SkillReporter	Nesusci Anne Wireless SkilReporter					
Installing Resusci Anne Wireless SkillReporter	Adobe Flash Player as GUI for F Wireless SkillReporter	esusci Anne Solaria Installation Complete	5				
Resusci Anne Wireless SkilReporter is being installed.	The application uses Adobe Flash Player for the graph installation dialog will appear after the installation. Plear	al user interface (GUI). A Flash Player follow the inituctions shown. If you	successfully installed.				
Please wait	already have Adobe Flash Player 11 installed on your m	chine, then exit the installation dialog. Click "Cloce" to exit.					
Please use Windows Update to check for any critical updates to the .NET Framework.							
Cancel	Cancel	Beck Next >	Cancel < Back Dose				
			Re Laerdal Medical Resusci Anne Wireless SkillReporter S				
		Adobe Flash Player 11 ActiveX	Completed				
順に「Next」をクリックし	ノ進めてください。	Please wait while Windows configures Adobe Flash Player 11 ActiveX	Show details				
必要に心し (Adobe I	Flash Player200						
アプリケーションもインス	トールされます。						
			anth				
			Cancel Nullsoft Install System V2.46 < Back Cose South States Sta				
			implement efficiently				
			The Incentration				

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