


▶ PALS Actions Overview > Legend

Action buttons (round buttons)

Clicking a round button initiates an action. Clicking this button, for example, checks the child's carotid pulse.

Menu buttons (triangular buttons)

A triangular menu button doesn't perform an action in itself, but opens a menu to the right of the button where related actions are available. The red outline around this button indicates that this menu is active.



The screenshot shows a grid of buttons categorized into: RESPONSE, AIRWAY, BREATHING, CIRCULATION, EXAMINATION & VITAL SIGNS, DRUGS & MANAGE IV/IO, TESTS & DIAGNOSTICS, and MISC. The 'CIRCULATION' category is expanded, showing buttons for: Switch providers, Check carotid pulse, Capillary refill time, CPR (active, red outline), CPR backboard, Check brachial pulse, Valsalva maneuver, Blood pressure, Length-based tape, Ice to face (active, green checkmark), Auscultate, and ECG & Defibrillator. To the right, a 'CPR' menu is open, showing options: Start CPR 15:2, Start CPR 30:2, and Start Continuous CPR. It also displays 'OXYGEN RATE 0 L/MIN' and 'VENTILATION RATE 0 /MIN'.

Active action button

The green check mark indicates that the given action is active—in this case ice has been placed on the child's face. To stop an active action simply click it again.

▶ Contents

Response	2
Airway	3
Breathing	7
Circulation	11
Examination & Vital Signs	14
Drugs & Manage IV/IO	17
Tests & Diagnostics	18
Misc	19
Other Actions	20

▶ Response

Tap and shout

Checks for response by tapping and shouting to determine the level of consciousness.

Look for breathing

Determines if the child is breathing by observing chest movement.

Sternal rub

Applies a painful stimulus by performing a sternal rub.
Not available for infant

Pinch infant/child

Applies a painful stimulus by pinching the child.

The screenshot shows a medical software interface with several tabs: RESPONSE, AIRWAY, BREATHING, CIRCULATION, EXAMINATION & VITAL SIGNS, DRUGS & MANAGE IV/IO, TESTS & DIAGNOSTICS, and MISC. The 'RESPONSE' tab is active, displaying a grid of icons for 'Tap and shout', 'Look for breathing', 'Check carotid pulse', 'Sternal rub', 'Pinch child', and 'Ask child: "Are you ok?"'. A central 'ASK PARENT QUESTIONS' section features a woman's face and icons for 'Current illness', 'Medical history', and 'Family history'. A 'QUESTIONS: CURRENT ILLNESS' menu is open, listing questions such as 'What happened?', 'When did the problems start?', 'Has there been any pain?', 'Have you noticed any change in consciousness?', 'Has he had fever?', 'Has he been coughing or wheezing or had other respiratory symptoms?', 'Has he been eating and drinking normally?', 'Has there been vomiting, diarrhea, or constipation?', and 'Has he urinated today?'.

Check carotid pulse

Performs a pulse check at the carotid artery.

Ask child: "Are you OK?"

Asks if the child is OK.
Not available for infant

Menu: Family history

Asks the parent questions in regard to family history.

Menu: Medical history

Asks the parent questions in regard to medical history.

Menu: Current illness

Asks the parent questions in regard to current illness.


Airway

Head tilt-chin lift
 Performs a head tilt-chin lift. This relieves airway obstruction caused by the tongue blocking the airway.

Suction
 Performs suction of the upper airway and airway adjuncts to remove secretions.

Jaw thrust
 Performs a jaw thrust. This relieves airway obstruction caused by the tongue blocking the airway.

Menu: Basic airways
 Presents options for basic airway adjuncts.



Length-based tape
 Uses a color-coded length-based resuscitation tape to determine correct drug dosages, equipment size, and defibrillation doses.

Menu: Auscultate
 Presents options for auscultation.

Menu: Advanced airways
 Presents options for advanced airway adjuncts.

Airway > Basic airways

Menu: Basic airways


Presents options for basic airway adjuncts.

Oropharyngeal airway


Inserts an oropharyngeal airway to improve patency of the upper airway.

Nasopharyngeal airway

Inserts a nasopharyngeal airway to improve patency of the upper airway.



Airway > Advanced airways



Menu: Airway expert

Presents options for the airway expert.

Menu: Endotracheal tube

Presents options for the endotracheal tube.

Menu: Laryngeal mask airway

Presents options for the laryngeal mask airway.

Menu: Advanced airways

Presents options for advanced airway adjuncts.

Airway > Advanced airways > Airway expert

Menu: Airway expert


Presents options for the airway expert.

Check tube

Checks the position of the endotracheal tube.

Colorimetric CO₂ detector

Detects exhaled CO₂ by colorimetry.



Withdraw the tube slightly

Adjusts the insertion depth of the endotracheal tube. Should be used to correct right main bronchus intubation.

Tube holder

Secures the placement of the tube with a commercial tube holder.


Intubation by Airway Expert

Airway expert intubates the child. Students will have to confirm and secure the device themselves.

Attach capnometer

Attaches a dynamic CO₂ detector/capnograph.

Airway > Advanced airways > Endotracheal tube



Menu: Endotracheal tube

Presents options for the endotracheal tube.

Insert tube


Inserts endotracheal tube to improve patency of the airway. The user will have to select the correct size and whether to use a cuffed or uncuffed device.

► **Airway > Advanced airways > Laryngeal mask airway**

Menu: Laryngeal mask airway
Presents options for the laryngeal mask airway.

Check tube
Checks the position of the endotracheal tube.

Colorimetric CO₂ detector
Detects exhaled CO₂ by colorimetry.




Insert mask airway
Inserts laryngeal mask airway to improve patency of the airway.

Attach capnometer
Attaches a dynamic CO₂ detector/capnograph.

Tube holder
Secures the placement of the tube with a commercial tube holder.

► **Airway > Auscultate**



Menu: Auscultate
Presents options for auscultation.

5-point auscultation
Listens to abdomen and both upper and lower lung zones to confirm proper airway device placement.

Two ventilations
Gives 2 ventilations with a self-inflating bag.

▶ Breathing

Look for breathing


Determines if the child is breathing by observing chest movement.

Attach pulse oximeter

Measures the oxygen saturation of the peripheral arterial blood and the peripheral pulse rate.

Menu: Ventilation

Presents options for ventilation.



Menu: Auscultate

Presents options for auscultation.
See page 6.

Menu: Breathing procedures

Presents options for treating pneumothorax.

Menu: Supplementary oxygen











Presents options for different oxygen supplements.

▶ Breathing > Ventilation

Menu: Ventilation
Presents options for ventilation.

Attach capnometer
Attaches a dynamic CO₂ detector/capnograph.

Attach PEEP valve
Applies PEEP with standard settings that cannot be changed.

AIRWAY	BREATHING	CIRCULATION	EXAMINATION & VITAL SIGNS	DRUGS & MANAGE IV/IO	TESTS & DIAGNOSTICS	MISC
<p>THING</p> <div style="border: 1px solid gray; padding: 5px; margin-bottom: 5px;">  <p>Ventilation</p> </div> <div style="border: 1px solid gray; padding: 5px; margin-bottom: 5px;">  <p>Supplementary oxygen</p> </div> <div style="border: 1px solid gray; padding: 5px;">  <p>Breathing procedures</p> </div>	<p>Start continuous ventilations</p> 	<p>Attach Capnometer</p> 	<p>Attach PEEP valve</p> 	<p>Warm and humidify oxygen</p> 	<p>Colorimetric CO₂ detector</p> 	<p>OXYGEN RATE 0 L/MIN</p>  <p>VENTILATION RATE 0 /MIN</p> 

Start continuous ventilations
Gives continuous ventilations with a self-inflating bag.

Warm and humidify oxygen
Reuse but change to child.

Colorimetric CO₂ detector
Detects exhaled CO₂ by colorimetry.

► Breathing > Supplementary oxygen

Menu: Supplementary oxygen

Presents options for different oxygen supplements.

Attach nasal cannula


Places a nasal cannula on the child, allowing for increased oxygen supply.

Attach CPAP

Places a CPAP mask on the child with standard settings that cannot be changed.

Attach T-piece

Places a T-piece on the endotracheal tube, allowing for spontaneous breathing and oxygen supply.



Attach nonrebreather
Places a nonbreathing mask on the child, allowing for increased oxygen supply.

Attach simple oxygen mask
Places a simple oxygen mask on the child, allowing for increased oxygen supply.

Attach capnometer
Attaches a dynamic CO₂ detector/capnograph.

Warm and humidify oxygen
Reuse but change to child.

Attach PEEP valve
Applies PEEP with standard settings that cannot be changed.


► Breathing > Breathing procedures

Needle decompression

Performs needle decompression on the right side.

Insert chest tube

Inserts a chest tube on the right side.



Menu: Breathing procedures

Presents options for treating pneumothorax.

Check chest tube

Checks if the chest tube is functioning correctly.

Circulation

Switch providers
Causes the rescuers to switch positions to prevent rescuer fatigue and maintain an adequate compression rate and depth.


CPR backboard
Places a rigid backboard to improve the quality of chest compressions.

Check carotid pulse
Performs a pulse check at the carotid artery.

Capillary refill time
Checks the capillary refill time of the child.

Menu: CPR
Presents options for CPR.

Valsalva maneuver
The child performs the Valsalva maneuver to increase vagal tone.
Not available for infant.



Length-based tape
Uses a color-coded length-based resuscitation tape to determine correct drug dosages, equipment size, and defibrillation doses.

Check brachial pulse
Performs a pulse check at the brachial artery.

Ice to face
Covers the child's forehead and eyes without obstructing the airway to elicit a vagal response that might terminate a supraventricular tachycardia.

Menu: Auscultate
Presents options for auscultation.
See page 6.


Menu: ECG & defibrillator
Presents options for various electrophysiologic interventions and cardiac monitoring options.

Menu: Blood pressure
Presents options for measuring the blood pressure of the child.

► Circulation > CPR

Menu: CPR

Presents options for CPR.



Start CPR 15:2

Performs CPR with cycles of 15 chest compressions and 2 ventilations. (Press button again to stop compressions.)


Start CPR 30:2

Performs CPR with cycles of 30 chest compressions and 2 ventilations. (Press button again to stop compressions.)

Start continuous CPR

Performs chest compressions at 100 compressions per minute continuously without pauses for ventilation. Meanwhile, the student should set the ventilation and oxygen rates delivered by the other provider (preferably 8 to 10 breaths per minute with high-flow oxygen).

► Circulation > Blood pressure



Manual blood pressure

Measures the blood pressure of the child manually.

NIBP interval

Sets the time interval for repeated automatic NIBP measurement.

Set up automatic NIBP

Initiates automatic noninvasive blood pressure (NIBP) measurement.


Remeasure NIBP

Triggers remeasurement of the NIBP.

Menu: Blood pressure

Presents options for measuring the blood pressure of the child.

► **Circulation > ECG & defibrillator**



Attach pads
Applies the AED pads.

Attach 3-lead ECG
Attaches a 3-lead ECG monitor.

Stand clear
Gives a verbal warning to stay clear of the child. Any rescuers in contact with the victim will stop touching the victim.

Obtain 12-lead ECG
Presents recorded 12-lead ECGs of the child.

Check electrodes
Checks the placement of the ECG electrodes.

Menu: ECG & defibrillator
Presents options for various electrophysiologic interventions and cardiac monitoring options.

Change ECG display lead
Changes the ECG lead displayed on the monitor.

Examination & Vital Signs

Check skin

Checks the condition of the child by inspecting clinical signs of the skin.

Examine ears


Examines the ears with an otoscopic procedure. Might reveal relevant information about possible etiology of an infection.

Capillary refill time

Checks the capillary refill time of the child.

Menu: Vital signs

Presents options for checking the vital signs in one place.



Length-based tape

Uses a color-coded length-based resuscitation tape to determine correct drug dosages, equipment size, and defibrillation doses.

Menu: Head-to-toe examination

Allows for full-body clinical examination, including findings from inspection and palpation of the different regions of the body.

Menu: Examine eyes

Presents options for examining the eyes of the child.

Menu: Auscultate

Presents options for auscultation.

See page 6.

Examination & Vital Signs > Vital signs

Look for breathing


Determines if the child is breathing by observing chest movement.

Attach pulse oximeter

Measures the oxygen saturation of the peripheral arterial blood and the peripheral pulse rate.

Menu: Blood pressure

Presents options for measuring the blood pressure of the child.



Measure temperature

Measures the temperature of the child.

Measure glucose

Checks the capillary glucose level.


Menu: Vital signs

Presents options for checking the vital signs in one place.

Check carotid pulse

Performs a pulse check at the carotid artery.

Examination & Vital Signs > Vital signs > Blood pressure



Menu: Blood pressure

Presents options for measuring the blood pressure of the child.

NIBP interval

Sets the time interval for repeated automatic NIBP measurement.

Manual blood pressure

Measures the blood pressure of the child manually.

Remeasure NIBP

Triggers remeasurement of the NIBP.

Set up automatic NIBP

Initiates automatic noninvasive blood pressure (NIBP) measurement.


▶ Examination & Vital Signs > Examine eyes

Test eye movement

Checks the eye movements of the child.

Test pupil responsiveness

Checks the pupillary light reflex of the child.




Menu: Examine eyes

Presents options for examining the eyes of the child.

Eyelash and corneal reflex

Checks the ciliary and/or corneal reflex of the child.

▶ Examination & Vital Signs > Head-to-toe examination



Menu: Head-to-toe examination

Allows for full-body clinical examination, including findings from inspection and palpation of the different regions of the body.

▶ Drugs & Manage IV/IO

Menu: Administer drugs and fluids

Presents the various drugs and fluids available for treating the child.

Menu: Manage IV/IO

Presents options for placing intravenous (IV) and/or intraosseous (IO) cannulas at different sites.

The screenshot shows a software interface for medical administration. At the top, there are tabs for 'RESPONSE', 'AIRWAY', 'BREATHING', 'CIRCULATION', 'EXAMINATION & VITAL SIGNS', 'DRUGS & MANAGE IV/IO', 'TESTS & DIAGNOSTICS', and 'MISC'. The 'DRUGS & MANAGE IV/IO' tab is active. On the left, there are three icons: a bottle and syringe labeled 'Administer drugs and fluids', a syringe and needle labeled 'Manage IV/IO', and a color-coded tape labeled 'Length-based tape'. The main area is divided into three sections: 'SELECT DRUG' with a list of options under 'FLUIDS & ELECTROLYTES' (Calcium chloride is selected), 'SELECT ROUTE' with 'IV/IO push' selected, and 'SELECT DOSE' for 'CALCIUM CHLORIDE' showing '000.0 mg' and a 'Slow injection' checkbox. An 'Administer' button is at the bottom right.

Length-based tape

Uses a color-coded length-based resuscitation tape to determine correct drug dosages, equipment size, and defibrillation doses.

Tests & Diagnostics

View lab results

Allows the student to study the results of the blood tests ordered.

Arterial blood gases

Performs an arterial blood gas test. The results are returned in the medical record after a short while.

Not available for infant

Venous blood sample


Collects venous blood samples. The results are returned in the medical record after a short while.

Urine dipstick and culture

Performs a urine dipstick test, reports the immediate findings, and sends the sample for culturing.

X-ray

Performs on-site chest x-ray.



Capillary blood gas

Performs a capillary blood gas test. The results are returned in the medical record after a short while.

Measure glucose

Checks the capillary glucose level.

Lumbar puncture

Performs a lumbar puncture, reports the immediate findings, and sends the sample for culturing.

Obtain 12-lead ECG

Presents recorded 12-lead ECGs of the child.

Venous blood gas

Performs a venous blood gas test. The results are returned in the medical record after a short while.

Blood culture

Obtains a blood culture. Results are not returned during the simulation.

View diagnostics results

Presents results of the diagnostics performed during simulation.

Misc

RESPONSE	AIRWAY	BREATHING	CIRCULATION	EXAMINATION & VITAL SIGNS	DRUGS & MANAGE IV/IO	TESTS & DIAGNOSTICS	MISC
<p>POSITION PATIENT</p> <p>Supine</p> <p>Sit up</p>		<p>EXPOSURE</p> <p>Blanket</p> <p>Ice to face</p> <p>Therapeutic hypothermia</p>	<p>Measure temperature</p>	<p>VARIOUS CATHETERS</p> <p>Insert nasogastric tube</p> <p>Insert urinary catheter</p> <p>Check catheter output</p>		<p>HELP</p> <p>Call MET</p> <p>Consult Cardiology</p> <p>Consult Neurosurgery</p>	<p>Open handbook</p>

Supine
Places the child in a supine position.

Blanket
Places a blanket on the child.

Ice to face
Covers the child's forehead and eyes without obstructing the airway to elicit a vagal response that might terminate a supraventricular tachycardia.

Measure temperature
Measures the temperature of the child.

Insert nasogastric tube
Inserts a nasogastric tube.

Open Handbook
Pauses the simulation and shows the *2010 Handbook of Emergency Cardiovascular Care for Healthcare Providers*.

Sit up
Places the child in a sitting position.

Therapeutic hypothermia
Begins preparation for post-cardiac arrest therapeutic hypothermia.

Insert urinary catheter
Places a urinary catheter in the bladder of the child.

Check catheter output
Checks the urinary output from the child.

Call MET
Activates the medical emergency team. This gives the learner the option to emphasize the need to call for help; however, there will be no response to this action during the simulation.

Consult Cardiology
Calls for cardiologic consultation. This gives the learner the option to emphasize the need to consult relevant expertise; however, there will be no response to this action during the simulation.

Consult Neurology
Calls for neurologic consultation. This gives the learner the option to emphasize the need to consult relevant expertise; however, there will be no response to this action during the simulation.

Other Actions

Call code

Activates the emergency response system, resulting in emergence of 2 additional rescuers to assist with CPR later.


Medical record

Opens the medical record to review the child's history, clinical signs, and lab findings obtained during the case.



Transfer patient

Ends the scenario by transferring the child. Presents options for ending the simulation.



Discharge patient

Ends the scenario by discharging the child from the hospital.

Transfer patient

Ends the scenario by transferring the child. Presents options for ending the simulation.

Stop resuscitation efforts

Ends the scenario by aborting all interventions and withholding further treatment.