EN

FR

DE

ES

IT

BR

NL

PL


JA

ZH

KO

# SimJunior

Directions for Use

<b>Introduction</b>	<b>2</b>
– Regulatory Information	2
– General Simulator Care	3
– SimJunior Overview Picture	5
– SimJunior Features	6
<b>Simulation Software</b>	<b>8</b>
– Instructor PC	8
– Software License	8
– Laerdal Simulation Home	8
– LLEAP	8
– Voice Conference Application	8
– Patient Monitor	8
– Session Viewer and SimView Server	8
– Other Applications	8
– SimDesigner	8
– Help Files	8
– Web Downloads	8
<b>Simulator Setup</b>	<b>9</b>
– Preparing for Simulations	9
– Preparing the IO Leg	9
– Administering IV Fluids	9
– Draining Excess IV Fluid	10
– Cleaning the IV Arm	10
– Installing Defibrillation Adapter Plates	10
– Connecting the Blood Pressure Cuff	10
– Calibrating the Blood Pressure Cuff	10
– Calibrating the Blood Pressure Cuff in SimPad	10
– Changing the Pupils	11
– Using an External Compressor	11
<b>Maintenance</b>	<b>12</b>
– Maintenance after each Simulation Session	12
– Opening the Torso	13
– Replacing the Simulator Battery	13
– Charging the Simulator Battery	14
– Replacing the Chest Rise Bladder	14
– Replacing the Lung Bladder	14
– Replacing Simulator Skins	15
– Replacing the Right IV Arm	15
– Replacing the Right IO Leg	15
<b>Troubleshooting</b>	<b>16</b>
<b>Specifications</b>	<b>19</b>
<b>Spare Parts and Accessories</b>	<b>21</b>

## Introducing SimJunior™ Simulator

SimJunior is a realistic, life-size pediatric patient simulator designed by Laerdal in collaboration with the American Academy of Pediatrics. SimJunior can be used for the education of individuals or teams of healthcare professionals from multiple medical disciplines.

The SimJunior Simulator comes completely assembled. There is no need to assemble parts of the Simulator.

For details on preparing SimJunior for simulation, go to “Simulator Setup” section.

## SimJunior Help

### Directions for Use (DFU)

Step-by-step instructions and illustrations for using the SimJunior Simulator.

### Original Manufacturer’s User Manuals

All separate user manuals and labeling from original manufacturers should be followed. The SimJunior Directions for Use does not replace or supersede those from the original manufacturer.

### SimJunior Quick Setup Guide

Step-by-step guide for setup of the SimJunior simulation system.

### SimPad Directions for Use (DFU)

SimPad Directions for Use provides instructions for using the Link Box and the SimPad with SimJunior.

### LLEAP Help Files (Laerdal Learning Application)

The LLEAP Help files are accessible from the Laerdal Simulation Home, see *Laerdal Simulation Home*.

Help topics include:

- LLEAP
- Session Viewer
- Voice Conference Application
- SimDesigner
- Trend Editor

### Technical Assistance

For technical assistance, contact your local Laerdal Technical Service Center.

### Web Downloads

Visit [www.laerdal.com/downloads](http://www.laerdal.com/downloads) to download the latest Directions for Use and newer versions of the simulation software.

## SimCenter Content

SimCenter offers easy access to validated content from worldwide simulation experts so you can take advantage of their experience and get the most from every learning opportunity. You will find a wide range of scenarios that have been developed specifically for the SimJunior. Visit [www.mysimcenter.com](http://www.mysimcenter.com) to learn how you can fully optimize your learning experience.

## Regulatory Information

### Federal Communications Commission Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

### Caution

Changes or modifications not expressly approved by the party responsible for compliance could void the user’s authority to operate the equipment.

**CE** The product is in compliance with the essential requirements of Council Directive 2004/108/EC on electromagnetic compatibility (EMC). The product is in compliance with Council Directive 2011/65/EU on restriction of the use of certain hazardous substances (RoHS).

**UK CA** UKCA (U.K. Conformity Assessed): This product complies with the requirements of U.K. legislation of Electromagnetic Compatibility Regulations (S.I. 2016/1101). This product complies with the requirements of U.K. legislation of Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations (RoHS) Regulations (S.I. 2012/3032).

 Dispose of in accordance with local requirements and regulations.

 Li-ion batteries should be recycled.

## Disclaimer

Use of the SimJunior patient simulation system to train personnel should be undertaken under supervision of suitably trained technical or medical personnel with an understanding of educational principles as well as recognized medical protocols. As with all Simulators or other such training devices, there may be approximations, variations and inaccuracies in anatomical features and the physiological modeling. This being the case, Laerdal does not guarantee that all features are completely accurate.

## Global Warranty

See the Laerdal Global Warranty Booklet, or see [www.laerdal.com](http://www.laerdal.com).

Country of origin - SimJunior is made in USA.

Laerdal Medical  
P.O. Box 38  
226 FM 116  
Gatesville, Texas 76528, USA

## General Simulator Care

Follow the instructions below to maintain optimum performance and longevity of the Simulator and its components.

### General Care

To maintain Simulator skins, wash hands before use and place the Simulator on a clean surface. Laerdal recommends to:

- Use gloves during simulation scenarios. Avoid using colored plastic gloves, as they may cause discoloration of the Simulator skin.
- Do not use felt-tipped markers, ink pens, acetone, iodine, or other staining medications near the Simulator. Take care not to place the Simulator on newsprint or colored paper. Staining may be permanent.
- Clean Simulator skins with mild soap and water.

Do not attempt to perform the following techniques on this Simulator due to the inability to properly sanitize the airway:

- Mouth-to-mouth/mouth-to-mask ventilation
- Insertion of simulated vomit or fluids for suctioning
- If a training session involves the use of fluids in the IV arm, drain the arm immediately following the training session.
- Only use Laerdal Airway Lubricant, and apply sparingly
- Rinse, clean, and dry Simulator component modules
- Fold the torso skin back and powder the inside of the torso skin to decrease friction. Do not spill powder into Simulator chest cavity.

### Environment

In cold conditions, wait until the Simulator has reached room temperature before starting up.

To avoid overheating and reduce wear:

- When using in temperatures above 40°C (104°F), always allow the Simulator to cool down between training sessions.
- When using in a bed, Simulator should not be covered with heavy bedding which prevents heat transfer from the Simulator.

## General Simulator Handling

Take the following precautions to avoid personal injury or damage to the product:

- Introduce fluids into the Simulator only as directed in this document
- Lubricate the oral and nasal airways with the lubricant provided prior to inserting any instrument, tube, or airway device. Also lubricate instruments and tubes prior to use.
- Do not introduce humidified air into the system during ventilation
- Do not use the Simulator if the internal tubing and cabling is disconnected
- Never use the SimJunior Simulator outdoors in wet conditions, as this may pose a shock hazard or damage the Simulator.
- Using a defibrillator in temperatures over 35° C (95° F) may cause overheating

 **Warning:** Do not use automated chest compression machines on the Simulator.

 **Caution Latex:** This product contains Natural Rubber latex, which may cause allergic reactions when in contact with humans.

## Defibrillation Hazards

SimJunior allows for defibrillation.

A conventional defibrillator may be used on SimJunior. During defibrillation, the defibrillator and Simulator may present a shock hazard. All standard safety precautions must be taken when using a defibrillator on the Simulator:

- Read and follow all safety and operating instructions provided with your defibrillator and associated equipment.
- Follow defibrillation protocol by avoiding contact between the external paddles and any of the electrode sites while defibrillating. A load of up to 300 joules can be delivered.
- Failure to follow safety measures could result in injury or death.

For more information, consult your defibrillator's original user manual.

 **Warnings:**

- When performing defibrillation, use the defibrillator connectors or the zap plate mounted on the Simulator's chest. Do not use the ECG connectors which are designed exclusively for ECG monitoring. Defibrillation on the ECG connectors will damage the internal electronics of the Simulator and may cause personal injury.
- Do not defibrillate the Simulator when it is turned OFF or if it is not functioning normally.

To prevent overheating during defibrillation, do not exceed a defibrillation sequence of 3 shocks in 45 seconds followed by 1 minute of CPR. After 30 minutes, cease all shocking for at least 15 minutes before starting a new sequence. Do not repeat this for more than a 4 hour period. Also, do not provide more than 2 x 300 defibrillator discharges per minute.

## Warnings:

- The Simulator must not come into contact with electrically conductive surfaces or objects during defibrillation.
- Avoid use in all flammable environments. For example, high levels of pure oxygen should be avoided during defibrillation. Ensure good ventilation if concentrated oxygen is used near the Simulator.
- The Simulator torso must always be kept dry.
- Allow the Simulator to acclimate before defibrillating. Sudden changes in temperature may result in condensation collecting on electronic components, which could pose a shock hazard.
- To prevent torso skin electrode pitting, do not apply conductive gel or conductive defibrillation pads intended for patient use.
- Pressing down too hard on the defibrillation connectors during defibrillation may also cause arcing and pitting.
- Do not defibrillate the Simulator if the torso skin is not in place.

## Mechanical or Electrical Hazards

Do not use the SimJunior Simulator if:

- Limbs are not attached to the torso
- Skins are torn or not properly fastened
- Internal or external cables, tubes or connectors are damaged
- There is fluid leakage in or on the Simulator
- There are unusual sounds indicating air leakage or mechanical damage
- There are signs of electrical malfunction, such as an unresponsive Simulator or unusual smell or smoke

## Warning:

- Avoid pinch hazards - Do not use the Simulator without the external skins

## Battery Use and Maintenance

- Always use batteries approved to power the SimJunior; Link Box, and SimPad.
- Ensure that the batteries are properly installed. Inserting and connecting batteries incorrectly could cause a short circuit.

## Warnings:

- Dispose of batteries according to local regulations.
- The external battery charger is for indoor use only.
- The Simulator battery should only be charged in temperatures ranging from 0 °C - 40°C (32 °F - 104 °F)
- Do not mistreat, disassemble, or attempt to repair the battery. Do not use the batteries if they are visibly damaged, malfunctioning, or appear to be leaking.

- Take extreme care to avoid direct contact with electrical, hot or smoking parts. In case of a leaking battery, disconnect and remove the battery when it is judged safe to do so.
- Exposure to fluids pose an explosion hazard.
- On every 30th charge cycle, drain the battery completely before recharging. To drain the battery, run the Simulator on battery power until automatic shutdown.
- Only replace with a Laerdal SimJunior battery.


## Storage and Transportation

- Never store fully charged batteries for longer than a month.
- The Simulator battery can be transported in the Simulator during air freight.
- When transporting spare batteries, contact the airline or freight company for the latest transport regulations.

## Warning:

- The SimJunior Simulator and accessories are heavy when combined in the carrying case. Always ensure that SimJunior is firmly secured during transportation and storage to prevent personal injury or damage to the product.

# SimJunior Simulator

## SimJunior Simulator

SimJunior facilitates interactive training of life-saving skills and responds to clinical intervention by instructor control and/or pre-programmed scenarios for effective practice of diagnosis and treatment of a patient.

With spontaneous breathing, airway control, voice, sounds, ECG, and other clinical features, SimJunior is a fully functional pediatric simulator. SimJunior allows observation and recognition of most vital signs, which enables the instructor to assess the student's skills based on a realistic clinical situation.

## SimJunior Features

- Realistic airway for simulation of difficult airway management, oral, and nasal intubation
- Observable breathing
- Cardiac features, including defibrillation and cardioversion
- Eyes with interchangeable pupils (normal, dilated, or constricted)
- Convulsions to simulate seizures
- Chest compressions
- Vascular access
- Normal and abnormal heart, breath, and bowel sounds
- Automatic Simulation Control based on pre-programmed and validated patient scenarios

Simulations can run autonomously using scenarios on the SimPad or with LLEAP (PC). The development of the patient's condition is pre-programmed and automatically responds according to the participant interventions.

## Overall Dimensions

Length / Width (Simulator only): 45 in x 9.84 in (120 cm x 42.5 cm)  
Weight (Simulator only): 25 lbs (11.36 kg)

## SimJunior Clothing

*SimJunior* comes with custom designed clothing with Velcro openings for easy removal. Washing instructions are found on the clothing label.

- Shirt
- Shorts
- Boxer Shorts

## General Clinical Features

### Airway Features

The airway is anatomically correct to the trachea. The airway can be manipulated by:

- Head tilt/chin lift
- Jaw thrust with articulated jaw
- Cricoid pressure and manipulation
- Suctioning (oral and nasopharyngeal)

The Simulator may be ventilated by normal and emergency methods:

- Bag-mask ventilation
- Orotracheal intubation
- Nasotracheal intubation

Prior to using airway devices, lubricate with Laerdal Airway Lubricant.

The following equipment or methods are suitable to secure the Simulator's airway:

- Laryngeal mask airways (size #2.5)
- Endotracheal tube intubation (size ID 4.5 cuffed, 5.5 uncuffed)

Use of a malleable stylet is recommended. Make sure it does not extend beyond the ET tube.

The following Simulator conditions indicate incorrect tube placement:

- Right main stem intubation – unilateral chest rise
- Stomach distention
- Lack of chest sounds (see Breathing section)

Simulator features may be configured to present various airway scenarios:

- Tongue edema - normal, medium, maximum levels
- Lungs open or closed

## Breathing Features

SimJunior can simulate spontaneous breathing with visible chest rise and fall and variable breathing rates. The breathing is generated by an enclosed air compressor in the Simulator's right thigh.

- Bilateral chest rise and fall with spontaneous breathing
- Unilateral chest rise and fall with right mainstem intubation during ventilations
- Unilateral and bilateral lung sounds
- Normal and abnormal breath sounds
- Variable respiration rate (0-60 breaths per minute)
- Anterior auscultation sites (4)

The left lung and right lung can be closed independently or together to create a partial or complete airway obstruction. The SimJunior Simulator can also be used with assisted ventilations.

**Note:** Lungs are not intended for use with PEEP-valves.

## Circulatory Features


### Cardiac:

- Extensive ECG library, pulses from 0-200
- Heart sounds - anterior location (1)
- ECG rhythm monitoring (4-connector, 3-lead ECG)
- 12-lead ECG display (Using LLEAP only)
- Pacing
- Defibrillation and cardioversion using live defibrillators


### Defibrillation:

- The number of shocks required for automatic conversion are set in each simulation Patient Case

### Defibrillation Studs:

### 3-Lead ECG Studs:

### Blood Pressure and Pulses:

- BP measured manually by auscultation of Korotkoff sounds
- Bilateral carotid pulse
- Central pulses can be set to Normal, Weak, Absent
- Radial/Brachial pulses can be set to Normal, Weak, Absent
- Carotid, brachial and radial pulses synchronized with ECG
- Pulse strength variable with BP
- Pulse palpation is detected and logged

## CPR Features

- Compliant with 2010 Guidelines
- Compressions generate palpable pulses, blood pressure wave form, and ECG artifacts
- Realistic compression depth and resistance
- Detection of depth and rate of compressions

## Vascular Features

### IV Access:

IV access is possible on the right arm and hand. The IV arm can be set up for IV insertion, infusion, and bolus into:

- Peripheral veins of the forearm
- Antecubital fossa and dorsum of the hand

For use of an optional catheter, see your local Laerdal representative.

### Intraosseous Access (IO):

Access for IO infusion is possible:

- Tibial Tuberosity (right)
- Medial Malleolus

**Note:** SimJunior does not have intra-muscular (IM) injection sites. Do not attempt to perform IM injections.

## Sound Settings for SimJunior

SimJunior supports multiple sites for auscultation and recognition of heart, lung, and bowel sounds. SimJunior also allows you to set vocal sounds.

For more information on setting sounds for SimJunior, see SimPad Directions for Use or LLEAP Help Files.

## Instructor PC

The Instructor PC uses LLEAP software to manage simulation sessions. LLEAP, as well as other useful simulation related applications, is accessible from Laerdal Simulation Home.

## Software License

The Instructor PC is shipped from Laerdal Medical with a pre-activated license for LLEAP installed.

Changes or updates of the computer hardware (e.g. new hard drive or mother board) may render the license invalid. Please contact your local Laerdal support for assistance for re-activating the license.

## Laerdal Simulation Home

Laerdal Simulation Home is an application from where LLEAP and other Laerdal programs related to patient simulation can be found and started. The help files can also be opened from here. Laerdal Simulation Home is located in the Laerdal Medical folder under the Windows start menu (Windows 7).

Software used in a simulation session can be divided in the following main applications:

- LLEAP (Laerdal Learning Application)
- Voice Conference Application
- Patient Monitor
- SimView Server or Session Viewer

In addition SimDesigner and other applications are used for designing or preparing a simulation.

## LLEAP

The Laerdal Learning Application (LLEAP) is the instructor's application from where the simulation session is run, controlled, and monitored. LLEAP can be operated in Automatic or Manual mode. Automatic mode is used for pre-programmed scenarios while Manual mode allows the instructor full manual control over the simulation session. Running simulations in Manual Mode generally requires some medical expertise to create clinically sound simulations.

## Voice Conference Application (VCA)

The VCA software is needed to transmit all vocal sounds used during simulation. It allows the instructor to communicate through the simulator during the session. VCA can also be used to communicate with other instructors on a network, and create separate channels that only members can talk and listen to.

## Patient Monitor

The Patient Monitor application emulates a typical hospital patient monitor. It is the learner's console and can be set up and controlled by the instructor, as well as by the learner, through on-screen touch menus.

## Session Viewer and SimView Server

Session Viewer and SimView Server are applications that record video and patient monitor screen captures during simulation, in addition to providing an interface to debrief your session. After a session is ended, log files generated in LLEAP are transferred and merged with the video files in Session Viewer or SimView Server for the debrief.

Session Viewer typically runs locally on the same computer used for LLEAP, SimView Server runs on a dedicated server in the local network. During the first start-up of LLEAP you are prompted to select a debriefing system available on your computer or local network. This can be changed later.

## Other Applications

There are also other programs that are used in conjunction with the simulation sessions, for example License Manager for handling program licenses and Simulator Firmware & Network Fixer for updating the firmware of the simulators or troubleshooting network problems.

## SimDesigner

SimDesigner is an application for configuring your own pre-programmed scenarios. It can also be used to analyze and print out a graphical representation of a scenario.

SimDesigner must be installed to allow conversion of legacy instructor application files to LLEAP compatible file formats.

## Help Files

For more information on the LLEAP software and other applications, see the Help files:

1. Click Windows <Start> <All programs>.
2. Select <Laerdal Simulation Home>.
3. Select <Help>.

## Web Downloads

Visit [www.laerdal.com/downloads](http://www.laerdal.com/downloads) to download the latest User Guide and Software.

## Preparing for Simulations

After the SimJunior Simulator is prepared and all applicable interactive devices (SimPad or Instructor PC and Patient Monitor) are configured and made ready for use, final preparations may need to be performed before SimJunior is ready to run simulations.

Final simulation preparations include the following activities:

- Preparing the IO leg
- Administering IV Fluids
- Draining Excess IV Fluid
- Cleaning the IV Arm
- Connecting Defibrillation Adapter Plates
- Connecting the Blood Pressure Cuff
- Changing the Pupils
- Connecting an External Compressor
- Changing Settings for the Internal Compressor

## Preparing the IO Leg

The Simulator's right leg is designed for practicing intraosseous infusion skills. Anatomy of the right leg includes knee, tibia, tibial tuberosity site, and medial malleolus site. The right leg comes with IV tubing and one IV bag for use as a reservoir during infusion.

Additional equipment needed:

- IV bag and administration set
- Simulated blood concentrate
- 35cc syringe and 16-gauge intraosseous needle

To fill the IO Leg and prepare it for use:

1. Raise the leg vertically.

2. Open the plug at the heel.
3. With a syringe, fill the bone with approximately 240 to 250 cc of simulated blood.
4. Attach tubing connected to an IV reservoir bag. (Use roller clamp to open or close the line as needed in order to relieve pressure build-up.)
5. Close off tubing that runs into the reservoir bag so that simulated blood will remain in the bone for aspiration.
6. When the reservoir bag is filled, discontinue infusion and replace with empty bag.

Before storing the Simulator:

- Detach the IV line and reservoir from the foot
- Allow all fluid to drain
- Store leg unplugged to allow it to air dry

Note: Do not use saline solution. It may cause deterioration of the leg. Use distilled or de-ionized water to avoid buildup of minerals.

## Administering IV Fluids

The right IV Arm supports training for IV drug administration, IV insertion, infusion, and bolus into peripheral veins of forearm, antecubital fossa, and dorsum of the hand.

Use 20 to 22 gauge needle for IV simulation. To prevent clogging of the IV system, use only distilled or de-ionized water to simulate IV drugs.

## Draining Excess IV Fluid

The IV fluid system is an open system. IV fluids are drained as they are administered.

Before each session:


1. Attach the IV overflow tube (in the Simulator's right arm) to an overflow container (such as an IV bag).
2. Allow excess fluid to drain into a container during the simulation.

## Cleaning the IV Arm

Clean the IV arm after each session or day of use by flushing the IV Arm with 60% isopropanol or 70% ethanol.

## Installing Defibrillation Adapter Plates

The Simulator torso is fitted with two stud connectors for defibrillator cables. Before attempting defibrillation, you must unscrew and remove the defibrillator stud covers.

Two defibrillation adapter plates are packaged with the Simulator. The defibrillator adapter plates must be fitted in place before using a live defibrillator with defibrillation paddles or adhesive pads. Screw the adapter plates firmly into place.

## During Defibrillation

A conventional defibrillator (not included) may be used on the SimJunior. During live defibrillation, the defibrillator and Simulator may present a shock hazard. All standard safety precautions must be taken when using the defibrillator on the Simulator.

**Caution:** The Simulator must not be in contact with electrically conductive surfaces or objects during defibrillation.

## Warnings:

- Defibrillation must be performed on the defibrillator connectors only.
- Do not press too hard over the defibrillation adapters as this may cause arcing and pitting.
- Do not defibrillate the Simulator without the torso skin in place.

## Connecting the Blood Pressure Cuff

The Simulator is delivered with a customized Blood Pressure Cuff. Connect the tube to the white BP connector on back of the Link Box before use.

**Note:** The speaker for the Simulator's blood pressure is located in the left antecubital fossa, circled below. The radial and brachial pulses are located in the rectangular areas.

## Calibrating the Blood Pressure Cuff


1. Select <Tools> from the menu
2. Select <Simulator Setup> from the Tools-menu
3. Select <Calibrate Blood Pressure> from the sub-menu
4. Follow the on screen wizard instructions to perform the calibration

## Changing the Pupils

SimJunior is delivered with normal pupils installed in the eyes. A separate kit containing plastic pupil inserts (constricted and dilated) comes with SimJunior:

To change the pupils:

1. Carefully open the eyelids wide, taking care not to tear the Simulator face skin.
2. Using the suction cup tool provided in the kit or with the edge of your fingernail, carefully remove the pupil from the eye.
3. Replace the pupil with the desired insert, using the suction cup tool or by gently pressing in place.

## Using an External Compressor

The external compressor connects to the clear tube found in the cable bundle exiting the right side of the Simulator. When using an external compressor you can turn off the internal compressor from SimPad or LLEAP.

For more information on external compressors and regulator panels compatible with SimJunior, contact your local Laerdal representative.

### Turning on or off the internal compressor in SimPad

1. Start **<Manual Mode>**, select a theme and start the session.
2. Select **<Miscellaneous>** from the bottom menu.
3. Select **<Compressor>** from the sub-menu.
4. Make the desired changes in the dialog and click OK.

### Turning on or off the internal compressor in LLEAP

1. Select **<Tools>** from the menu.
2. Select **<Simulator Setup>** from the Tools-menu.
3. Select **<Internal compressor>** from the sub-menu.
4. Make the desired changes in the dialog and click OK.

### Maintenance after each Simulation Session

The following preventive measures are required to ensure longevity of the SimJunior Simulator:

- Power off the SimPad, Link Box and PCs
- Charge batteries if necessary
- Flush the IV-arm system and drain the IO leg (leave the plug in the bottom of the IO leg open)
- After each session where the IV-arm has been used, flush the IV arm with 60% Isopropanol alcohol or 70% ethanol
- Wipe the skin with a moist cloth
- Remove wet clothes or linens
- Remove tape residue with a moist cloth or rubbing alcohol
- Based on the use of the Simulator, replace modules that are consumed or damaged
- Replace cables or connectors showing visible damage
- Servicing should be performed by qualified personnel at regular intervals

Always perform a service:

- If liquids have been spilled in or on the Simulator
- After use in dusty environments

**Note:** For information on and spare parts, see section: Spare Parts and Accessories.

### Before Storage or Shipping

- Flush the IV arm with 60% Isopropanol alcohol or 70% ethanol. Open the drain valve on the bottom of the Simulator's right foot.

## Opening the Torso

Open the Simulator torso to perform procedures such as:

- Attaching or replacing limbs
- Replacing the Simulator Battery
- Replacing the Chest Rise Bladder and Lung Bladder
- Replacing the Torso Skin
- For general inspection

**Warning:** Open the torso from left to right, gently unhook the skin. Do not pull hard to lift the Torso Skin, this could damage the chest cables.

1. Unhook Torso Skin holders on left side of the torso and each shoulder:

2. Fold the Torso Skin over to the right side.
3. Disconnect the Defibrillation connector cable.

4. The Defibrillation connector cable must be disconnected at the right side of the Simulator below the chest form.

**Note:** DO NOT disconnect the tubes and cables connecting the stomach foam to the Simulator.


5. Replace the Torso Skin by performing steps 1- 4 in reverse.

## Replacing the Simulator Battery


1. Open the Simulator torso.

2. Slide the battery out of the battery strap.

3. Pull the battery power connection up from inside the Simulator. Unhook the battery power connection.

4. Lift the battery out of the Simulator.
5. Replace the charged battery or install a new Laerdal approved SimJunior battery.
6. Connect the battery to the Simulator.
7. Replace the Torso Skin.

## Charging the Simulator Battery

The Simulator battery charges itself as it is operated. If necessary, the Simulator battery can be charged using an optional external power adapter:

1. Power down the Simulator via the Link Box.
2. Remove the Simulator battery from the Simulator.
3. Connect an approved external battery charger (with a power cord that meets local specifications) to the bottom of the Simulator battery. Connect the battery charger's power cord into a wall outlet.



4. When the battery has finished charging, replace the battery in the Simulator.

## Replacing the Chest Rise Bladder

The Chest Rise Bladder is located in the lower end of the chest plate, centered under the lungs.

1. Open the Torso Skin to expose the chest plate. See *Opening the Torso*.
2. Disconnect the clear tubing from the Lung Bladder and remove the Lung Bladder. For more details, see *Replacing the Lung Bladder*.

3. Disconnect the clear tubing from the Chest Rise Bladder.
4. Discard the old bladder.
5. Insert new bladder.
6. Re-connect the tubing to the new bladder.
7. Close the torso and replace the Torso Skin.

## Replacing the Lung Bladder

If leaking occurs, the Lung Bladder (in the chest cavity) should be replaced.

1. Open the Torso Skin.
2. Disconnect the clear tubing from under each side of the lung bladder. Refer to step 2 illustration in *Replacing the Chest Rise Bladder*.
3. Disconnect the Lung Bladder clip from the top center.

4. Gently pull the top center portion of the Lung Bladder up from its socket.
5. Reverse the process to install the new Lung Bladder.

## Replacing Simulator Skins

The Simulator skins may need to be replaced if they are torn, perforated, or stained.

### Torso Skin:

To replace the Simulator torso skin, see *Opening the Torso*.


### Arm Skin:

To replace the Simulator Arm Skin:

1. Remove the old Simulator Arm Skin by working it down over the arm and hand or \*cut the skin.


**Warning:** \*Do not puncture or cut the mandrel arm.

2. Discard old skin.
3. Lubricate inside of new Arm Skin with mild liquid detergent mixed with water.
4. Lubricate (or lather) mandrel with a mix of mild liquid soap and water.
5. Slide Simulator hand into skin.
6. Work skin over fingers (as with a glove).
7. Work skin up, over the arm mandrel.


## Replacing the Right IV Arm

1. Open the torso, see *Opening the Torso*.
2. Remove old arm by unscrewing the bolt holding the arm to the Simulator's shoulder:

3. Ensure that the shoulder screw is loose enough to allow the arm axle to slide out easily.
4. Discard the old arm according to specified guidelines.

## To attach the new arm:

1. Ensure that the shoulder screw is loose enough to allow the arm axle to slide easily into place.
2. Connect the arm to the corresponding connection points in the torso.
3. Carefully push the arm axle into the shoulder bracket, so that the axle is flush with the inside of the bracket.
4. Tighten the shoulder screw with an Allen wrench.

## Replacing the Right IO Leg

To detach the right lower leg:

1. Ensure that any fluid in the leg has been drained.
2. Remove old IO Leg by unscrewing the bolt holding it to the thigh.
3. Slide the IO lower Leg out from the thigh.
4. Discard the old IO Leg according to specified guidelines.

Reverse the above process to attach the new IO Leg.

## System Setup

### Problem

- Lost data or total system failure (General System failure).

### Possible Solution

- Should system shutdown occur or all data is lost or corrupted; please contact your local Laerdal Service Center.

## LLEAP

- See Software Help Files > LLEAP > Troubleshooting Guide (Help is accessible from Laerdal Simulation Home)

## Laerdal Patient Monitor

- See Software Help Files > Laerdal Patient Monitor > Troubleshooting Guide. (Help is accessible from Laerdal Simulation Home)

## Debriefing and video recording

- See section for Troubleshooting in SimView Help (accessible from a link at the bottom of the screen in SimView) or Troubleshooting Guide from SessionViewer Help (accessible from Laerdal Simulation Home)

## Simulator

### Problem

- Unpredictable behavior.

### Possible Solutions

- Simulator malfunctions may be caused by loose cables, tubes or connectors. Open the torso and check if any items appear to have become disconnected or are leaking. See section: Simulator Setup - Opening the Torso.
- In case of fluid leakage, power off the Simulator and contact Laerdal Technical Service.

## Airway Contamination

### Problem

- Simulator airways have become contaminated from mouth-to-mouth rescue breathing.

### Possible Solution

- Clean the outside of the Simulator with Simulator wipes. Clean the inside of the oral cavity with Simulator wipes. Change the lung bladder; see section: Maintenance - Replacing Lung Bladder.

## Chest Movement

### Problem

- No chest rise on Simulator.

### Possible Solutions

- Check that power to Simulator is ON.
- In SimPad or LLEAP, check that
  - awRR is not set to zero
  - Chest Rise is not set to OFF
  - ECG is set to a perfusing rhythm.
  - Airway complications like maximum air resistance or laryngospasm are not set.
  - Verify that the Integrated Compressor is enabled in LLEAP. See section: Turning the Internal Compressor Off.
- Check that any external compressed air source is switched off and that the air tube is disconnected from the Simulator.
- The internal compressor may be overheated. Wait approximately 20 minutes for it to cool down.
- Chest rise is set to bilateral (for example if ET-tube is inserted too far into the bronchia).
- Chest rise bladder is leaking or tubing to chest rise bladder is twisted kinked or disconnected. Replace chest rise bladder if it is leaking, see Maintenance: *Replacing the Chest Rise Bladder*.
- Check air tubing for leakage; check that all connections are intact. Replacing tubing section if leaking.
- Shallow chest movements and the internal compressor runs continuously. The internal compressor may be worn. Contact Laerdal Technical Service.

## Lungs

### Problem

- Lungs not functioning properly

### Possible Solutions

- Open the torso and chest plate. Check that the lungs are free to expand and are not restricted by any cables.
- Check that the Lung Bladder is properly connected, and that the tubes are not twisted.
- Check that the Lung Bladder is in a horizontal position and inserted correctly.

- Check for flaws or tearing in the lung bladder.
- Check that there are no obstructions inside the Simulator airways which may block air flow.
- If there is no change when adjusting lung compliance, contact Laerdal Technical Service.
- If there is no change in lung resistance, contact Laerdal Technical Service.

## Mechanical Noise during Auscultation

In LLEAP click <Auscultation focus>.

## Pulses (Radial and Brachial)

### Problem

- Cannot feel pulses

### Possible Solutions

- Skin may be too tight over pulse units – re-adjust skin and reboot.
- Ensure blood pressure is set to at least 60 systolic.
- Ensure compressor is ON.

**Note:** If the radial pulse is set below 60, peripheral pulses may disappear.

## Simulator Shutdown

### Problem

- The Simulator is unresponsive.

### Solution

1. Power off the Link Box.
2. Turn the SimPad off, or if you are using LLEAP turn the instructor PC off.
3. Restart everything.

If the Simulator is still unresponsive, contact your Laerdal representative.

## Overall Dimensions

Length / Width (Simulator only): 45 in x 9.84 in (120 cm x 42.5 cm)  
 Weight (Simulator only): 25 lbs (11.36 kg)

## Simulator Power

External power: Input voltage 9VDC, 1.5A  
 Internal battery: 7.4V, 4.4Ah, Lithium-Ion

Only use approved SimJunior external power supply and batteries.

## Air Pressure

Internal air tank: Max 6 psi  
 External air connection: Max 16 psi

## Temperature Limits

Operating temperatures: +4 °C to 40 °C (39 °F to 104 °F)  
 Storage temperatures: -15 °C to 50 °C (5 °F to 122 °F)

## Environment - Simulator only

Relative humidity: 20% -90% (non-condensing)  
 DO NOT use outdoors in wet conditions.  
 Not tested with salt spray.

## RF Communication

Operation range: 10 m (30 ft) max.

## Material Chart for Simulator

Clothes:	Cotton, Nylon
Skins and airways:	PVC
External hard plastics:	PVC, ABS
Inner plastics:	Silicone, TPU, PVC, Nitrile ABS, POM, Nylon + GF Epoxy-Polyurethane
Metal components:	Aluminum, Brass, Steel

## Minimum Computer Requirements

- Intel i-core generation 3 or newer
- Score over 3000 PassMark - CPU Mark
- 4 GB RAM
- 120GB hard disk space
- 1366x768 or better
- Headset w/microphone
- USB ports – 2

## Minimum Software Requirements

- Windows 7 or Windows 8
- 100% DPI

## Cleaning Fluids

To clean the Simulator use one of the following:

- 60% Isopropanol alcohol
- 70% Ethanol
- Mild solution of liquid soap and water

## IV Fluids

Only use distilled or deionized water to simulate IV and IO fluids.

## Spare Parts and Accessories

For latest version of Spare Parts and Accessories, visit [www.laerdal.com](http://www.laerdal.com)

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© 2022 Laerdal Medical AS. All rights reserved.
Manufacturer: Laerdal Medical Corporation