CPRmeter™ with Q-CPR® Technology

Directions for Use

ENGLISH	Directions for Use	5
DEUTSCH	Gebrauchsanweisung	31
ITALIANO	Istruzioni per l'uso	55
ESPAÑOL	Instrucciones de utilización	79
NORSK	Bruksanvisning	103
SVENSKA	Bruksanvisning	127
FRANÇAIS	Mode d'emploi	151
NEDERLANDS	Handleiding	175
DANSK	Brugsanvisning	199
SUOMI	Käyttöohje	223

CPRmeter™ Directions for Use

Contents

1 - Description and Indications for Use	7
2 - Setup	8
3 - Use during an Emergency	10
4 - Compression Feedback	12
5 - Debriefing after Use	15
6 - Maintenance and Cleaning	17
7 - Troubleshooting Guide	20
8 - Accessories, Consumables and Spare Parts	21
9 - Glossary of Symbols	23
10 - Additional Warnings and Cautions	24
11 - Recommendations	25
12 - Specifications	26
13 - Environmental Considerations	28
14 - Electromagnetic Conformity	29

About this edition

The information in this *Directions for Use* applies to the model 801-00140 CPRmeter™. This information is subject to change. Please contact your local Laerdal representative for information on revisions.

Printed in Norway

Copyright

© 2013 Laerdal Medical AS. All rights reserved.

Device Manufacturer

Laerdal Medical AS, Tanke Svilandsgate 30, PO Box 377, 4002 Stavanger, Norway, (+47) 51511700.

Distributed in the USA by Laerdal Medical Corporation, 167 Myers Corners Rd, Wappingers Falls, NY 12590 (877) 523-7325

The CPRmeter™ with Q-CPR® technology is protected under U.S. Patents 6,306,107, 6,351,671, 6,390,996, 7,074,199, 7,108,665, 7,118,542, 7,122,014, 7,220,235 and other foreign patents. Other U.S. and foreign patents pending. The design of CPRmeter™ is protected under design registrations/design patents US609,813; US628,212; EU997416, EU997424, NO81223, NO81219, AU323899, AU324044 and other foreign design registrations.

CPRmeter™ and Q-CPR® are trademarks or registered trademark of Laerdal Medical AS. Energizer® is a registered trademark of Eveready Battery Company, Inc. Panasonic® is a registered trademark of Matsushita Electric Industrial Co., Ltd. Duracell® is a registered trademark of The Gillette Company or its affiliates. SD™ and microSD™ are trademarks of SD-3C LLC.


CPRmeter Overview

Front View


(this surface up)

Rear View

(apply this surface to patient)

Side View

1 Description and Indications for Use

The CPRmeter with Q-CPR® technology is a small, lightweight device powered by a replaceable battery. The CPRmeter is intended for use by responders who have been trained in CPR and use of the CPRmeter.

The CPRmeter is used as a guide in administering cardiopulmonary resuscitation (CPR) to a suspected sudden cardiac arrest (SCA) victim at least 8 years old. If in doubt about the appropriateness for use, perform CPR without using the CPRmeter.

When attached to the bare chest of a suspected victim of SCA, the CPRmeter provides real-time feedback on CPR compressions in accordance with current CPR guidelines. It displays CPR feedback indicators for depth, release, and rate of chest compressions. It also counts the number of compressions in a series, and provides notification of lack of expected CPR activity.

Rx only

Caution: Federal law (USA) restricts the CPRmeter to sale by or on the order of a physician.

⚠ WARNING: The CPRmeter is not intended for use on SCA victims under 8 years old.

💬 IMPORTANT NOTE: CPR cannot assure survival, no matter how well it is performed. In some patients, the underlying problem causing the cardiac arrest is not survivable despite any available care.

2 Setup


Check the contents of the CPRmeter box. It should contain:

- the CPRmeter
- Red sleeve
- set of 3 CPRmeter Patient Adhesives in a resealable pouch
- Directions for Use
- Laerdal Global Warranty and WEEE-statement leaflets
- microSD™ Memory Card with SD™ Card Adapter
- 3V Lithium 123 (Li/MnO₂) battery

 Remove the protective film from the CPRmeter display screen before first use.

Inserting a new Battery

1. Unscrew the Rear Cover using a flathead screwdriver (not supplied) or a coin, and lift out the cover.
2. Insert a new battery into the battery well, in the orientation as directed by the battery symbol inside the well.
3. Ensure that the Rear Cover's vent membrane is not soiled or damaged. Refer to *Chapter 6, Replacing the Rear Cover* for instructions.
4. Position the Rear Cover on the CPRmeter. Thread the screws through the washers and into the cover; and tighten the screws.


The CPRmeter's status light will flash green every 5 seconds to indicate that it has passed its self-test and the battery is OK for use.

If for any reason the CPRmeter's status indicator light is orange or does not flash green, refer to *Chapter 7, Troubleshooting Guide* for instructions.

Inserting a microSD™ Memory Card


Use of a microSD™ Memory Card with the CPRmeter is optional. Refer to *Chapter 5, Q-CPR® Review* for further details.

1. Remove the Rear Cover and Battery.
2. Insert a microSD™ card into the card slot inside the battery well, in the orientation as directed by the card symbol inside the well.
3. Re-insert the battery and attach the Rear Cover as described in the previous section.

Applying a new Patient Adhesive

1. Check the label on the CPRmeter patient adhesives resealable pouch to be sure that they are within their expiration date.
2. Open the package of CPRmeter patient adhesives and peel one from the white liner strip to expose the adhesive surface on its underside.
3. Align the bottom of the patient adhesive with the yellow patient adhesive mount area on the Rear Cover of the CPRmeter. Make sure the channel on the adhesive is directly over the vent membrane. Press the patient adhesive into place.
4. Do not peel off the green liner from the attached adhesive until you are ready to apply the CPRmeter to a patient in an emergency. Return remaining patient adhesives to their resealable pouch.




Storing the CPRmeter between Use

Use the CPRmeter red sleeve to shield the display screen from scratches and to protect the patient adhesive from damage. Ensure that the On/Off button cannot be inadvertently activated during storage. Perform CPRmeter Routine Maintenance (Section 6) at least once every six months while the CPRmeter is in storage.

3 Use During an Emergency

Attach the CPRmeter to the Patient's Chest

1. Make certain the patient is on a firm surface.
2. Remove clothing from the patient's chest.
3. Remove the CPRmeter from its red sleeve.
4. Turn on the CPRmeter by pressing the On / Off button once.
5. If the patient's chest is wet, dry it before attaching the CPRmeter.
6. Peel off the green liner from the CPRmeter patient adhesive to expose the white adhesive surface.
7. Position the CPRmeter so that the compression area is on the lower half of the sternum (breastbone), on the centerline of the bare chest, as illustrated on the front of the CPRmeter.
8. If the CPRmeter moves during use, re-position it to the center of the chest, as described above.
9. If difficulty is encountered in applying the CPRmeter, do not delay initiation of CPR. Remove the CPRmeter and begin compressions.


⚠ WARNING: Do not use the CPRmeter in conjunction with any mechanical or automated compression device.

⚠ WARNING: Do not use the CPRmeter on top of defibrillation pads, unless the manufacturer of the defibrillator and the defibrillation pads has explicitly stated that the CPRmeter can be used in such manner.


Perform CPR. Follow Feedback.

Using standard CPR technique, place the heel of one hand directly onto the center of the compression area of the attached CPRmeter. Place your other hand on top of the first.

You should be able to see the display area of the CPRmeter to look for feedback.


1. Provide chest compressions according to your CPR protocol.
2. When compressions are first detected by the CPRmeter; the display changes to enlarged compression indicators, as shown:

3. Follow feedback provided by the indicators on the CPRmeter display. The indicators are described in the following section. If no CPRmeter feedback is observed, remove the CPRmeter and begin chest compressions.

⚠ WARNING: When the CPRmeter is used together with a defibrillator, make sure to follow the defibrillator manufacturer's instructions. Stop compressions, remove hands from the CPRmeter and remain clear of all patient contact during defibrillation or when otherwise required, in accordance with a proper defibrillation protocol.

4 Compression Feedback


The display of the CPRmeter provides real-time feedback to the responder as described below.

Good Compressions

Each compression performed is represented on the CPRmeter display by a moving white compression depth indicator bar. When the CPRmeter detects that a compression meets the targets for depth or release, the target briefly illuminates.


If the CPRmeter detects that the rate of compressions is within the target area, the "speedometer" needle on the compression rate indicator points to the green target area and the green target briefly illuminates.

Compression Depth

If the CPRmeter detects a compression that fails to meet the proper minimum compression depth target (50 mm / 2 inches – for a patient on a hard surface), the compression depth target does not illuminate.

If 4 consecutive compressions fail to meet the compression depth target, the display shows a yellow arrow pointing to the compression depth target.


If the CPR meter detects a compression that exceeds 60 mm (2.36"), the CPRmeter shows the depth indicator below the target area. If a specific CPR event requires CPR to be performed on a patient lying on a mattress, slide a backboard under the patient and compensate for the mattress softness by ensuring that for each chest compression the area below the compression depth target lights up.

⚠ WARNING: When performing CPR on a patient lying on a mattress, a backboard must be used to limit the amount of compressed depth which is absorbed by the mattress. Depending on characteristics of the mattress, backboard and patient, the depth compensation does not guarantee that the patient chest is compressed by 50 mm / 2 inches.

Incomplete Release

If the CPRmeter detects that pressure is not fully released between compressions, the release target does not illuminate. If 4 consecutive compressions fail to meet the release target, the display shows a yellow arrow pointing to the compression release target. Release pressure completely between compressions.

IMPORTANT NOTE: Release pressure completely between compressions.


Compression Rate

If the CPRmeter detects that the compression rate is greater than the target rate, the needle on the compression rate indicator points to the right of the green area. If the CPRmeter detects that the compression rate is less than the target rate, the needle on the compression rate indicator points to the left of the green area.

Compression Counter

When compressions start, the CPRmeter display shows a numerical compression counter. The counter digits turn solid white between 25 and 30 compressions, as a guide to the responder when performing a cycle of 30 compressions. Beyond 30 compressions, the counter digits flash solid white for every tenth compression*. The compression counter is reset after 3 seconds without a compression.

Compression Counter

* When performing continuous compression CPR on an intubated patient, the compression counter may be used to time delivery of ventilations. At a compression rate of 100 per minute: If a ventilation is delivered for every 10 continuous compressions, ventilation rate should be 10 per minute.

Inactivity

If compressions are stopped during a CPR event, after 3 seconds the CPRmeter displays an inactivity timer which counts the seconds since the last compression.

After 20 seconds since the last compression, the inactivity timer starts flashing.

After 1 minute, the CPRmeter display fades down to conserve battery power. The display is restored when a new compression is delivered.

After 10 minutes of inactivity, the CPRmeter is automatically turned off. To restart the CPRmeter, press the On/Off button.

Turning Off the CPRmeter


To turn the CPRmeter off, press and hold the On / Off button for at least 1 second.

5 Debriefing After Use

Q-CPR® Quick Review


The CPRmeter can display CPR performance statistics for the last CPR event. When the CPRmeter is turned on, press the On / Off button once to activate Q-CPR Quick Review. The statistics are shown over two displays.

- Percentage of compressions with release force within the Compression Release Target.
- Percentage of compressions with compression depth which have reached the Compression Depth Target.
- Percentage of compressions with compression rate within the Compression Rate Target.

Press the On / Off button once to cycle to the next display.

- Duration of the CPR event (minutes: seconds).
- Percentage of the CPR event duration where chest compressions have been delivered.

Press the On / Off button once to cycle between the displays.

To turn the CPRmeter off, press and hold the On / Off button for at least 1 second.

The CPRmeter reverts to Compression Feedback mode if a compression is delivered.

The CPR event statistics are stored when the CPRmeter is turned off. When turned on again, the statistics from the stored CPR event can be reviewed, as described above.

When the CPRmeter is used in a new CPR event, the preceding event's statistics are deleted and the new event's statistics are stored in Q-CPR Quick Review.

CPR performance statistics are only calculated if at least 10 compressions have been delivered.

Q-CPR® Review

When a microSD™ memory card has been inserted into the CPRmeter prior to use, data from CPR events are stored on the card. The data can only be viewed on a PC with the Laerdal Q-CPR Review software program, version 3.1 or newer, installed.

Q-CPR Review software is available for download by visiting: www.laerdal.com/downloads


To transfer CPRmeter data to a PC with a SD™ Memory Card reader:

1. Remove the microSD™ card from the CPRmeter and insert it into the provided SD™ Card Adapter.
2. Insert the SD™ Card Adapter into the SD™ card reader.
3. Start the Q-CPR Review software on the PC and follow instructions.

Remember to re-insert a microSD™ card into the CPRmeter before further CPR event data recording.

If the 'remove memory card' symbol is shown on the CPRmeter display, the memory card is full and should be replaced to prevent data from subsequent CPR events from being lost.

The CPRmeter's compression feedback and Q-CPR Quick Review features will function as normal even if the Memory Card is full.

6 Maintenance and Cleaning

Routine Maintenance


On a routine basis, at least once every six months, check the following:

1. Verify that the CPRmeter status light is flashing green. If the status light is not flashing green every 5 seconds: See *Chapter 7, Troubleshooting Guide*.
2. Check that the CPRmeter has a patient adhesive in place and that the liner remains on it. Replace the patient adhesive at least every 2 years if it is not used.
3. Replace the battery in accordance with Battery Monitoring (below), and at least every 2 years.
4. Replace the Patient Adhesive every time the Rear Cover is opened.

Battery Monitoring

The CPRmeter continuously monitors the power of its battery. If the remaining power is less than that required for an entire CPR event, the following visual indicators signal that the battery should be replaced before the next use:

- The green status light does NOT flash green every 5 seconds when the CPRmeter is turned off.
- A small battery-low icon is shown in the upper right corner of the display when the CPRmeter is being turned on.
- A large battery-low icon is shown on the display when the CPRmeter is being turned off.

IMPORTANT NOTE: If the remaining battery power during use becomes too low to sustain further operation, the battery-low icon is shown for 10 seconds and then the CPRmeter turns itself off.

WARNING: Do not interrupt CPR to replace the battery. Continue CPR without feedback from the CPRmeter.

Battery Replacement

1. Remove the CPRmeter Patient Adhesive from the back of the CPRmeter.
2. Unscrew the Rear Cover using a flathead screwdriver (not supplied) or a coin, and lift out the cover.
3. Remove the old battery and discard.
4. Insert a new battery into the battery well, in the orientation as directed by the battery symbol inside the well. See *Chapter 12, Specifications*, for battery specifications.
5. Ensure that the Rear Cover's vent membrane is not soiled or broken.
6. Position the Rear Cover on the CPRmeter. Thread the screws through the

washers and into the cover, and tighten the screws.

7. Apply a new CPRmeter patient adhesive to the CPRmeter as described in *Chapter 2, Setup*.

After each Use

After use on a patient, the CPRmeter may be contaminated and should be handled appropriately.

1. Place the contaminated CPRmeter in a plastic bag until it can be cleaned. Do not insert a contaminated CPRmeter into the red sleeve.
2. If it is visibly soiled, wipe the CPRmeter with a soft cloth or paper towel to remove as much contamination as possible.
3. Remove the CPRmeter patient adhesive from the back of the CPRmeter.
4. Clean the CPRmeter as described under *Cleaning*.
5. Examine the vent membrane in the center of the yellow Rear Cover. If the vent membrane is soiled or damaged, see *Replacing the Rear Cover*.
6. Check the exterior of the CPRmeter for signs of damage. Contact Laerdal to arrange for replacement if needed.
7. Apply a new CPRmeter patient adhesive to the CPRmeter as described in *Chapter 2, Setup*.

Cleaning

If the CPRmeter has been used in a training situation, it may be wiped using an alcohol wipe (70% ethanol solution).

If the CPRmeter has been used in a clinical situation, clean it as follows:

1. Clean and scrub the exterior using a mild detergent and either a soft cloth or toothbrush until the surfaces are visually clean.
2. Wipe the exterior with a soft cloth dampened in lukewarm water.
3. Clean the exterior using a 0.55% solution of ortho-phthalaldehyde. Spray the solution on to cover all exterior surfaces, and allow to sit for a minimum of 10 minutes. An alternative cleaning agent is isopropyl alcohol (70% solution).
4. Wipe the exterior with a clean soft cloth dipped in water. Allow to dry completely.

⚠ CAUTION: This cleaning procedure ensures that the CPRmeter is clean; however it does not ensure disinfection.

⚠ CAUTION: Do not immerse the CPRmeter in water, hold it under running water, or allow moisture to penetrate it. Do not sterilize the CPRmeter.

Replacing the Rear Cover

The vent membrane allows the CPRmeter to adapt to minute air pressure changes while protecting it from fluid entry. If the vent membrane in the center of the yellow Rear Cover is soiled, contaminated, scratched or damaged, the Rear Cover must be replaced. Do not attempt to clean the vent membrane.

⚠ WARNING: A soiled or damaged vent membrane may lead to inaccurate feedback or fluid entry into the CPRmeter. If there are signs of fluid entry, remove the CPRmeter from service and contact your Laerdal representative for advice.

If the Rear Cover is damaged or does not fit flush with the CPRmeter housing, or the Rear Cover screws have sharp edges or the washers under each of the Rear Cover screws are missing, the Rear Cover must be replaced. See *Chapter 8, Accessory, Consumables and Spare Parts*.

1. Unscrew the Rear Cover using a flathead screwdriver (not supplied) or a coin, and lift out the cover.
2. Remove the replacement Rear Cover from its packaging and check to be sure it is not damaged.
3. Position the Rear Cover on the CPRmeter. Thread the screws through the washers and into the cover, and tighten the screws.

Customer Service Indicator

The Customer Service Indicator (right) appears on the CPRmeter only at shutdown, after 500 000 chest compressions. Contact your local Laerdal representative for further instructions.

7 Troubleshooting Guide

Problem	Suggested Action
The CPRmeter display is dark.	Make sure the CPRmeter is turned on.
CPRmeter status light does not flash green.	Replace the battery. Ensure that the battery is inserted in the correct orientation.
CPRmeter status light is orange (solid or flashing) and the CPRmeter display is dark.	Remove the CPRmeter from service. Contact Laerdal for technical support.
CPRmeter does not adhere to patient's chest.	If there is no other responder present, provide CPR using the CPRmeter even if it does not adhere to the patient's chest. Be sure to keep the CPRmeter properly positioned.
	If there is another responder present to provide CPR, remove the CPRmeter and replace the CPRmeter patient adhesive.
	Dry the patient's chest if wet and re-apply the CPRmeter with a minimum interruption to CPR.
User is unfamiliar with the CPRmeter; or thinks there is a problem with the CPRmeter.	Ignore feedback from the CPRmeter and continue with CPR. At an appropriate time after the CPR event, inspect the vent membrane for soiling or damage. Apply a new patient adhesive and let a trained user perform chest compressions with the CPRmeter on a CPR manikin. If the CPRmeter does not appear to function as expected, remove it from service and contact Laerdal or approved Laerdal representative for technical support.

8 Accessories, Consumables and Spare Parts

Contact Laerdal to order the following replacements for the CPRmeter:

- 10 replacement packs of CPRmeter patient adhesives (each pack holds 3)
[REF: 801-10850]

- Replacement Rear Cover with screws and rubber washers
[REF: 801-10750] (spare recommended)

- CPRmeter Red Sleeve
[REF: 801-10150]

- CPRmeter Hard Case
[REF: 801-10550]
The accessory CPRmeter Hard Case for holding the CPRmeter can be attached to a belt or a clip hook

- Reusable Silicone Cover kit, includes one pack of 3 large patient adhesives
[REF: 801-10650]

- 10 replacements packs of CPRmeter large patient adhesives (each pack holds 3)
[REF: 801-10950]
- CPRmeter battery (each pack holds 5)
[REF: 801-10350]
- microSD card kit (each pack holds 5)
[REF: 801-10450]


Reusable Silicone Cover Kit

The CPRmeter can be fitted with the accessory Reusable Silicone Cover and the large patient adhesive (single-use). This increases the contact area to the patient.

Applying the Reusable Silicone Cover and large patient adhesive

1. Check that the silicone cover has been cleaned and disinfected after last use.
2. Insert the CPRmeter into the sleeve.
3. Check the label on the CPRmeter large patient adhesives resealable pouch to be sure that they are within their expiration date.
4. Open the package of CPRmeter large patient adhesives and peel one from the white liner strip to expose the adhesive surface on its underside.
5. Align the bottom of the patient adhesive with the yellow patient adhesive mount area on the Rear Cover of the CPRmeter and the adjacent Silicone Cover's bottom surface. Make sure the channel on the adhesive is directly over the vent membrane. Press the large patient adhesive into place.
6. Do not peel off the green liner from the attached large patient adhesive until you are ready to apply the CPRmeter to a patient in an emergency. Return remaining patient adhesives to their resealable pouch.

IMPORTANT NOTE: Do not delay CPR. If the Reusable Silicone Cover is generally used with the CPRmeter, the cover and the Large Patient Adhesives should be applied to the CPRmeter in advance.

Use and maintenance

The CPRmeter use or maintenance is not changed by use of the Silicone Cover. Clean the Silicone Cover separately from the CPRmeter. Disinfect the Silicone Cover using the following method:

1. Clean and scrub the cover with a small soft brush, using a mild detergent until the cover is visibly clean.
2. Rinse in lukewarm water.
3. Submerge the cover in a 0.55% solution of ortho-phthalaldehyde in accordance with manufacturer's specifications.
4. Rinse again and dry.

The Silicone Cover can also be autoclaved in distilled water at 136 °C (210 °F) and 2.0 kg/cm² for 10-20 minutes.

WARNING: The CPRmeter with the Silicone Cover can not be considered as disinfected or sterilized.

9 Glossary of Symbols

The following symbols appear on the CPRmeter; accessories and packaging:

Symbol	Definition
	The product is in compliance with the essential requirements of Council Directive 93/42/EEC as amended by Council Directive 2007/47/EC.
	Compliance with applicable U.S and Canadian safety standards has been certified by Canadian Standards Association.
	These CPRmeter patient adhesives are disposable and are for single patient use only. Do not re-use. Re-use will lead to increased risk of cross contamination, and/or degradation of adhesive performance.
	Defibrillation protection. The CPRmeter is defibrillation protected, type BF patient connection.
	Manufacturer:
	Dispose of in accordance with your country's requirements.
	Reference order number:
	The CPRmeter meets IEC 60529 class IP55.
	Serial number:
	Expiration date for patient adhesives, formatted YYYY-MM (year-month).
	The patient adhesives, CPRmeter or other parts do not contain natural rubber latex.
	Warning / Caution symbol
	Store the CPRmeter patient adhesives between the temperatures shown. See Chapter 12, Specifications.
	Contains number of CPRmeter patient adhesives shown as "#."
	Refer to Directions for use.
	Consult Directions for Use
	Lift here to peel off the patient adhesive liner and apply to patient's bare chest.
	Not for use on children under 8 years old
	microSD memory card Battery – 3V Lithium 123 type
Rx only	Federal law (US) restricts this device to sale by or on the order of a physician

10 Additional Warnings and Cautions

A Warning identifies conditions, hazards, or unsafe practices that can result in serious personal injury or death. A Caution identifies conditions, hazards, or unsafe practices that can result in minor personal injury or damage to the CPRmeter.

- WARNING: The CPRmeter is not intended for use in a moving environment, such as an ambulance. If used during patient transport, the CPRmeter may provide inaccurate feedback. If CPR is indicated in a moving environment, do not rely on the CPRmeter depth feedback during such conditions. It is not necessary to remove the device from the patient.
- WARNING: Do not practice by using the CPRmeter on a person. The CPRmeter may be used with a training manikin or simply on a compliant surface for practice.
- WARNING: Properly performed CPR may result in fracturing of the patient's ribs.* If rib integrity has been compromised, continue to provide CPR in accordance with your local protocol.
- WARNING: Properly performed CPR may result in chest injuries* e.g. external chest wall bruising or abrasion.
- WARNING: Do not rely on CPRmeter feedback during aircraft ascent and descent, as its accuracy is reduced in such conditions.
- CAUTION: Do not apply the CPRmeter to an open wound or recent incision site.
- CAUTION: The CPRmeter is designed to be used only with Laerdal-approved accessories. The CPRmeter may perform improperly if non-approved accessories are used. Do not attempt to modify the CPRmeter in any way.
- CAUTION: Use only model 801-10850 Patient Adhesives with the CPRmeter. If the Reusable Silicone Cover is applied to the CPRmeter, use only model 801-10950 Large Patient Adhesives.

*Black CJ, Busuttill A, Robertson C. Chest wall injuries following cardiopulmonary resuscitation. *Resuscitation*. 2004;63:339 –343.

11 Recommendations

To prolong the life of the display, avoid storing the CPRmeter where it is exposed to direct sunlight when not in use.

Responders should receive training, including regular refresher training, in use of the CPRmeter. When training with the CPRmeter on a CPR manikin, disable or ignore feedback from the manikin.

12 Specifications

CPR Targets

Category	Specification
Compression Depth Target	> 50 mm (2") <i>Depth accuracy: ±10%</i>
Compression Release Target Force	< 2.5 kg (5.5 lbs) Force accuracy: +1.5 kg, -2.0 kg (+3.3 lbs, - 4.4 lbs)
Compression Rate Target	100 to 120/min ± 3/min

CPRmeter [REF 801-00140]

The CPRmeter meets the performance requirements of IEC 60601-1, 2nd and 3rd edition.

Category	Specification
Dimensions	154 mm × 64 mm × 28 mm (6.1" × 2.5" × 1.1")
Weight	227 g (7.3 oz)
Battery	3V Lithium 123 (Li/MnO ₂) The CPRmeter has been tested with Energizer® batteries. Use only Energizer®, Panasonic® or Duracell® batteries.
microSD™ Memory Card	The CPRmeter has been tested with 2Gb microSD™ memory cards. Always test a new card before using it in data collection.
Operating temperature	0°C to 40°C (32°F to 104°F) (Fault condition: If the battery has a short circuit at maximum ambient temperature, the unit surface temperature may reach 60°C (140°F). In this condition the unit will be inoperable.)
Storage temperature	-20°C to 70°C (-4°F to 158°F)
Operating Relative humidity	5% to 95%
Storage Relative humidity	5% to 75%
Operating/storage atmospheric pressure	1014 to 572 hPa (101 to 57 kPa)
IP protection class per ISO/IEC 60529	IP 55

Electromagnetic Compatibility	Meets IEC 60601-1-2 and RTCA/DO-160E.
Recovery time	Recovery time after defibrillation: 0 s

CPRmeter Patient Adhesives [REF 801-10850]

Category	Specification
Dimensions	39 mm × 90 mm (1.5" × 3.5")
Temperature and Relative Humidity	Storage temperature: -20°C to 70°C (-4°F to 158°F) Relative humidity 0% to 75% Operating temperature: 0°C to 50°C (32°F to 122°F) Relative humidity 0% to 95%.
Material	Foam pad with biocompatible adhesive on each side.
Shelf Life	2 years when applied to the CPRmeter or 4 years in unopened package. Do not exceed the expiration date on the packaging.

CPRmeter Large Patient Adhesives [REF 801-10950]

Category	Specification
Dimensions	64 mm × 128 mm (2.5" × 5").

Reusable Silicone Cover [REF 801-10650]

Category	Specification
Dimensions	66 mm × 156 mm × 31 mm (2.6" × 6.1" × 1.2")
Weight	36.5 g (1.3 oz)
Material	Silicone

Warranty

The Laerdal CPRmeter has a one-year Limited Warranty. Refer to the Laerdal Global Warranty for terms and conditions.

13 Environmental Considerations

Product	Information
CPRmeter	The CPRmeter contains electronic components. Dispose of it at an appropriate recycling facility in accordance with local regulations.
CPRmeter Patient Adhesive	The used patient adhesive may be contaminated with body tissue, fluid, or blood. Dispose of it as infectious waste.

14 Electromagnetic Conformity

Guidance and manufacturer's declaration: The CPRmeter is intended for use in the electromagnetic environment specified in the tables below. The user of the CPRmeter should assure that it is used in such an environment.

Electromagnetic Emissions

Emissions Test	Compliance	Electromagnetic environment - Guidance
RF CISPR 11	Group 1 Class	The CPRmeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The CPRmeter is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Electromagnetic Immunity

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	There are no special requirements with respect to electrostatic discharge.
Power Frequency (50/60/400 Hz) Magnetic Field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial/hospital environment. There are no special requirements for non-commercial/non-hospital environments.
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	10 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the CPRmeter, than is absolutely necessary.†,‡The recommended separation distances for various transmitters and the CPRmeter are shown in the following table. Interference may occur in the vicinity of equipment marked with the following symbol: 

† The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,660 MHz to 40,700 MHz.

‡ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CPRmeter is used exceeds the applicable RF compliance level above, the CPRmeter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the CPRmeter.

Recommended separation distances between portable and mobile RF communications equipment and the CPRmeter

The CPRmeter is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CPRmeter can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CPRmeter as recommended below, according to the maximum output power of the communications equipment.

Electromagnetic Emissions

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (M)	
	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0,12	0,23
0.1	0,38	0,72
1	1,2	2,3
10	3,8	7,28
100	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter; where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1. At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2. The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,660 MHz to 40,700 MHz.

NOTE 3. An additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.

NOTE 4. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

NOTE 5. Transmitters/antennas of this power-level are most likely mounted on an emergency vehicle chassis. The distances cited here are for open field. For an external antenna, the separation distance is most likely shorter.

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CPRmeter™ e Q-CPR® sono marchi commerciali registrati o depositati di Laerdal Medical AS. Energizer® is a registered trademark of Eveready Battery Company, Inc. Panasonic® è un marchio commerciale registrato di Matsushita Electric Industrial Co., Ltd. Duracell® è un marchio commerciale registrato di The Gillette Company o delle sue affiliate. SD™ e microSD™ sono marchi commerciali di SD-3C LLC.