EN

FR

DE

ES

IT

NL

NO

SV

DA

FI

ZH

JA

KO

Premature Anne & Premature Anne Task Trainer

Important Product Information

www.laerdal.com

Laerdal

helping save lives

Read these instructions thoroughly.

Observe all warnings, precautions and instructions on the product, in the User Guide and in this Important Product Information booklet.

Retain this booklet for future reference.

Warnings and Cautions

A Warning states a condition, hazard, or unsafe practice that can result in serious personal injury or death.

A Caution states a condition, hazard, or unsafe practice that can result in minor personal injury or damage to the manikin.

Umbilical Cord

Cautions

- Do not use real umbilical cords during simulation with the Premature Anne.
- Do not insert the simulated umbilical cord further in than 1 cm (due to risk of spillage into the internal parts).

Use of Fluids During Simulation

Cautions

- Do not apply any fluids directly into the mouth, airways and stomach of the manikin.
- Only the following fluids can be used during simulation:
 - Cat. No. 250-21050 Laerdal Airway Lubricant. Apply sparingly on the ET tube and laryngoscope prior to intubation.
 - Cat. No. 377-18850/377-19950 Laerdal Birthing Simulation Lubricant. Apply on the front side of the manikin only.
 - Cat. No. 240-00250 Laerdal Umbilical Blood. Insert 3 ml into the umbilical reservoir prior to simulation.
 - Water (for simulated volume expander applied through the umbilical cord - max 2 x 7.5 ml).

During Use

Cautions

- Do not overfill the umbilical reservoir (due to risk of spillage into the internal parts).
- Do not use tools exceeding the sizes and dimensions specified in the User Guide (due to risk of damage to the product).
- Do not use the manikin within 2.5m of a patient.
- It is recommended that a "soft" 2.5mm nasal ET tube be used for the practice of nasal intubation.
- Do not use the accessories supplied with this product on real patients. The accessories are intended for training only and not for clinical use.

Cautions

Do not attempt to perform the following on the manikin:

- Mouth-to-mouth or mouth-to-mask ventilation.
- Insertion of simulated vomit or fluids for suctioning.
- Insertion of simulated surfactant or nutrition.
- Insertion of fluids into the peripheral IV ports.
- Defibrillation.

Use of Medical Devices

Cautions

Do not use medical devices during simulation with the manikin unless the medical devices have been subjected to the local cleaning and disinfection procedures both prior to and immediately after simulation (due to risk of cross-contamination and infections).

Care, Maintenance and Cleaning

- Use mild soap and water for cleaning of the skin using a dampened soft cloth.
- In case of simulation scenarios involving the use of fluids in the umbilical reservoir, drain the reservoir immediately following the simulation.
- In case of simulation scenarios involving the use of birthing lubricant on the front of the manikin, clean skin immediately following the simulation.
- Do not submerge the product.

These can discolour the manikin:

- Pigments from lipstick and pens
- Latex gloves
- Using clothing other than that provided with the Premature Anne

Premature Anne - Battery and Charging

Cautions

- Do not attempt to disassemble the product. There are no user serviceable parts in Premature Anne.
- Do not use the battery in temperatures above what is specified in this IPI.
- Do not heat or incinerate.
- Do not crush the battery.
- Do not short circuit the battery terminals.
- Do not immerse in water.
- Use only the AC charger supplied with the product.
- Do not use the mini software upgrade USB to charge the manikin. The mini cable is only for software updates with SimPad PLUS.
- Do not charge Premature Anne when it is located inside an incubator.
- Do not charge Premature Anne in areas with a high oxygen concentration.

Federal Communications Commission (FCC) Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- *This device may not cause harmful interference, and*
- *This device must accept any interference received, including interference that may cause undesired operation.*

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- *Reorient or relocate the receiving antenna*
- *Increase the separation between the equipment and the receiver*
- *Connect the equipment into an outlet on a circuit different from that to which the receiver is connected*
- *Consult the dealer or an experienced radio/TV technician for help.*

The user is cautioned that any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

The radio module is located in the skull and the FCC ID is referenced on the manikin label.

FCC ID 2AAQS-ISP091201

Industry Canada Rules

This device complies with RSS-210 of the Industry Canada Rules. Operation is subject to the following two conditions:

- *This device may not cause harmful interference.*
- *This device must accept any interference received, including interference that may cause undesired operation.*

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The user is cautioned that any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

The radio module is located in the skull and the IC ID is referenced on the manikin label.
IC ID 11306A-ISP091201

Canadian ICES-003 Statement

This Class B digital apparatus complies with Canadian ICES-003.

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Europe

This product is in compliance with the essential requirements of Council Directive 2014/53/EU on Radio Equipment (RED) and Council Directive 2011/65/EU on Restriction on the use of certain hazardous substance (RoHS).

Japan

MIC certification R 001-A06167

Waste Handling

This appliance is marked according to the European directive 2012/19/EC on Waste Electrical and Electronic Equipment (WEEE).

By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

The symbol on the product, or on the documents accompanying the product, indicates that this appliance may not be treated as household waste. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment. Disposal must be carried out in accordance with local environmental regulations for waste disposal.

For more detailed information about treatment, recovery and recycling of this product, please contact your local city office, your household waste disposal service or Laerdal representative.

Symbol Glossary

	CE mark
	Japanese MIC Certification
	Dispose of in accordance with your country's recommendations.
	Australian Radiocommunications and EMC Compliance Mark
	Manufacturer
	Date of Manufacture
	Reference Number
	Serial Number
	Warning/Caution

Premature Anne

Temperature and Environment	
Operating temperature	15 °C to 50 °C (59 °F to 122 °F)
Storage temperature	-15 °C to 50 °C (5 °F to 122 °F)
Operating humidity	10% to 95% relative humidity
Storage humidity	0% to 90% relative humidity
Dimensions	
Length - extended joints (manikin only)	Approximately 30 cm (12 in)
Weight (manikin only)	Approximately 435 g (0.9 lbs)
Volume of umbilical reservoir (with inserted umbilical)	Maximum 18 ml (3 ml simulated blood, 2 x 7.5 ml simulated volume expander)
RF Communication	
Bluetooth Low Energy (BLE)	2.4GHz
BLE Operation range	< 8m indoors (free sight)

Battery and Charging	
Battery type, volt- age and capacity	NiMH button cell battery, 3.6V, 600mAh
Battery run time (fresh battery)	Approx. 5 hours at 22°C (using cyanosis 1/3 of the time). The run time is HIGHLY dependent on the use of cyanosis)
Battery charging	USB micro B input, 5V, 250mA max
Charging time	Approx. 4 hours
Recommended Temperature Range for charging	15 °C to 35 °C (59 °F to 95 °F)

Premature Anne Task Trainer

Temperature and Environment	
Operating temperature	15 °C to 50 °C (59 °F to 122 °F)
Storage temperature	-15 °C to 50 °C (5 °F to 122 °F)
Operating humidity	10% to 95% relative humidity
Storage humidity	0% to 90% relative humidity
Dimensions	
Length - extended joints (manikin only)	Approximately 30 cm (12 in)
Weight (manikin only)	Approximately 435 g (0.9 lbs)
Volume of umbilical reservoir (with inserted umbilical)	Maximum 18 ml (3 ml simulated blood, 2 x 7.5 ml simulated volume expander)

Warranty

Premature Anne and Premature Anne Task Trainer have a one year limited warranty. Refer to the Laerdal Global Warranty for terms and conditions. For more information visit www.laerdal.com.

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- This device may not cause harmful interference, and
- This device must accept any interference received, including interference that may cause undesired operation.

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- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and the receiver
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected
- Consult the dealer or an experienced radio/TV technician for help.

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- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

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© 2019 Laerdal Medical AS. All rights reserved.