EN

FR

DE

ES

IT

BR

NL

PL

JA

ZH

KO

SimMan

Important Product Information

CAUTIONS AND WARNINGS

Disclaimer

Use of the SimMan Patient Simulator system to train personnel should be undertaken under supervision of suitably trained medical personnel with an understanding of educational principles as well as recognized medical protocols.

As with all Patient Simulators or other such training devices there may be approximations, variations and inaccuracies in anatomical features and the physiological modeling. This being the case, Laerdal Medical does not guarantee that all features are completely accurate.

General Simulator Handling

It is important to follow the instructions below, as well as other available User Information, in order to maintain optimum performance and longevity of the simulator components.

Cautions

- Do not provide artificial respiration to the Patient Simulator using oxygen enriched air or flammable gases.
- Do not introduce fluids into or onto the Patient Simulator (except as directed in the *User Guide*) as this may damage the Patient Simulator and its components.
- Do not introduce humidified air into the system during ventilation.
- Never perform mouth-to-mouth or mouth-to-nose rescue breathing on the patient simulator. The Patient Simulator's airways are not designed for cleaning or disinfection.
- Do not use the Patient Simulator if the internal tubing and cabling is disconnected.
- Do not store or use this product outside its specified operating and storage conditions.
- Do not use automated chest compression machines on the Patient Simulator.

Defibrillation Hazards

SimMan allows for defibrillation in accordance with 2015 international guidelines for CPR.

A conventional defibrillator may be used on the SimMan. During live defibrillation, the defibrillator and Patient Simulator may present a shock hazard. All standard safety precautions must be taken when using a defibrillator on the Patient Simulator. For more information, consult your defibrillator's User Guide.

The ECG connectors are designed exclusively for ECG monitoring and must not be used for defibrillation. Defibrillation on the ECG connectors will damage the internal electronics of the Patient Simulator and may cause personal injury.

Warnings

- Defibrillation must be performed on the defibrillator connectors only.
- The Patient Simulator must not come into contact with electrically conductive surfaces or objects during defibrillation.
- Do not defibrillate the patient simulator in a flammable or oxygen enriched atmosphere.
- The Patient Simulator torso must always be kept dry. Allow the Patient Simulator to acclimate before defibrillating. Sudden changes in temperature (moving the Patient Simulator from a cold environment to a warm environment and vice versa) may result in condensation collecting on the base board and pose a shock hazard.

Battery

Warnings

- The external battery charger is for indoor use only.
- The batteries should only be charged and used at the specified temperature ranges.
- Inserting and connecting batteries incorrectly, short circuiting or exposure to fluids pose an explosion hazard.
- Do not mistreat, disassemble or attempt to repair the battery.
- Do not use the batteries if they are visibly damaged, malfunction or appear to leak electrolyte.
- Take extreme care to avoid direct contact with electrolyte, hot or smoking parts. In case of the above, disconnect and remove the battery when it is judged safe to do so.

For other battery related information, consult User Guide.

REGULATORY INFORMATION

Federal Communications Commission (FCC) Statement

FCC-ID: QHQ-212-00001 (SimMan 3G only)

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1 This device may not cause harmful interference, and
- 2 This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna

- Increase the separation between the equipment and the receiver
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected
- Consult the dealer or an experienced radio/TV technician for help.

Canadian ICES-003 Statement

This Class B digital apparatus complies with Canadian ICES-003.

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Caution

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

CE EU

The product is in compliance with the essential requirements of Council Directive 1999/5/EC on radio and telecommunications terminal equipment (R&TTE). The product is in compliance with Council Directive 2011/65/EU on restriction of the use of certain hazardous substances (RoHS).

WASTE HANDLING

Dispose of in accordance with local requirements and regulations.

Li-ion

SimMan is operated on a Lithium Ion (Li-Ion) battery. Li-ion batteries should be recycled or disposed of in accordance with local regulations.

SPECIFICATIONS**Size and Weight**

Dimensions (Patient Simulator only)	1800 mm (l) x 550 mm (b) chest (5.90 ft x 1.80 ft)
Weight (Patient Simulator only)	38.5 kg (85 lbs)
Weight (with clothes)	40 kg (88 lbs)

Patient Simulator Power

External power	Input voltage 24VDC, 6.25A
Internal batteries (two)	Each 14.8V, 4.6Ah, Lithium-Ion
Only use the SimMan external power supply and batteries	

Air & CO₂ Pressure

Internal air tank	Max 0.85 bar
External air connection	Max 1.4 bar
External CO ₂ to Patient Simulator	Max 1.4 bar

Temperature Limits

Operating temperatures	+4 °C to 40 °C (39 °F to 104 °F)
Storage temperatures	-15 °C to 50 °C (5 °F to 122 °F)

Environment - Patient Simulator Only

Relative humidity	20% -90% (non-condensing)
DO NOT use outdoors in wet conditions	
Not tested with salt spray	

RF Communication

WLAN frequency ranges	- 2.4 GHz WLAN, channels 1-11 - 5 GHz WLAN, channels 36, 40, 44 and 48
WLAN operation range	100 m (300 ft) outdoors
RFID frequency ranges (SimMan 3G only)	13.56 MHz
RFID operation range (SimMan 3G only)	< 0.2 m

Material Chart for Patient Simulator

Clothes	Cotton, Nylon
Skins and airways	PVC (DEHP free)
External hard plastics	PP, PA, PC, PC/PET
Inner plastics	Silicone, TPU, TPE, PVC, Nitrile PA, PA+GF, PC, ABS, POM, HDPE, PET
Metal components	Aluminum, Brass, Steel

Recommended Hardware Specifications

Intel i-core generation 3 or newer
Score over 3000 PassMark - CPU Mark
4 GB RAM
120 GB hard disk space
1366x768 or better
Headset w/microphone
USB ports – 2

Minimum Software Requirements

Windows 7 or Windows 8
100% DPI

Acceptable Fluids for Patient Simulator

To clean the simulator's fluid and blood system use one of the following:	
Distilled water or De-ionized water	
Cleaning Fluids	- 60-70% Isopropanol

Simulated IV Fluids

Only use purified water to simulate IV fluids:	
Simulated IV fluids	- Distilled Water - Deionized Water

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© 2016 Laerdal Medical AS. All rights reserved.
SimMan® and Laerdal® are registered trademarks of Laerdal Medical AS.