EN

FR

DE

ES

IT

BR

NL

PL

JA

ZH

KO

Nursing Anne Simulator & Nursing Anne Simulator Geriatric

Important Product Information
重要产品信息

www.laerdal.com

Laerdal

helping save lives

English

Nursing Anne Simulator/ Nursing Anne Simulator Geriatric

Read these instructions thoroughly. Observe all warnings, precautions and instructions on the product, in the User Guide and in this Important Product Information booklet. Retain this booklet for future reference.

Warnings and Cautions

A Warning states a condition, hazard, or unsafe practice that can result in serious personal injury or death.

A Caution states a condition, hazard, or unsafe practice that can result in minor personal injury or damage to the product.

Note

A note states important information about the product or its operation.

Disclaimer

Use of the Nursing Anne Simulator system to train personnel should be undertaken under supervision of suitably trained technical or medical personnel with an understanding of educational principles as well as recognized medical protocols.

As with all patient simulators or other such training devices, there may be approximations, variations and inaccuracies in anatomical features and physiological modeling. Laerdal does not guarantee that all features are completely accurate.

General Simulator Handling

Follow all instructions available in the User Guide, Quick Setup Guide and this document, in order to maintain optimum performance and longevity of the simulator components.

General Operation

Cautions

- Avoid pinch hazards – Do not use the patient simulator without the external skins.
- Avoid contact with sharp objects or surfaces and excessive pulling on skin.
- Do not leave the patient simulator unattended when sitting up unassisted.

This poses a fall hazard, risking damage to simulator or equipment in area and injury to personnel.

- Do not introduce fluids into or onto the patient simulator (except as directed in the User Guide) as this may damage it and its components.
- Do not use the patient simulator if the internal tubing and cabling is disconnected.
- Do not use the patient simulator in aircrafts when they are airborne.
- Do not use outdoors in wet conditions. Not tested with salt spray.
- Do not perform mouth-to-mouth resuscitation.
- Use minimum two people when lifting the manikin.
- Do not perform chest compressions while spontaneous breathing is enabled as this can damage internal parts.
- Do not use mineral oil, MEK, or acetone as a cleaning agents.
- Avoid using alcohol on the hard, polycarbonate underbody surfaces near mounting screws.

Patient Simulator Power - Battery

Warnings

- Only use the designated Laerdal approved battery in Nursing Anne Simulator.
- Do not use the Laerdal Li-Ion Battery for other purposes than specified.
- Use only the AC charger adapter supplied with the product.
- Do not use the battery in temperatures above what is specified in this Important Product Information.
- Do not heat or incinerate. Do not crush the battery.
- Do not short circuit the battery contact.
- Do not immerse in water.
- The external battery charger is for indoor use only.
- Inserting and connecting battery incorrectly, short circuiting or exposure to fluids pose an explosion hazard.

Nursing Anne Simulator/ Nursing Anne Simulator Geriatric

- Do not mistreat, disassemble or attempt to repair the battery.
- Do not use the battery if it is visibly damaged, malfunctions or appears to leak electrolyte.
- Take extreme care to avoid direct contact with electrolyte, hot or smoking parts. In case of the above, disconnect and remove the battery when it is judged safe to do so.
- Do not apply conductive gel or conductive defibrillation pads intended for patient use to prevent chest skin pitting.
- Do not use cables or connectors with visible damage.

ESD

Warning

Take care when cleaning the silicone skins with flammable chemicals as sparks or fire may occur.

Oxygen


Warning

Take care when using oxygen together with the manikin. There is a risk of fire.

Patient Simulator - External Power

Warnings

- Only connect the Patient Simulator to a Laerdal approved external power supply dedicated to use on Nursing Anne Simulator.
- Do not use the Patient Simulator in close proximity to real patients when it is connected to the external power supply.
- Patient Simulator should be connected to an electrical outlet that is properly grounded.

Air Tank

Warnings

- Do not disassemble or attempt to repair the internal parts of the right thigh.
- Do not puncture.

Defibrillation

Warnings

- Observe normal safety precautions for use of defibrillators.
- Defibrillation must be performed using ShockLink only. Refer to ShockLink Important Product Information.
- Do not attach medical device defibrillator electrodes or place defibrillator paddles on any part of the simulator.
- Only defibrillate the patient using ShockLink when it is fully assembled, dry and undamaged.
- Do not defibrillate the patient simulator when it is turned off.
- Do not defibrillate the Patient Simulator in a flammable or oxygen-enriched atmosphere. This is a fire-hazard.

CAN CPU Module (Component of simulator)

Cautions

- The device for operation in the band 5150-5250 MHz is only for indoor use to reduce the potential for harmful interference to co-channel mobile satellite systems;
- The maximum antenna gain permitted for devices in the bands 5250-5350 MHz and 5470-5725 MHz shall comply with the e.i.r.p. limit; and
- The maximum antenna gain permitted for devices in the band 5725-5825 MHz shall comply with the e.i.r.p. limits specified for point-to-point and non point-to-point operation as appropriate.
- The worst-case tilt angle(s) necessary to remain compliant with the e.i.r.p. elevation mask requirement set forth in Section 6.2.2(3) shall be clearly indicated.
- Users should also be advised that high-power radars are allocated as primary users (i.e. priority users) of the bands 5250-5350 MHz and 5650-5850 MHz and that these radars could cause interference and/or damage to LE-LAN devices.

USA and Canada

FCC

Federal Communications Commission Statement and Industry Canada Statements

This device complies with part 15 of the FCC rules and RSS-210 of the Industry Canada rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. this device must accept any interference received, including interference that may cause undesired operation.

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This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.

- Consult the dealer or an experienced radio/TV technician for help.

Contains FCC ID: QHQ-2011480

Contains IC ID: 2063-2011480

The term "IC" before the equipment certification number only signifies that the Industry Canada technical specifications were met.

Caution: Changes or modifications not expressly approved by Laerdal Medical could void the user's authority to operate the equipment.

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Radiation Exposure Statement

This equipment complies with IC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance 20 cm between the radiator & your body.

EU

This product is in compliance with the essential requirements of Council Directive 2014/53/EU on Radio Equipment (RED) and Council Directive 2011/65/EU on restriction of the use of certain hazardous substances (RoHS).

Korea

R-R-Lm1-320-05050

China

CMIIT ID: 2019AJ3157

Japan

MIC ID: R201-171040

Symbol Glossary

	CE Mark
	Waste Electrical and Electronic Equipment symbol
	Li-Ion batteries recycling symbol
	Australian Regulatory Compliance Mark
	Korea Certification (KC) Mark
	MIC Technical Conformity Mark (Japan)
	Manufacturer
	Date of Manufacture
	Reference Number
	Serial Number

Waste Handling

Dispose of in accordance with your country's recommendations.

This appliance is marked according to the European directive 2012/19/EC on Waste Electrical and Electronic Equipment (WEEE).

By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

The symbol on the product, or on the documents accompanying the product, indicates that this appliance may not be treated as household waste. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment. Disposal must be carried out in accordance with local environmental regulations for waste disposal.

For more detailed information about treatment, recovery and recycling of this product, please contact your local city office, your household waste disposal service or Laerdal representative.

Battery

Nursing Anne Simulator is operated on a Lithium Ion (Li-Ion battery). Li-Ion batteries should be recycled or disposed of in accordance with local regulations.

Warranty

Refer to the Laerdal Global Warranty for terms and conditions. For more information visit www.laerdal.com

Size and Weight	
Dimensions (Patient Simulator only)	170 cm (5' 7")
Weight (Patient Simulator only)	31.75 kg (70 lbs)
Patient Simulator Power	
External Power	12VDC, 3.3A
Internal Battery	14.52 V, 5.2 Ah, 76Wh
Battery Time	6 hours (default settings)
Temperature Limits	
Operating Temperature	4 °C to 40 °C (39 °F to 104 °F)
Storage Temperature	-15 °C to 50 °C (5 °F to 122 °F)
Battery Charging	19 °C to 25 °C (66 °F to 77 °F)
Environment - Patient Simulator Only	
Relative Humidity	20% - 90% (non-condensing)
RF Communication	
WLAN Frequency Ranges	2.4 GHz WLAN 5 GHz WLAN
WLAN Operation Range	100 m (300 ft) outdoors
Pressure Data (Compressor)	
Working Pressure	1.5 bar
Working Temperature	0 °C to 40 °C (32 °F to 104 °F)
Capacity of the vessel (V in l)	0.3 l
Blood Pressure	
Pressure Range	0-50 mmHg
Accuracy	+/-4 mmHg
Chest Compressions	
Maximum Compression Depth	65 mm

Material Chart for Patient Simulator	
Clothes	Cotton, Polyester
Skins and airways	Silicone, PVC (DEHP free)
External hard plastics	PP, PA, PC, PC/PET, PVC, ABS
Inner plastics	Silicone, TPU, TPE, PVC, Nitrile PA, PA+GF, PC, ABS, POM, HDPE, PET
Metal components	Aluminium, Brass, Steel

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USA and Canada

FCC

Federal Communications Commission
Statement and Industry Canada Statements

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This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

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© 2021 Laerdal Medical AS. All rights reserved.