EN

FR

DE

ES

BR

JA

ZH

KO

Laerdal-SonoSim

Ultrasound Solution

&

Procedure Trainer

Important Product Information

www.laerdal.com

Laerdal

helping save lives

This Important Product Information covers Laerdal-SonoSim Ultrasound Solutions and Laerdal-SonoSim Procedure Trainer. Unless otherwise specified the information applies to all product configurations.

Read these instructions thoroughly. Observe all warnings, cautions and instructions in the User Guide and in this Important Product Information booklet. Retain this booklet for future reference.

Warnings and Cautions

A Warning states a condition, hazard, or unsafe practice that can result in serious personal injury or death.

A Caution states a condition, hazard, or unsafe practice that can result in minor personal injury or damage to the product.

Notes

A note states important information about the product or its operation.

Laerdal-SonoSim Ultrasound Solution

LS Tags

Warnings

LS Tags must not be used on anyone using a pacemaker, implantable cardioverter defibrillator, or other electronic medical device. Components of the technology may interfere with such medical devices.

Cautions

- *The LS Tags located in the chest, gravid, and non-gravid skins should not be used during manual or automatic birthing simulations.*
- *Do not use ultrasound gel.*
- *Avoid exposing the LS Tags and LS Probe to any liquids.*
- *Do not cut or puncture the LS chest, gravid, and non-gravid skins.*

Laerdal-SonoSim Procedure Trainer

Needle & Syringe

Caution

The tip of the Needle & Syringe is sharp. Take care when handling.

Warranty

Refer to the Laerdal Global Warranty for terms and conditions. For more information visit www.laerdal.com.

Federal Communications Commission Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Changes or modifications not expressly approved by Laerdal Medical could void the user's authority to operate the equipment.

Industry Canada Statement

This device complies with RSS-210 of the Industry Canada Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Caution: Changes or modifications not expressly approved by Laerdal Medical could void the user's authority to operate the equipment.

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The term "IC" before the equipment certification number only signifies that the Industry Canada technical specifications were met.

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USA and Canada

Laerdal-SonoSim Ultrasound Solution

LS Probe contains FCC IC: 2AEMESSLS
 and IC Certificate: 20197-SSLS

Laerdal-SonoSim Procedure Trainer

Needle & Syringe contains
 FCC ID: 2AEME-1002040 and
 IC Certificate: 20197-1002040

EU

CE:

This product is in compliance with the essential requirements of Council Directive, 2014/53/EU on Radio Equipment (RED) and Council Directive 2011/65/EU on restriction of the use of certain hazardous substances (RoHS).

Waste Handling

 Dispose of in accordance with your country's recommendations

This appliance is marked according to the European directive 2012/19/EC on Waste Electrical and Electronic Equipment (WEEE). By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

The symbol on the product, or on the documents accompanying the product, indicates that this appliance may not be treated as household waste. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment. Disposal must be carried out in accordance with local environmental regulations for waste disposal.

For more detailed information about treatment, recovery and recycling of this product, please contact your local city office, your household waste disposal service or Laerdal representative.

Symbol Glossary

	CE mark
	WEEE symbol
	Australian Radiocommunications and EMC Compliance Mark
	Warning/Caution
	Note
	Manufacturer
	Reference Number
	Serial Number

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© 2019 Laerdal Medical AS. All rights reserved.