Stifneck

User Guide

98010033 Stifneck Baby No-Neck (Qty. 1)
98020033 Stifneck Pediatric (Qty. 1)
98030033 Stifneck No-Neck (Qty. 1)
98040033 Stifneck Short (Qty. 1)
98050033 Stifneck Regular (Qty. 1)
98060033 Stifneck Tall (Qty. 1)

Manufactured in China for:
LAERDAL MEDICAL AS, Norway
P.O. Box 377, N-4002 Stavanger, Norway
Tel. +47 51 51 17 00,
Fax +47 51 52 35 57

ENGLISH
ČESKÝ
DANSK
DEUTSCH
ESPAÑOL
FRANÇAIS
ITALIANO
NEDERLANDS
NORSK
POLSKI
PORTUGUÊS
Русский
SVENSK
SUOMI
中文
한국어

Laerdal

helping save lives

Contents

ENGLISH	4
ČESKÝ	5
DANSK	6
DEUTSCH	7
ESPAÑOL	8
FRANÇAIS	9
ITALIANO	10
NEDERLANDS	11
NORSK	12
POLSKI	13
PORTUGUÊS	14
Русский	15
SVENSK	16
SUOMI	17
中文	19
한국어	20
REGULATORY INFORMATION	21

Intended Use

This product is intended to support cervical spine in a neutral position during transportation, in combination with other cervical and full body immobilization devices (CSIDs and FBSIDs).

Important


- In cases of a suspected spinal cord injury, proper neck immobilization is only one part of total immobilization.
- It is imperative that the patient be properly immobilized to prevent any movement of the spinal column.
- Should any serious malfunction, undesirable incident with, or deterioration in the functionality or performance of the device occur, contact Laerdal promptly. The competent authority where the incident took place and/or the device was used should also be notified.

Warnings

- Stifneck Collars should only be used by qualified healthcare providers.
- Improper use may result in injury or permanent disability. Always follow local spinal immobilization protocol.
- Use the correct size collar. An oversized collar may hyperextend a patient's cervical spine. An undersized collar may not provide appropriate stability.

User Instructions


1. Measure the Patient


Align the head to neutral or "eyes forward" position unless contraindicated by local protocol.

2. Match the collar size to the patient:


Plastic edge to hole or to rivet if assembled.


3. Insert the rivet into the hole.

4. Preform the collar


5. Apply the collar while manually maintaining neutral head position

Place the chin support under the chin.

Pull the back of the collar snug while holding the front in place, and then fasten.

For a supine patient, slide the rear panel behind the neck before placing the chin support.

Technical Specifications

Storage:
Do not store collar in folded position.
Store in flat position.

Storage Temperature Range:
-34 °C to 52 °C
(-30 °F to 125 °F)

Operating Temperature Range:
-18 °C to 43 °C
(0°F to 110 °F)

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