EN

FR

DE

ES

IT

BR

NL

PL

JA

ZH

KO

SimNewB

Important Product Information
重要产品信息

www.laerdal.com

Laerdal

helping save lives

This Important Product Information covers SimNewB. Unless otherwise specified the information applies to all product configurations.

Read these instructions thoroughly. Observe all warnings, cautions and instructions in the User Guide and in this Important Product Information booklet. Retain this booklet for future reference.

Warnings and Cautions

A Warning states a condition, hazard, or unsafe practice that can result in serious personal injury or death.

A Caution states a condition, hazard, or unsafe practice that can result in minor personal injury or damage to the product.

Notes

A note states important information about the product or its operation.

General Simulator Handling

Warnings

- Do not use medical devices during simulation with SimNewB unless the medical devices have been subjected to the local cleaning and disinfection procedures both prior to and immediately after the simulation (risk of cross-contamination and infections).
- Do not use SimNewB closer than 2.5 m vicinity to a patient.
- Do not provide artificial respiration to the Patient Simulator using oxygen enriched air or flammable gases.
- Do not use or charge the simulator in the presence of oxygen sources over 21%, or in the presence of flammable gases.
- Do not carry out chest compression during spontaneous breathing of the simulator.
- Do not disconnect the legs at the hip.
- Do not use real umbilical cords during simulation with SimNewB.

Cautions

- Do not lift SimNewB by the chest.
- Do not perform peripheral IV procedures
- Do not apply counterpressure to the arms during use of automated movements.
- Do not introduce fluids into or onto the patient simulator (except as directed in the User Guide) as this may damage it and its components.
- Do not store or use this product outside its specified operating and storage conditions.
- Do not use the Patient Simulator in aircrafts when they are airborne.
- Do not store the simulator with fluid in the Umbilical reservoir and IO system. Use a syringe to remove remains of any injected water from the tubing and components prior to storage.

Defibrillation

Warning

Do not perform defibrillation and pacing.

Airway Management

Cautions

- Never perform mouth-to-mouth or mouth-to-nose ventilation on the Patient Simulator. The airways are not designed for cleaning or disinfection.
- Do not put biological materials or other fluids into the simulator's airways or stomach. They are not designed for cleaning or disinfection.
- Do not introduce humidified air into the system during ventilation.
- Use ONLY silicone lubricant provided by Laerdal Medical. (Cat. No. 410-00050)
- Do NOT use Laerdal Training Lubricant (Cat.No. 250-21050/ 252090)
- Using other lubricant can damage the airway and/or simulator.
- Lubricate laryngoscope blade and tubes prior to insertion into the airway. Non-lubricated instruments and tubes may cause damage to the airway.

Patient Simulator Power - Battery

Warnings

- Use only power supply that comes with the product.
- Do not charge the simulator when it is placed on a heating mattress and/or under a heater.
- Do not charge the simulator when it is covered in blankets.
- Do not use the battery in temperatures above what is specified in this Important Product Information.
- Do not heat or incinerate. Do not crush the battery.
- Do not short circuit the battery contact.
- Do not immerse in water.

Notes

- In case of high battery temperature during simultaneous operation and charging, the charging would pause automatically to protect the battery. Charging will automatically resume when the battery temperature has decreased.
- The Patient Simulator contains electronic components. All components, including the battery must be recycled and disposed of in compliance with applicable laws and regulations.

For other battery related information, consult User Guide.

Patient Simulator - External Power

Warning

Only connect the Patient Simulator to a Laerdal approved external power supply dedicated to use on SimNewB.

CPU Module (component of simulator)

Cautions

- The device for operation in the band 5150-5250 MHz is only for indoor use to reduce the potential for harmful interference to co-channel mobile satellite systems;
- The maximum antenna gain permitted for devices in the bands 5250-5350 MHz and 5470-5725 MHz shall comply with the e.i.r.p. limit; and
- The maximum antenna gain permitted for devices in the band 5725-5825 MHz shall comply with the e.i.r.p. limits specified for point-to-point and non point-to-point operation as appropriate.
- The worst-case tilt angle(s) necessary to remain compliant with the e.i.r.p. elevation mask requirement set forth in Section 6.2.2(3) shall be clearly indicated.
- Users should also be advised that high-power radars are allocated as primary users (i.e. priority users) of the bands 5250-5350 MHz and 5650-5850 MHz and that these radars could cause interference and/or damage to LE-LAN devices.

FCC

Federal Communications Commission Statement and Industry Canada State- ments

This device complies with part 15 of the FCC rules and RSS-210 of the Industry Canada rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. this device must accept any interference received, including interference that may cause undesired operation.

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SimNewB contains IC: 20263-2011480

Radiation Exposure Statement

This equipment complies with IC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance 20 cm between the radiator & your body.

Japan

MIC certification R201-171040

Korea

R-R-Lm1-296

China

CMIIT ID: 2020AJ8314

EU

CE: This product is in compliance with the essential requirements of Council Directive 2014/53/EU on Radio Equipment (RED) and Council Directive 2011/65/EU on restriction of the use of certain hazardous substances (RoHS).

Waste Handling

Dispose of in accordance with your country's recommendations.

This appliance is marked according to the European directive 2012/19/EC on Waste Electrical and Electronic Equipment (WEEE).

By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

The symbol on the product, or on the documents accompanying the product, indicates that this appliance may not be treated as household waste. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment. Disposal must be carried out in accordance with local environmental regulations for waste disposal.

For more detailed information about treatment, recovery and recycling of this product, please contact your local city office, your household waste disposal service or Laerdal representative.

Symbol Glossary

Symbol	Definition
	CE mark
	WEEE symbol
	Australia Radiocommunications and EMC Compliance Mark
	MIC Technical Conformity Mark (Japan)
	Korean Certification (KC) Mark
	Warning / Caution
	Note
	Manufacturer
	Date of Manufacture
	Reference Number
	Serial Number

Warranty

Refer to the Laerdal Global Warranty for terms and conditions. For more information visit www.laerdal.com.

Size and Weight	
Length (patient simulator only)	52 cm (20,4 in)
Weight (patient simulator only)	2.3 kg (4.4 lbs)
Patient Simulator Power	
External Power	Input voltage 12VDC, 6,25 A, 75W
Internal Battery (one)	12V, 2000mAh, 24Wh, NiMH
Battery Time	≈ 5 hours
Temperature Limits	
Operating temperatures	5 °C to 40 °C (39 °F to 104 °F)
Storage temperatures	-15 °C to 50 °C (5 °F to 122 °F)
Battery Charging temperatures	5 °C to 35 °C (39 °F to 95 °F)
Environment - Patient Simulator Only	
Relative humidity	10% - 90% (non-condensing)
DO NOT use outdoors in wet conditions	
Not tested with salt spray	
RF Communication	
WLAN Frequency Bluetooth Frequency	2.4 GHz and 5 GHz 2.4 GHz
Operation range	10 m (30 ft) max.
Material Chart for Patient Simulator	
Clothes	Cotton, Nylon
Skins and airways	Silicone
External hard plastics	PVC, ABS
Inner plastics	Silicone, PVC, ABS, Nylon, Polyurethane
Metal components	Aluminum, Brass, Stainless steel, copper
Cleaning Fluids	
Use mild solution of liquid soap and water to clean the Simulator	

Compatible Adjunct Sizes	
ET Tube	3.5mm (uncuffed) - to mark 9
LMA	Size 1
Laryngoscope	Size 0 or 1 Miller Blade
OG / NG tube	8 Fr
Suction catheter	10-14 Fr
UVC / UAC catheter	5 Fr / 3.5 Fr
IO	≤ 14 Gauge
Needle decompression (pneumothorax)	18 – 22 Gauge
Chest drain	10 – 12 Fr

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FCC**Federal Communications Commission Statement and Industry Canada Statements**

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This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.

- Increase the separation between the equipment and receiver.
- Consult the dealer or an experienced radio/TV technician for help.

The term "IC" before the equipment certification number only signifies that the Industry Canada technical specifications were met.

Caution: Changes or modifications not expressly approved by Laerdal Medical could void the user's authority to operate the equipment.

SimNewB contains FCC ID: QHQ-2011480

Canada**ICES-003 Statement**

This Class B digital apparatus complies with Canadian ICES-003.

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© 2022 Laerdal Medical AS. All rights reserved.