# Resusci Anne QCPR RQI

Important Product Information  
重要产品信息

EN

FR

DE

ES

IT

BR

NL

NO

SV

DK

FI

JA

ZH

KO

RU

PL

[www.laerdal.com](http://www.laerdal.com)

**Laerdal**  
helping save lives

Read these instructions thoroughly. Observe all warnings, precautions and instructions on the product, in the User Guide and in this Important Product Information booklet. Retain this booklet for future reference.

 **Warnings and Cautions**

A Warning states a condition, hazard, or unsafe practice that can result in serious personal injury or death.

A Caution states a condition, hazard, or unsafe practice that can result in minor personal injury or damage to the product.

 **Note**

A note states important information about the product or its operation.

**General**

 **Caution**

*Use of automatic chest compression machines may damage the manikin.*

**Resusci Anne QCPR AED**

 **Warnings**

- *Observe all standard safety precautions for the use of defibrillators.*
- *Do not allow the manikin to contact electrically conductive surfaces or objects during defibrillation.*
- *Do not defibrillate in a flammable or oxygen-enriched atmosphere.*
- *Only defibrillate using ShockLink as described in ShockLink instructions.*

 **Cautions**

- *When removing or replacing the chest skin, do not pull or damage the wires connecting the chest skin to the battery box.*
- *To prevent skin pitting on the Resusci Anne QCPR AED manikin, do not apply conductive gel or conductive defibrillation pads intended for patient use.*

 **Note**

*Use training pads with a white connector (Cat. No. 94 50 90).*

### **Federal Communications Commission (FCC) and Industry Canada (IC) Statements**

This device complies with part 15 of the FCC Rules and RSS-210 of IC rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Consult the dealer or an experienced radio/TV technician for help.

**Caution:** Changes or modifications not expressly approved by Laerdal Medical could void the user's authority to operate the equipment.

The term "IC" before the equipment certification number only signifies that the Industry Canada technical specifications were met.

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This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

1. This device may not cause interference.
2. This device must accept any interference, including interference that may cause undesired operation of the device.

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Resusci Anne QCPR contains FCC ID:  
QH20-10494

Contains IC Certification Number:  
IC 20263 20-10494

### **Canada**

This Class B digital apparatus complies with Canadian ICES-003.

## **Japan**

MIC certification 012-180007

## **Korea**

R-CMM-Lm1-173-00160

## **China**

Basic Torso: CMIIT ID: 2018DJ4593

AED Torso: CMIIT ID: 2018DJ4596

## **EU**

CE: This product is in compliance with the essential requirements of Council Directive 2014/53/EU on Radio Equipment (RED), and Council Directive 2011/65/EU on restriction of the use of certain hazardous substances (RoHS).

## **UK**

UKCA (U.K. Conformity Assessed): This product complies with the requirements of U.K. legislation of Radio Equipment Regulations (S.I. 2017/1206).

UKCA (U.K. Conformity Assessed): This product complies with the requirements of U.K. legislation of Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations (RoHS) Regulations (S.I. 2012/3032).

## **Waste Handling**

Dispose of in accordance with your country's recommendations.

This appliance is marked according to the European directive 2012/19/EC on Waste Electrical and Electronic Equipment (WEEE).

By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

The symbol on the product, or on the documents accompanying the product, indicates that this appliance may not be treated as household waste. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment. Disposal must be carried out in accordance with local environmental regulations for waste disposal.

For more detailed information about treatment, recovery and recycling of this product, please contact your local city office, your household waste disposal service or Laerdal representative.

## Symbol Glossary

|                                                                                   |                                                        |
|-----------------------------------------------------------------------------------|--------------------------------------------------------|
|   | CE Mark                                                |
|   | UKCA Mark                                              |
|   | Australian Radiocommunications and EMC Compliance Mark |
|   | MIC Technical Conformity Mark (Japan)                  |
|   | Manufacturer                                           |
|   | Date of Manufacture                                    |
|   | WEEE Symbol                                            |
|   | Reference number                                       |
|   | Serial Number                                          |
|  | Warning / Caution symbol                               |

## Specifications

| Resusci Anne QCPR, QCPR AED |                                                  |
|-----------------------------|--------------------------------------------------|
| Dimensions                  | 157 cm x 52 cm x 25 cm<br>(61.8" x 20.5" x 9.8") |
| Weight                      | 79.2 lbs (36 kg)                                 |
| Operating temperature       | 0 °C to +40 °C<br>(32 °F to 104 °F)              |
| Humidity                    | < 95% relative humidity                          |
| Storage temperature         | -15 °C to +50 °C<br>(5 °F to +122 °F)            |
| Electronics                 |                                                  |
| Maximum Output Power        | -2.7dBm                                          |
| Frequency Range             | 2402 MHz to 2480 MHz                             |

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### **Federal Communications Commission (FCC) and Industry Canada (IC) Statements**

This device complies with part 15 of the FCC Rules and RSS-210 of IC rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Changes or modifications not expressly approved by Laerdal Medical could void the user's authority to operate the equipment.

The term "IC" before the equipment certification number only signifies that the Industry Canada technical specifications were met.

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© 2021 Laerdal Medical AS. All rights reserved.