EN

Laerdal BaXstrap Spineboard - User Guide

[REF] 982500, 982599, 982600 and 982699

Brand Name: Laerdal BaXstrap Spineboard
Common Name: Backboard/Spineboard/Longboard
Description: Full-body length spinal immobilization device

Intended Use: A base structure to be used with other adjunct cervical spine and head immobilization devices to facilitate in-line, neutral immobilization and transport of adult and pediatric patients.

Indications for Use: For use with patients having known or suspected spine or spinal cord injuries when patient immobilization and transportation is indicated.

Contraindications: None known

- Notes**
- Always follow local protocol for patient movement and immobilization techniques.*
- Not made with natural rubber latex.*
- The product is X-ray, CT, and MRI compatible.*
- Should any serious malfunction, undesirable incident with, or deterioration in the functionality or performance of the device occur, contact Laerdal promptly. The competent authority where the incident took place and/or the device was used should also be notified.*

Warnings and Cautions

Read and thoroughly understand these directions. These directions are not intended as, and do not in any way constitute, medical treatment advice, instruction in or teaching of medical practice, or imply or convey any rights, medical direction or authority. This and other spinal immobilization adjunct devices should only be used by personnel who have received adequate training. Inspect the BaXstrap prior to initial use and before each patient use. Pay particular attention to the underside of the head end. If the product develops punctures, broken pins or cracks then discontinuse use immediately. Ensure that the patient is secured using the straps. The length of the non-adjustable strap can be shortened by wrapping it around the handle of the BaXstrap.

Directions for Use

- Immobilize the cervical spine using a proper Cervical Immobilization Device
- Position the patient on the BaXstrap using appropriate number of personnel
- Properly fasten/strap the patient to the BaXstrap to protect against patient movement in any direction (up, down, left or right).
- Secure the head to the BaXstrap using a proper Head Immobilization Device.
- Assure in-line neutral immobilization of the head, neck, torso and pelvis.

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Cleaning: Non-sterile, non-critical infection risk. Withstands cleaning in soap and water (1:16) and 1:100 household bleach: water.

Technical Specifications	
Operating and storage temperature	-34 °C to +52 °C (-29 °F to +126 °F)
Operating and storage humidity	5-95% RH non-condensing
Shelf-life	3 years
Expected service-life	5 years <p>300 cycles of use Note: This is based on an average patient size of 150 kg and the defined reprocessing method.</p>
Maximum load	200 kg
Symbol Glossary	
	This medical device complies with the general safety and performance requirements of Regulation (EU) 2017/745 for medical devices.
	Medical Device

The product has a 5 year limited warranty. See the enclosed "Laerdal Global Warranty" for terms and conditions. The warranty is also available at www.laerdal.com.

CE

MD
Medical Device

The product has a 5 year limited warranty. See the enclosed "Laerdal Global Warranty" for terms and conditions. The warranty is also available at www.laerdal.com.

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CE

The product is accompanied by a limited 5 year warranty. See the enclosed "Laerdal Global Warranty" for terms and conditions. The warranty is also available at www.laerdal.com.

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© 2023 Laerdal Medical AS. All rights reserved.
BaXstrap is a trademark of Laerdal Medical AS.
Formed Plastics Inc. 207 Stonehenge Lane Carle Place, NY 11514, USA