
- EN
- FR
- ES
- IT
- BR
- NL
- PL
- JA
- ZH
- KO

SimMom

Guidelines on Best Practices on PPH Functionality

GUIDELINES ON BEST PRACTICES ON PPH FUNCTIONALITY

Introduction

During Postpartum Hemorrhage (PPH) simulation scenarios, some parts of the simulator inside the pelvis can be damaged by exposition to large amounts of fluids. To avoid this, Laerdal Medical AS recommends:

- The installation of a new power hatch plug, to use during every PPH simulation.
- Guidelines on how to set up the simulator before starting PPH simulation.

Installing the power hatch plug

This new accessory is intended to stop fluid from entering the SimMom pelvis during Postpartum Hemorrhage simulation scenarios.

1. Open the power hatch located on the side of the pelvis.

2. Insert the power hatch plug in the power hatch.

3. Press to keep it in place.

4. Close the power hatch.

5. Check that the power hatch can close properly with the power hatch plug inside.

 Note: After every PPH simulation, sponge out the blood gathered inside the power hatch.

Setting up the simulator for PPH simulation

Caution:

- Do not open the service panel located on the back side of the pelvis. The service panel should only be removed and installed by a qualified Laerdal Service Technician.
- Do not perform a PPH simulation when the simulator's head is down. Make sure the head and upper body of the simulator are raised during PPH simulation to ensure the correct flow of blood through the simulator.

To avoid excessive pooling of fluids in the bed, spread an underpad under the simulator's pelvis.

 Note: Sponge out the excess of blood accumulated in the birthing canal before turning the simulator on its side

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