EN

ES

NeoBeat

Newborn Heart Rate Meter

User Guide - US version

NeoBeat 532-01001

NeoBeat Mini 531-01001

www.laerdal.com

Laerdal

helping save lives

Clinical Indications	4
Clinical Indications	5
Important Information	6
Overview	8
Charging	10
Clinical Use	12
Clinical Use – Display	14
Reprocessing	18
Installation	21
Service and Warranty	24
Troubleshooting	25
Specifications	28
Technical Description	31

Clinical Indications

Device Description

NeoBeat is a heart rate measurement device using dry-electrode technology, with integrated digital display.

Indication for Use

NeoBeat and NeoBeat Mini are indicated to continuously measure and display the heart rate of neonates using dry electrodes on the torso during transition, stabilization and/or resuscitation. The devices are intended to be used in healthcare facilities.

NeoBeat is intended for use on newborns approximately 1.5 – 5 kg.

NeoBeat Mini is intended for use on newborns approximately 0.5 – 2 kg.

Intended Users

NeoBeat is intended to be used by healthcare professionals involved in newborn care.

Clinical Benefits

The intended clinical benefit of NeoBeat is the positive impact on newborn assessment by providing rapid and accurate presentation of the newborn's heart rate in real time.

Clinical Outcome

Desired outcome of NeoBeat is the presentation of the heart rate in real time.

Known Side Effects

None known.

Contraindications

None known.

Limitation of the Device

NeoBeat does not provide any type of alarms and has no monitoring capabilities. Clinical personnel must be present during the use of NeoBeat at all times. NeoBeat cannot be used unattended.

Important Information

The information in this User Guide applies to both NeoBeat and NeoBeat Mini.

Prior to first use, read the User Guide completely to become familiar with the operation and maintenance of NeoBeat. Read all Cautions and Warnings before using NeoBeat.

Warnings and Cautions

A Warning states a condition, hazard, or unsafe practice that can result in serious personal injury or death.

A Caution states a condition, hazard, or unsafe practice that can result in minor personal injury or damage to the product.

Notes

A Note states important information about the product or its operation.

Warnings

- *Decisions on when to start or end resuscitation efforts should not be made based on the output of this device alone.*
- *If an error has occurred, if the device provides no heart rate or if you do not trust the output, continue therapy without device.*
- *Handling of the newborn may cause false heart rate readings even when the patient has no heart rate (asystole). The device alone should not be used to confirm asystole/stillbirth.*
- *Care should be taken with newborns with exposed internal organs. Ensure placement on area of torso with intact skin.*
- *To avoid potential skin damage and infection, take care when placing NeoBeat Mini on preterms.*
- *Heart rate may not be detected correctly in the case of severe arrhythmia.*

Important Information

Cautions

- Excessive patient handling and movement may cause missing or erroneous heart rate readings.
- Do not misuse the device; e.g. using it on adults or children, or exposing the device to hard surface impacts.
- Do not use the product if it is damaged or cracked.
- NeoBeat is not intended to be used together with a defibrillator or HF surgical equipment. Remove NeoBeat if intending to use such devices on the newborn.
- Not intended to be sterilized as it may damage the device

Notes

- The device detects and displays the electrical heart rate, which in some conditions (e.g. pulseless electrical activity), may not reflect presence of circulatory pulse.
- Another person touching the electrodes simultaneously with the patient may interfere with the heart rate measurement.
- The device may indicate patient contact if one of the charging pins is touched while holding one ECG electrode.
- If the device is placed on a conductive surface, e.g. a metal tray, it may not go to standby, thus draining the battery.
- The electrodes are intended for use on moist newborn skin, and may reduce the heart rate accuracy when used with the thicker and drier skin of e.g. an adult's fingers.
- NeoBeat may not display a heartrate when used in close proximity to a Wireless Power Transfer charger.
- NeoBeat does not store, display, or output an ECG signal.


Before Use

Clean and disinfect NeoBeat as described in Reprocessing.

Charge NeoBeat as described in Charging.

Overview


NeoBeat *

Note

ECG electrode consists of metal disc and the black material (conductive TPU).


* Patient Applied Part

Charging

First time use

To initiate the device, NeoBeat must be charged for up to 3 hours before first use. Place it on the Charging Stand to charge the battery.


Notes

- *Maintenance charging: If NeoBeat is not stored on the Charging Stand, recharge it at least monthly.*
- *The USB port on the Charging Stand can be used to connect one (and only one) additional Charging Stand.*
- *It may take up to 30 minutes before the device indicates charging if the battery is very depleted.*
- *If a device does not start, perform reset procedure (p. 42) and leave it charging for another 30 minutes.*
- *NeoBeat is activated by motion. When not in use, store it on the Charging Stand to avoid unnecessary battery use.*
- *Device cannot be used clinically while charging*

Clinical Use

Remove NeoBeat from the Charging Stand and place it around the newborn's torso. NeoBeat automatically turns on when it detects motion.



✗ Do not over-bend.

⚠ Caution

Do not bend the device more than necessary for application to the newborn. Excessive bending may cause the device to become permanently deformed and damaged.

2

Ensure that both electrodes have good contact with patient's skin.

Notes

- *NeoBeat's position on the torso can be shifted so that its placement does not interfere with other therapies, e.g. chest compressions, palpation, auscultation, umbilical access.*
- *NeoBeat can also be repositioned on the torso to improve signal quality.*


Activated - No contact

NeoBeat is activated, but there is no or inadequate patient contact. If no patient contact, the device will return to standby automatically after 10 seconds.

Note

Skin contact is necessary to measure heart rate. If the patient skin is too dry, add a drop of water under the electrodes to achieve contact or improve performance of the device.

Initial calculation

NeoBeat has patient contact and is calculating the heart rate.

Heart rate unknown

Heart rate cannot be detected. This may be due to poor contact of the ECG electrodes, or lack of detectable heart rate despite having good contact.

Note

Reposition device if no heart rate can be obtained. Moving the device closer to the heart may give a stronger signal.

Clear detection

Heart rate is detected and there is good signal quality.

Weak detection

A typically transient condition where heart rate is detected, but with reduced signal quality. The heart rate displayed may be less accurate during the weak detection display.

If this display persists, ensure that the electrodes of NeoBeat are in full contact with the skin of the newborn's torso. Adding a drop of water under each electrode will improve signal quality if the skin is too dry. Moving the electrodes closer to the heart may give a stronger signal.

Undetectable due to motion

There is too much motion to detect heart rate.

When excessive movement (e.g. stimulation) is detected, the device will not show a heart rate.

Low battery

From when low battery is first indicated, the device will have approximately 30–60 minutes remaining run time. Recharge NeoBeat after use and reprocessing, by placing it on the Charging Stand.

Error display

NeoBeat has detected an internal technical error. See troubleshooting section for details.

Use other methods for assessing newborn heart rate, e.g. by auscultation or palpation.

Reprocessing

Clean and disinfect NeoBeat after each patient use to minimize the risk of cross-contamination.

Warning

Do not place a used NeoBeat back onto the charging stand before it has been cleaned and disinfected.

Cleaning

1. Rinse NeoBeat under running lukewarm water for 15 seconds.
2. Clean all surfaces of NeoBeat using a cloth dampened with lukewarm (30 to 40 °C, 86 to 104 °F) tap water and mild dishwashing detergent for 30 seconds. To remove difficult soil, use a bristled brush (e.g. toothbrush) dipped in the cleaning solution. Clean for a minimum of 2 minutes, ensuring that all soil has been removed. Visually inspect NeoBeat for soil. Repeat this cleaning step in the case that residual soil is still visible on NeoBeat.
3. Rinse NeoBeat under running lukewarm water, for 15 seconds. Repeat once.
4. Dry NeoBeat using a clean cloth or by air drying.

Disinfection

1. Wipe all surfaces of NeoBeat with a clean cloth soaked with 70% ethanol or 70% isopropanol for a minimum of 2 minutes.
2. Spray 70% ethanol or 70% isopropanol on all surfaces of NeoBeat. Ensure it remains wet for a minimum of 12 minutes. Repeat spraying as necessary to account for evaporation.
3. Allow to air dry.

Cautions

- *Effective disinfection is not possible without first performing a thorough cleaning.*
- *Care should be taken while handling the product between cleaning and disinfection.*
- *Do not submerge any of the product components in liquid (including ethanol or isopropanol).*
- *Do not use sodium hypochlorite (bleach) to disinfect.*
- *The reprocessing method is designed to ensure NeoBeat is adequately disinfected between uses. Any deviation will increase the risk of cross-contamination especially for newborns that may have compromised immune defense, such as a pre-term baby or in the case of outbreaks of highly transmissible pathogens.*

Reprocessing

Inspection

After reprocessing, inspect NeoBeat for cracks and damage with particular attention to the highlighted areas.

If there is any damage, remove the device from service. Otherwise, put the device back onto the Charging Stand.


Note



Other than the visual inspection after reprocessing, there is no need for regular functional testing of NeoBeat. NeoBeat runs a self-test when activated.

Charging Stand

When needed, the Charging Stand can be cleaned and disinfected using the same method as previously described.

NeoBeat Charging Stand can be placed on a tabletop. NeoBeat Charging Stand can optionally be mounted vertically to a wall, rail or post.

Caution

Do not mount or place the Charging Stand in close proximity to heat sources (e.g. directly below an infant warmer heat lamp).

Ensure that the mounting of NeoBeat Charging Stand to any surface is safe and is performed by a competent person.



Installation




Vertical Mount

Use double sided foam tape (not included) or screws to mount the Mounting Bracket.

Multiple mounting screw hole patterns:

9 × 9 mm, 21 × 21 mm, 75 × 75 mm.

Service and Warranty

NeoBeat does not have any serviceable parts.

NeoBeat has a one-year limited warranty. Refer to company website for terms and conditions.

Main Device and System:

532-01001 NeoBeat Main Device

531-01001 NeoBeat Mini Main Device

532-00001 NeoBeat System

531-00001 NeoBeat Mini System

Accessories

532-21050 Power Supply and Plug Kit

532-20050 NeoBeat Charging Stand

Symptom/Display	Possible cause	Possible solution
<p>Patient contact is not detected.</p>  <p>The image shows the NeoBeat device's display with the LAERDAL logo at the top and NeoBeat at the bottom. The central display area is mostly dark, with only a few scattered white dots, indicating that no patient contact has been detected.</p>	<p>Poor or no contact with skin.</p> <p>Or</p> <p>The skin is too dry for NeoBeat to detect patient contact.</p>	<p>Check for barriers or obstruction to skin contact, e.g. towel.</p> <p>Wet the electrodes with water and/or reposition the device.</p>
<p>Heart rate is not detected.</p>  <p>The image shows the NeoBeat device's display with the LAERDAL logo at the top and NeoBeat at the bottom. The central display area shows a pattern of white dots forming a heart rate waveform, but it is incomplete and sparse, indicating that a heart rate is not being detected.</p>	<p>The device is poorly positioned or not in direct skin contact with patient.</p> <p>Or</p> <p>Stimulation/movement/handling of patient temporarily generating too much disturbance.</p> <p>Or</p> <p>Heart rate is outside detectable range.</p>	<p>Reposition the device. Moving the device closer to the heart may give a stronger signal.</p> <p>If the problem continues, use alternative means of measuring newborn heart rate, e.g. a stethoscope.</p>

Troubleshooting

Symptom/Display	Error codes	Possible solution
<p>Critical technical error: Error code is shown in the display and error light activated.</p> 	<ul style="list-style-type: none">x01 - Program memoryx02 - Calibration memoryx04 - Data memoryx08 - RTC crystalx10 - Display driver communicationx20 - Accelerometer communicationx40 - Light sensor communicationx80 - ECG analog signal chain	<p>Continue clinical procedures without use of the device.</p> <p>Use alternative means of measuring newborn heart rate, e.g. a stethoscope.</p> <p>Place the device on the charging stand. If the problem persists, remove device from service.</p>

Troubleshooting

Symptom	Possible cause	Possible solution
NeoBeat does not turn on or it turns off during use.	Battery depleted.	<p>Continue standard clinical procedure without use of the device.</p> <p>After the procedure, reprocess and charge the device.</p> <p>It may take up to 30 minutes before the device indicates charging if the battery is very depleted.</p> <p>If a device does not start, perform reset procedure (p. 42) and leave it charging for another 30 minutes.</p>
Nonfunctional or damaged device detected during equipment inspection.		Remove the device from service. A replacement NeoBeat may be ordered.

Specifications

Environmental	
Temperature	Operating: 0 – 40 °C (32 – 104 °F)
	Storage / shipping: Short-term: -20 - 60 °C (-4 - 140 °F) Long-term: 15 - 25 °C (59 - 77 °F)
Atmospheric Pressure	620 – 1060 hPa (up to 4000 meters above sea level)
Relative Humidity	Operating: 15% - 90%, non-condensing
	Storage / Shipping: 15% - 90%, non-condensing

Heart Rate Meter	
Dimensions	NeoBeat: 83 x 87 x 40 mm (3.2 x 3.4 x 1.6 inches) NeoBeat Mini: 70 x 70 x 40 mm (2.8 x 2.8 x 1.6 inches)
Weight:	NeoBeat: 31 g (1.1 oz) NeoBeat Mini: 27 g (1 oz)

Materials	
NeoBeat Body	Polyamide Conductive TPU
NeoBeat Metal Electrodes	Stainless steel

Power	
Battery	Internal rechargeable lithium-ion button cell, 3.7V, 120 mAh Typical service life of battery: 3 – 6 years depending on use
Run time	>4 hours (full charge on new battery) >3 hours (full charge at expected end of battery service life depending on use and if properly maintained)
Charge time	Up to 3 hours (full charge of empty battery)
Power supply	Input 100 – 240V AC, 50 – 60 Hz, 0.3 A Output 5V DC, 1 A

Lifetime parameters	
Shelf life	3 years
Expected Service Life	NeoBeat - typically 100 cycles of reprocessing
	NeoBeat Mini - typically 100 cycles of reprocessing
	Ensure that the device is inspected after reprocessing in accordance with the Reprocessing Inspection on page 20.

Caution

Only use provided power supply, PSAI05R-050QL6-R, or an alternative 5 W USB power supply, 5 V DC, 1 A, that is IEC 60950-1, IEC 62368-1 or IEC 60601-1 certified.

Heart Rate Measurement	
Accuracy	Short term average $\pm 10\%$ or ± 5 bpm, whichever is greater, in the range 30 – 250 bpm and 0.2 - 5 mV QRS amplitude. No detectable heart rate is displayed as “-?-“

Specifications

Classification	
Ingress Protection	Heart rate meter: IP55 - Protected against ingress from dust and water jet spraying. All other components: IPX0 - Not protected against liquid ingress.
IEC 60601 - 1	Internally powered/class II equipment type BF

 Warning

Do not modify this equipment without authorization of the manufacturer.

 Caution

Changes or modifications not expressly approved by Laerdal Medical could void the user's authority to operate the equipment.

Technical Description

Federal Communications Commission (FCC) and Industry Canada (IC) Statement

This device complies with part 15 of the FCC Rules and Industry Canada's licence-exempt RSSs. Operation is subject to the following two conditions:

This device may not cause harmful interference, and this device must accept any interference received, including interference that may cause undesired operation.

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FCC ID: QHQ-20-09917

IC: 20263-2009917

Symbol Glossary	
	Medical Device
	This appliance is marked according to the European directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE).
	Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
	Protected against ingress from dust and water jet spraying
	Meets IEC type BF applied part leakage current requirements
	Manufacturer
	Consult User Guide
	The products is eligible to bear the CSA Mark with adjacent indicator US for US

Technical Description

Symbol Glossary	
	Temperature limitation
	Atmospheric pressure limitation
	Humidity limitation
XXX 	Machine readable Unique Device Identification (UDI). The last three digits of the UDI are printed above the machine readable UDI for easier distinction between devices.
	MR unsafe

Waste Handling

This appliance is marked according to the European directive 2012/19/EC on Waste Electrical and Electronic Equipment (WEEE). By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product. The symbol on the product, or on the documents accompanying the product, indicates that this appliance may not be treated as household waste. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment. Disposal must be carried out in accordance with local environmental regulations for waste disposal. For more detailed information about treatment, recovery and recycling of this product, please contact your local city office, your household waste disposal service or Laerdal representative.

Electromagnetic Conformity

NeoBeat is IEC 60601-1-2 compliant and intended for use in the following environments: Health care facilities except for near HF surgical equipment and the RF shielded room for magnetic resonance imaging.

No particular actions are required to maintain safety and performance with regard to electromagnetic disturbances for the expected service life.

NeoBeat was tested according to the recommendations of IEC TR 60601-4-2: Electromagnetic immunity: performance of medical electrical equipment and medical electrical systems.

Warnings

- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the NeoBeat, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical Description

Electromagnetic Emissions Tests

Emissions Test	Standard or test method	Compliance
Radiated RF emissions	CISPR 11	Group 1 Class B
Conducted RF emissions	CISPR 11	N/A
Harmonic distortion	IEC 61000-3-2	N/A
Voltage fluctuations/ flicker	IEC 61000-3-3	N/A

Electromagnetic Immunity Tests

Immunity Test	Standard or test method	Immunity Test Level
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air maximum recovery time following a TRANSIENT phenomenon: 2s
Radiated RF EM fields	IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz 80% AM at 2 Hz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	380-390 MHz: 27 V/m 430-470 MHz: 28 V/m 704-787 MHz: 9 V/m 800-960 MHz: 28 V/m 1700-1990 MHz: 28 V/m 2400-2470 MHz: 28 V/m 5100-5800 MHz: 9 V/m

Immunity Test	Standard or test method	Immunity Test Level
Rated power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz
Electrical fast transients / bursts	IEC 61000-4-4	N/A
Surges: Line-to-line	IEC 61000-4-5	N/A
Surges: Line-to-ground	IEC 61000-4-5	N/A
Conducted disturbances induced by RF fields	IEC 61000-4-6	N/A
Voltage dips	IEC 61000-4-11	N/A
Voltage interruptions	IEC 61000-4-11	N/A

Technical Description

Performance

ECG Performance Specification as required by IEC 60601-2-27	
Leads-off detection	For sensing leads-off, NeoBeat applies a 200 nA p-p 100 kHz square wave current through the two sense electrode (electrodes on arms)
Noise suppression	The right-leg-drive amplifier through the reference electrode (middle electrode) limits the current to be <5 μ A even in single fault condition
Tall T-wave rejection capability	1.2 mV
Heart rate averaging method	Normally, heart rate is $\frac{1}{2}$ *(mean of 12 most recent included RR intervals) + $\frac{1}{2}$ *(median of 12 most recent included RR intervals). During "weak detection" periods, heart rate is median of 12 most recent included RR intervals.
Display update rate	2 s
Heart rate meter accuracy and response to irregular rhythm	A1 Ventricular bigeminy: 80 bpm A2 Slow alternating ventricular bigeminy: 60 bpm A3 Rapid alternating ventricular bigeminy: 60 bpm A4 Bidirectional systoles: 60 bpm
Response time of heart rate meter to change in heart rate	HR change from 80 to 120 bpm: 6 s HR change from 80 to 40 bpm: 14 s
Measured accuracy of detected heart rate	$\leq \pm 1\%$ or ± 1 bpm, whichever is greater; when tested in accordance with Clause 201.12.1.101.15.


Clinical Accuracy

Data from clinical testing show that NeoBeat has an RMS (root mean square) heart rate accuracy of ± 3 bpm during Clear Detection (i.e., with good signal quality) and ± 6 bpm during Weak Detection (i.e., reduced signal quality). The estimates are based on 19 clinical cases, including resuscitation cases, representing over 4 hours of ECG data. The 19 cases were randomly selected from a large database containing newborn ECGs from four countries outside the United States. The database contains over 1,000 cases of ECGs collected from researchers using NeoBeat on both term and pre-term newborns within the first few minutes of birth, including resuscitation cases. The QRS complexes of the 19 cases were then annotated by two experts, an ICU physician and a scientific expert in ECG signal processing and analysis. Heart rate based on expert annotation was considered “ground truth”. Heart rate based on NeoBeat’s algorithm was then compared to the “ground truth” heart rate to determine accuracy.

Technical Description



Machine Readable Unique Device Identification (UDI)

The GS1 DataMatrix located on the outside of the electrode arm of NeoBeat contains the UDI of the device. This barcode can be read using apps on a smart phone or tablet.

<p>Identify the 3-digit short SN (serial number) and the GS1 2D datamatrix:</p>	<p>On a GS1 compliant scanner, or with a smartphone with a suitable App, capture the datamatrix.</p>	<p>For manual reading of the code, copy the numeric string that represents the UDI identifier:</p>	<p>For manual decoding of the serial number; the UDI code can be read as below example</p>
 <p>Example: Short SN printed on device</p>	 <p>Example</p>	 <p>Example</p>	<p>UDI code(example): 0107045420815342 134215312200033</p> <p>SN is the last 10 digits: 531 22 00033</p> <p>531 = NeoBeat Mini 532 = NeoBeat</p> <p>22 = manufacturing year</p> <p>00033 = unit production number; (033 is the short SN on the device)</p>

Information Available Electronically

Hold NeoBeat upside down and double tap firmly on the charging pin side to display the unique device identification (UDI), FCC ID, IC certification number and software version.

UDI 0107045420815342134215312200033

Technical Description

Bluetooth® Low Energy transmitter

Frequency band: 2.400 – 2.4835 GHz

Modulation: Gaussian frequency shift modulation

Maximum radio-frequency power transmitted: 1 mW

Effective radiated power: 0 dBm

Restart NeoBeat

NeoBeat can be restarted by placing NeoBeat on and off the charging stand 10 consecutive times and then leaving it on the charging stand.

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© 2024 Laerdal Medical AS. All rights reserved.
NeoBeat is protected by US and international registered patents and design rights. Laerdal® and NeoBeat® are registered trademarks of Laerdal Medical AS.