EN

FR

DE

ES

IT

BR

NL

PL

RU

JA

ZH

KO

SimMan ALS

Important Product Information
重要产品信息

www.laerdal.com

Laerdal

helping save lives

Read these instructions thoroughly. Observe all warnings, cautions and instructions in the User Guide and in this Important Product Information booklet. Retain this booklet for future reference.

Warnings and Cautions

A Warning states a condition, hazard, or unsafe practice that can result in serious personal injury or death.

A Caution states a condition, hazard, or unsafe practice that can result in minor personal injury or damage to the product.

Notes

A note states important information about the product or its operation.

Disclaimer

Use of the SimMan ALS Patient Simulator system to train personnel should be undertaken under supervision of suitably trained medical personnel with an understanding of educational principles as well as recognized medical protocols.

As with all Patient Simulators or other such training devices there may be approximations, variations and inaccuracies in anatomical features and the physiological modeling.

This being the case, Laerdal Medical does not guarantee that all features are completely accurate.

The medical equipment and simulated medical equipment included in the product might be modified and should be used for training purposes only.

General Simulator Handling

It is important to follow the instructions below, as well as other available User Information, in order to maintain optimum performance and longevity of the simulator components.

Warnings

- Do not provide artificial respiration to the Patient Simulator using oxygen enriched air or flammable gases.
- The transportation cases are heavy. Lifting it alone may lead to muscle strains.

Cautions

- Do not introduce fluids into or onto the patient simulator (except as directed in the User Guide) as this may damage it and its components.
- Do not put biological materials or other materials in the simulator's airways. They are not designed for cleaning or disinfection.
- Do not introduce humidified air into the system during ventilation.
- Never perform mouth-to-mouth or mouth-to-nose ventilation on the Patient Simulator. The airways are not designed for cleaning or disinfection.
- Do not use the Patient Simulator if the internal tubing and cabling is disconnected.
- Do not store or use this product outside its specified operating and storage conditions.
- Do not use automated chest compression machines on the Patient Simulator.
- Do not use the Patient Simulator in aircrafts when they are airborne.
- Do not store the simulator with fluid in the IV Arm System. Use a syringe to remove remains of any injected water from the tubing and components prior to storage.

English

SimMan ALS - Important Information

Note

When used in a humid environment, some drops of water may leak from the manikin's right lower leg during startup and shutdown. This is normal and does not require assistance.

Defibrillation

Warnings

- Do not attach medical device defibrillator electrodes or place defibrillator paddles on any part of the simulator.
- Do not defibrillate on the ECG connectors or other conductive parts.
- Do not defibrillate the Patient Simulator in a flammable or oxygen-enriched atmosphere. This is a fire hazard.

For more information refer to the Important Product Information booklet enclosed with Shocklink.

Note

Electrical arcing may occur during defibrillation. This happens when the pads are not stuck properly to the manikin, the gel can become overheated and produce a strong smell. This is not harmful to neither the user nor the product.

Patient Simulator Power - Battery

Warnings

- Only use the designated Laerdal approved battery in SimMan ALS.
- Do not use the Laerdal Li-Ion Battery for other purposes than specified.
- Use only the AC charger adapter supplied with the product.
- Do not use the battery in temperatures above what is specified in this Important Product Information.
- Do not heat or incinerate. Do not crush the battery.
- Do not short circuit the battery contact.
- Do not immerse in water.

Note

The Patient Simulator contains electronic components. All components, including the battery must be recycled and disposed of in compliance with applicable laws and regulations.

For other battery related information, consult User Guide.

Patient Simulator - External Power

Warnings

- Only connect the Patient Simulator to a Laerdal approved external power supply dedicated to use on SimMan ALS.
- Do not use the Patient Simulator close to real patients when it is connected to the external power supply.

Air Tank

Warnings

- Do not disassemble or attempt to repair the internal parts of the right thigh.
- Do not puncture.

Torso Skin

Warnings

- Ultrasound Live Scan (LS) tags must not be used on anyone using a pacemaker, implantable cardioverter defibrillator, or other electronic medical devices. Components of the technology may interfere with such medical devices
- Do not use ultrasound gel and avoid exposing the LS Tags to any liquids.

Care, Maintenance and Cleaning

- Keep the simulator clean and in a hygienic condition.
- Clean all skin parts regularly.
- Use warm soapy water or Laerdal Manikin Wipes.

Federal Communications Commission (FCC) and Industry Canada (IC) Statements

This device complies with part 15 of the FCC Rules and RSS-210 of IC rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Changes or modifications not expressly approved by Laerdal Medical could void the user's authority to operate the equipment.

The term "IC" before the equipment certification number only signifies that the Industry Canada technical specifications were met.

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This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

1. This device may not cause interference.
2. This device must accept any interference, including interference that may cause undesired operation of the device.

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This Class B digital apparatus complies with:

Link Box PLUS 2	
FCC ID:	QHQ-LB1YM
IC:	20263-LB1YM

EU

This product is in compliance with the essential requirements of Council Directive 2014/53/EU on Radio Equipment (RED). The product is in compliance with Council Directive 2011/65/EU on restriction of the use of certain hazardous substances (RoHS).

Waste Handling

Dispose of in accordance with your country's recommendations. This appliance is marked according to the European directive 2012/19/EC on Waste Electrical and Electronic Equipment (WEEE).

By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

The symbol on the product, or on the documents accompanying the product, indicates that this appliance may not be treated as household waste. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment. Disposal must be carried out in accordance with local environmental regulations for waste disposal.

For more detailed information about treatment, recovery and recycling of this product, please contact your local city office, your household waste disposal service or Laerdal representative.

Battery

SimMan is operated on a Lithium Ion (Li-Ion battery). Li-Ion batteries should be recycled or disposed of in accordance with local regulations.

Symbol Glossary

	CE mark
	Korean Certification (KC) Mark
	Australia Radiocommunications and EMC Compliance Mark
	Waste Electrical and Electronic Equipment symbol
	Warning/ Caution
	Note
	Manufacturer
	Date of Manufacture
	Reference Number
	Serial Number
	Li-Ion batteries recycling symbol

Specifications

Size and Weight	
Dimensions (Patient Simulator only)	1800 mm (l) x 550 mm (b) chest (5.90 ft x 1.80 ft)
Weight (Patient Simulator only)	33 kg (73 lbs)
Weight (with clothes)	34 kg (75 lbs)
Patient Simulator Power	
External Power	Input voltage 12VDC, 3,33A, 40W
Internal Battery (one)	Each 7.34V, 5700mAh, 42Wh, Lithium-Ion
Battery Time	2 hours
Temperature Limits	
Operating Temperature	0 °C to 40 °C (32 °F to 104 °F)
Storage Temperature	-15 °C to 50 °C (5 °F to 122 °F)
Battery Charging	10 °C to 25 °C (50 °F to 77 °F)
Environment - Patient Simulator Only	
Relative Humidity	20% -90% (non-condensing)
Do not use outdoors in wet conditions	
Not tested with salt spray	
RF Communication	
WLAN Frequency Ranges	- 2.4 GHz WLAN, channels 1-11 - 5 GHz WLAN, channels 36, 40, 44 and 48
WLAN Operation Range	100 m (300 ft) outdoors
Pressure Data (Compressor)	
Working Pressure	1.5 bar
Working Temperature	0 °C to 40 °C (32 °F to 104 °F)
Capacity of the vessel (V in l)	0.3 l

Material Chart for Patient Simulator	
Clothes	Cotton, Nylon
Skins and airways	PVC (DEHP free)
External hard plastics	PP, PA, PC, PC/PET
Inner plastics	Silicone, TPU, TPE, PVC, Nitrile PA, PA+GF, PC, ABS, POM, HDPE, PET
Metal components	Aluminium, Brass, Steel
Recommended Hardware Specifications	
https://laerdal.my.site.com/HelpCenter/s/article/Will-there-be-a-compatibility-matrix-to-help-identify-if-the-hardware-is-compatible-with-LLEAP	
	
Minimum Software Requirements	
SimPad ver: 9.0.	
LLEAP ver: 8.6.1.	
Simulated Medications and IV Fluids	
Only use purified water for simulated medications and IV fluids.	

Warranty

Refer to the Laerdal Global Warranty for terms and conditions.
For more information visit www.laerdal.com

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© 2024 Laerdal Medical AS. All rights reserved.
SimMan® and Laerdal® are registered trademarks of Laerdal Medical AS.