
# The BAG™ Bag Refill Valve Adapter

## User Guide

www.laerdal.com

### INSTRUCTIONS FOR USE

#### A Overview

1. Bag Refill Valve Adapter
2. Connector (ID 32 mm)
3. Silicone seal

#### B Preparation of The BAG

1. Disconnect The BAG's oxygen tube.

#### Note

The Oxygen tube requires a tool, e.g. a flat screwdriver blade, in order to be loosened from The BAG Resuscitator's Oxygen Inlet Connector. The blade so that the Oxygen tube's connector is lifted and allows partial detachment from the resuscitator's oxygen inlet connector.

- a) Orient the resuscitator so that the Oxygen Tube's connector is accessible.
- b) Insert the blade into the gap just below the Oxygen tube's connector.
- c) Twist the blade so that the Oxygen tube's connector is lifted and allows partial detachment from the resuscitator's oxygen inlet connector.
- d) Manipulate and pull the Oxygen Tube's connector so that it detaches.

#### C Attaching and Testing

1. Align the Bag Refill Valve Adapter with The BAG inlet valve oxygen ports.
2. Press The Bag Refill Valve Adapter firmly onto The BAG's inlet valve so that it is locked in place.
3. Check that the Bag Refill Valve Adapter is solidly connected.
4. Block The Bag Refill Valve Adapter's open connector with a hand.
5. Squeeze the resuscitator bag fully.
6. Observe that the resuscitator bag does not re-expand until the hand is removed. If it fails, remove the Bag Refill Valve Adapter and check the Silicone Seal's position.

#### Warning

If the bag refill valve adapter is not correctly connected to the resuscitator, it may lead to leakage of ambient air into the breathing gas and reduced oxygen content provided to the patient. Ensure that the bag refill valve is activated for each ventilation.

#### D Attaching and Testing with Bag Refill Valve

1. Connect the bag refill valve.
2. Test function by squeezing the resuscitator bag fully.
3. Observe if the resuscitator bag re-expands quickly with oxygen from the bag refill valve.

#### Warnings

- A resuscitator bag that re-expands slowly or not at all indicates that the oxygen supply from the bag refill valve is too low or non-existent. Remove the bag refill valve and proceed with ventilation with ambient air.
- Check that the bag refill valve is activated for each ventilation.

#### E Revert to using The BAG as normal

1. Remove the bag refill valve.
2. Follow the instructions on the Bag Refill Valve Adapter to remove it by squeezing the two areas marked with "release-indicators" areas on its side edge while pulling the adapter off the resuscitator.
3. Reattach The Oxygen Reservoir Bag and Oxygen Tube onto The BAG's inlet valve. Connect to an oxygen source as described in The BAG's user guide.

#### Using the Bag Refill Valve

1. Remove The BAG's Oxygen Reservoir Bag and Oxygen Tube (Fig. B).

#### Note





In case of needing to re-apply the original components, keep them protected from contamination.

1. Attach The Bag Refill Valve Adapter to The BAG's inlet valve (Fig. C).
2. Perform a function test and attach the bag refill valve (Fig. D).
3. Ventilate as clinically indicated with the Bag Refill Valve.

### SPECIFICATIONS

#### Standards

- The BRVA complies with the following standards:
  - ISO 10651-4:2002 - Lung ventilators Part 4: Particular requirements for operator-powered resuscitators
  - ISO 10993-1:2018 Biological evaluation of medical devices - Part 1

### ISO 18562-1:2017 Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 11: Evaluation and testing within a risk management process

### EN 1789:2020 - Medical valves and their equipment

Environmental conditions	
Storage	Extreme temperatures, short term: -40 °C to 70 °C (-40 °F to 140 °F) Non-condensing
Operating Temperature	-18 °C to 50 °C (0.4 °F to 122 °F)

#### Weight / Dimensions / Volume (approximate)

BRVA	Weight: 0,33 kg / 0,73 lbs Length: 460 mm / 18,1 inches Height: 135 mm / 5,3 inches Width: 130 mm / 5,1 inches
------	-------------------------------------------------------------------------------------------------------------------------

#### Disposal: In accordance with local protocol.

#### Materials

Valve housings	Polycarbonate (PC)
Seal	Silicone rubber

#### Symbol Glossary

CE 2460	This medical device complies with the particular requirements of Regulation (EU) 2017/745 for medical devices.
MD	Medical Device
UDI	Unique Device Identifier
LATEX	Not made with natural rubber latex
⊗	Do not re-use
1P	Single patient multiple use
ⓘ	Consult Operators Manual

#### Global 1-Year Limited warranty: www.laerdal.com

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### CLINICAL INDICATIONS

#### Device Description

The Bag Refill Valve Adapter (BRVA) is a non-sterile single patient use accessory for provision of oxygen from a Bag Refill Valve to Laerdal The BAG Resuscitator. It may be used for single-patient multiple use when kept free from contamination.

#### Indication for Use

The Bag Refill Valve Adapter is indicated for use when oxygen is pulled into the resuscitator from a bag refill valve (not provided by Laerdal).

#### Intended Use

The Bag Refill Valve Adapter is intended for attachment of a bag refill valve as an oxygen source into the resuscitator.

#### Intended Users

Healthcare professionals trained in resuscitation, delivering ventilatory support and in the use of manual resuscitators.

#### Clinical Benefits

Positive impact on clinical outcome, by respiratory support that reduces probability of adverse outcomes, such as morbidity and mortality caused by hypoxia.

#### Clinical Outcome

Desired outcome of ventilation is oxygenation of the patient, often evaluated using SpO<sub>2</sub>, EtCO<sub>2</sub>, blood gas analysis or other method of analysis.

#### Known Side Effects

Oxygen toxicity

#### Contraindications

None known.

### IMPORTANT INFORMATION

Read this User Guide and become familiar with the operation of the BRVA prior to use. Use the BRVA only as described in this User Guide.

#### Warnings

A Warning states a condition, hazard, or unsafe practice that can result in serious personal injury or death.

- The BRVA is for single patient use only. It is not designed for reprocessing. Replace the BRVA with a new one if it has come in contact with any fluids, mucus or vomit. Reuse with multiple patients will lead to risk of cross contamination.
- The BAG Resuscitator's Oxygen Reservoir Bag and Oxygen Tube must be removed before the BRVA can be fitted.
- Perform a function test before use.
- In case of a function of the bag refill valve is compromised, remove it immediately from the resuscitator to allow respiration with ambient air. Re-attach The BAG's oxygen tube and oxygen reservoir bag for enabling provision of oxygen as described in The BAG's user guide.

#### Cautions

A Caution states a condition, hazard, or unsafe practice that can result in minor personal injury or damage to the product.

- The BRVA should only be used by persons who have received adequate training in the use of resuscitators and bag refill valves.
- The use of third-party products and oxygen delivery devices with The BAG may change product performance. Please consult with the manufacturer of the third-party device to verify compatibility with The BAG and obtain information on possible performance changes.

#### Notes

Important information about the product or its operation.

- Should any serious malfunction, undesirable incident with, or deterioration in the functionality or performance of the device occur, contact Laerdal promptly. The competent authority where the incident took place and/or the device was used should also be notified.

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