EN

FR

DE

ES

IT

PT

NL

PL

JA

ZH

KO

MamaAnne

Important Product Information

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Laerdal

helping save lives

**limbs&
things**

Read these instructions thoroughly. Observe all warnings, cautions and instructions in the User Guide and in this Important Product Information booklet. Retain this booklet for future reference.

Warnings and Cautions

A Warning states a condition, hazard, or unsafe practice that can result in serious personal injury or death.

A Caution states a condition, hazard, or unsafe practice that can result in minor personal injury or damage to the product.

Notes

A note states important information about the product or its operation.

Disclaimer

Use of the MamaAnne system to train personnel should be undertaken under supervision of suitably trained technical or medical personnel with an understanding of educational principles as well as recognized medical protocols.

As with all Patient Simulators or other such training devices, there may be approximations, variations and inaccuracies in anatomical features and physiological modeling.

Laerdal does not guarantee that all features are completely accurate.

General Simulator Handling

It is important to follow the instructions below, as well as other available user information, in order to maintain optimum performance and longevity of the simulator components.

General Operation

Cautions

- Avoid pinch hazards – Do not use the Patient Simulator without the external skins.
- Avoid contact with sharp objects or surfaces and excessive pulling on skin.

- Removing jewelry is recommended to prevent damage to soft skins.
- Do not leave the Patient Simulator unattended when sitting up unassisted. This poses a fall hazard, risking damage to simulator or equipment in area and injury to personnel.
- Do not use the Patient Simulator in aircrafts when they are airborne.
- Do not use outdoors in wet conditions. Not tested with salt spray.
- Do not perform mouth-to-mouth resuscitation.
- Use minimum of two people when lifting the Patient Simulator.
- Do not perform chest compressions while spontaneous breathing is enabled as this can damage internal parts.

Fluids

Cautions

- Do not introduce fluids into or onto the Patient Simulator (except as directed in the User Guide) as this may damage it and its components.
- Follow cleaning protocols as instructed in the User Guide. Failure to do so may damage the product.
- Do not use the Patient Simulator if the internal tubing and cabling is disconnected.
- Do not use mineral oil, MEK, or acetone as cleaning agents.
- Do not rotate the Patient Simulator if fluid has leaked into the birthing chute or in the lower pelvis.

CPU Module (Component of simulator)

Cautions

- The device for operation in the band 5150-5250 MHz is only for indoor use to reduce the potential for harmful interference to co-channel mobile satellite systems.
- The maximum antenna gain permitted for devices in the bands 5250-5350 MHz and 5470-5725 MHz shall comply with the e.i.r.p. limit.
- The maximum antenna gain permitted for devices in the band 5725-5825 MHz shall comply with the e.i.r.p. limits specified for point-to-point and non point-to-point

operation as appropriate.

- The worst-case tilt angle(s) necessary to remain compliant with the e.i.r.p. elevation mask requirement set forth in Section 6.2.2(3) shall be clearly indicated.
- Users should also be advised that high-power radars are allocated as primary users (i.e. priority users) of the bands 5250-5350 MHz and 5650-5850 MHz and that these radars could cause interference and/or damage to LE-LAN devices.
- Do not interrupt or power down the system during a software update.

Battery

Warnings

- Only use the designated Laerdal approved battery in the Patient Simulator.
- Do not use the Laerdal Li-ion battery for other purposes than specified.
- Do not use the battery in temperatures above what is specified in this Important Product Information.
- Do not heat or incinerate. Do not crush the battery.
- Do not short circuit the battery contact.
- Do not immerse in water.
- Inserting and connecting battery incorrectly, short circuiting or exposure to fluids pose an explosion hazard.
- Do not mistreat, disassemble or attempt to repair the battery.
- Do not use the battery if it is visibly damaged, malfunctions or appears to leak electrolyte.
- Take extreme care to avoid direct contact with electrolyte, hot or smoking parts. In case of the above, disconnect and remove the battery when it is judged safe to do so

External Power

Warnings

- Only connect the Patient Simulator to a Laerdal approved external power supply dedicated to use on MamaAnne.
- Do not use the Patient Simulator in close proximity to real patients when it is connected to the external power supply.
- Patient Simulator should be connected to an

electrical outlet that is properly grounded.

- Use only the external power supply supplied with the product.
- The external power supply is for indoor use only.

Air Tank

Warnings

- Do not disassemble or attempt to repair the internal parts of the right thigh.
- Do not puncture.

Defibrillation

MamaAnne allows for defibrillation in accordance with AHA 2020 international guidelines for CPR when it's connected to ShockLink compatible defibrillators.

An approved conventional defibrillator may be used only with ShockLink on the Patient Simulator. During live defibrillation with ShockLink, the defibrillator and Patient Simulator may present a shock hazard. All standard safety precautions must be taken when using the defibrillator. For more information, consult your defibrillator's User Guide.

Warnings

- Only use defibrillators approved for use with Laerdal ShockLink. For more information, see <https://laerdal.com/products/skills-proficiency/defibrillation-cardiology/shocklink/>
- Observe normal safety precautions when using defibrillators.
- Defibrillation must be performed using ShockLink only. Refer to ShockLink Important Product Information.
- Do not attach medical device defibrillator electrodes or place defibrillator paddles on any part of the simulator.
- Do not attach medical device defibrillator electrodes or places defibrillator paddles on ECG monitoring leads.
- Do not defibrillate the Patient Simulator with ShockLink if it is disassembled, wet or damaged.
- Do not defibrillate the Patient Simulator with ShockLink when it is turned off.

- Do not defibrillate the Patient Simulator with ShockLink in a flammable or oxygen-enriched atmosphere. This is a fire-hazard.

ESD

Warning

- The production of sparks or fire may occur when cleaning the silicone skins with flammable chemicals.

Oxygen

Warning

- Oxygen should not be administered to the Patient Simulator. There is a risk of fire.

Birthing

Warning

- Do not operate the birthing mechanism without abdomen installed.
- Do not insert hands into the birthing mechanism past the soft birth canal while it is moving.
- Keep hands clear of the birthing chute proximal wall when the delivery mechanism is returning to home position.
- Do not insert fingers into the birthing mechanism chute to cover gap when the mechanism is moving.

Waste Handling

Recycle and dispose of the product in accordance with your country's recommendations.


By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

For more detailed information about treatment, recovery and recycling of this product, contact your local waste disposal service or local governmental office for waste management.

Battery

MamaAnne is operated on a Li-ion battery. Li-ion batteries should be recycled or disposed of in accordance with local regulations.

Symbol Glossary

	Waste Electrical and Electronic Equipment symbol
	CE Mark
	Australian Radiocommunications and EMC Compliance Mark
	MIC Technical Conformity Mark (Japan)
	Korean Certification (KC) Mark
	Li-ion batteries recycling symbol
	Manufacturer
	Date of Manufacture
	Reference Number
	Serial Number
	Use only at altitude 2000 meters or lower
	Use only in non-tropical climate conditions

Warranty

See Laerdal Global Warranty for terms and conditions. For more information visit laerdal.com.

Federal Communications Commission (FCC) and Industry Canada (IC) Statements

This device complies with part 15 of the FCC Rules and Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. this device must accept any interference received, including interference that may cause undesired operation.

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This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

The term "IC" before the equipment certification number only signifies that the Industry Canada technical specifications were met.

Caution: Changes or modifications not expressly approved by Laerdal Medical could void the user's authority to operate the equipment.

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Radiation Exposure Statement

This equipment complies with IC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance 20 cm between the radiator & your body.

EU

This product is in compliance with the essential requirements of Council Directive 2014/53/EU on Radio Equipment (RED).

The product is in compliance with Council Directive 2011/65/EU on restriction of the use of certain hazardous substances (RoHS).

The product is in compliance with Council Directive 2006/42/EC on machinery.

Size and weight	
Dimensions (Patient Simulator)	1800 mm (L) x 500 mm (W) x 370 mm (stomach depth) (71 in x 20 in x 14.5 in)
Weight (Patient Simulator)	44 kg (97 lbs)
Patient Simulator power	
External power	Input voltage 20VDC, 5A
Internal battery	Li-ion 14.6V, 5.2Ah, 76Wh
Air & CO ₂	
Internal air tank	0.9 bar
Environmental conditions	
Operating temperatures	5 °C to 40 °C (43 °F to 104 °F)
Storage and transport temperatures	-15 °C to 50 °C (5 °F to 122 °F)
Battery charging temperatures	19°C to 25 °C (66 °F to 77 °F)
Relative humidity	20% -90% (non-condensing humidity)
Not tested with salt spray	

Material chart for Patient Simulator	
Clothes	Cotton, Polyester
Airways	PVC (DEHP free)
Skins	Silicone
External hard plastics	ABS, PBT
Inner plastics	ABS, POM, PBT, PP, PC/ABS, PA, COP, LSR, Silicone
Metal components	Aluminum, Brass, Steel
Cleaning fluids	
Cleaning fluids (blood and fluid system)	1:1 mixture of distilled water with 70% isopropanol
Cleaning surfaces	Soap and water
Simulated IV fluids	
Simulated IV fluids	Distilled water
Blood	Distilled water
Urine	Distilled water

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© 2025 Laerdal Medical AS. All rights reserved.
MamaAnne® and Laerdal® are registered trademarks of Laerdal Medical AS.