

Little Anne

Advanced Airway Trainer

Important Product Information



Read these instructions thoroughly. Observe all warnings, cautions and instructions in the User Guide and in this Important Product Information booklet. Retain this booklet for future reference.

Warnings and Cautions

A Warning states a condition, hazard, or unsafe practice that can result in serious personal injury or death.

A Caution states a condition, hazard, or unsafe practice that can result in minor personal injury or damage to the product.

Notes

A note states important information about the product or its operation.

General guidelines

Warnings

- For training use only - May be operated by humans, but not for clinical use and not to be applied to any person (including actors or standardized patients).
- Do not perform mouth-to-mouth ventilation on airway management manikins that are not equipped with an airway filter.

Caution

- Manikin electronics have no IP rating. Do not expose any electronic components to liquids of any form.

Care, Maintenance and Cleaning

Cautions

- Clean the manikin and all accessories immediately after training and dry thoroughly. Failure to do so may lead to bacterial growth, mildew, or product damage.
- Do not submerge any parts.

Notes

- These can discolor the manikin:
 - Pigments from lipstick and pens
 - Latex gloves

- Using clothes not supplied by Laerdal
- Simulated blood
- Using bleach on external manikin surfaces

- Friction between the leaf spring and torso may cause small fragments to appear in and under the manikin. This is part of normal function and does not damage the product.

Follow the user guide instructions to keep the manikin clean and in a hygienic condition. Change the lung if there are visible signs of wear or damage to lung bag or metal patch.

Battery

Warnings

- Batteries must be replaced with correct type
- Do not use batteries with a nominal cell voltage higher than 1.5 V.
- Remove batteries if product is not to be used for an extended period

Notes

- If the device is used with rechargeable batteries, the "Low battery warning" deviates from normal behavior
- Internal charging is not supported

Following battery types are recommended:

- Non-rechargeable alkaline batteries, LR6 1.5 V AA
- Rechargeable NiMH batteries, HR6 1.2 V AA

Accessories for LAAS

Warning

- Stifneck Select Collar and Thomas Tube Holder are provided as accessories for training with the LAAS kit. Do not use in real patient care or on a standardized patient. Cleaning procedures are not included.

Federal Communications Commission (FCC)

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Changes or modifications not expressly approved by Laerdal Medical could void the user's authority to operate the equipment.

Little Anne contains FCC ID: QHQ-LA

Waste Handling

Recycle and dispose of the product in accordance with your country's recommendations.

By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

For more detailed information about treatment, recovery and recycling of this product, contact your local waste disposal service or local governmental office for waste management.

Specifications

| Size and Weight | |
|-----------------------|---|
| Dimensions | 64 cm x 37 cm x 21,5 cm (25.2" x 14.5" x 8.5") |
| Weight | 3 kg (± 0.2 kg) |
| Environment | |
| Operating temperature | 5 °C to 35 °C (41 °F to 95 °F) Degradation in product performance may be experienced when used at low temperatures (below 15 °C / 59 °F) |
| Storage temperature | -10 °C to 50 °C (14 °F to 122 °F) |
| Relative humidity | 5 % - 95 % relative humidity |

Warranty

Refer to the Laerdal Global Warranty for terms and conditions. For more information visit www.laerdal.com

© 2026 Laerdal Medical AS. All rights reserved.



Laerdal Medical AS,
P.O. Box 377
Tanke Svilandsgate 30,
4002 Stavanger, Norway
T: (+47) 51 51 17 00

20-23812 Rev C

laerdal.com



Laerdal
helping save lives