
MOYO

Fetal Heart Rate Monitor

User Guide

CONTENTS

- Clinical indications**..... 4
- Important information**..... 6
- Overview** 8
 - Device* 8
 - Accessories*..... 9
- Charging and setup**..... 11
 - Charging* 11
 - Set the device language, date and time*..... 12
- User interface**..... 13
 - Detailed*..... 13
 - Basic*..... 13
 - History*..... 13
 - Signal quality* 16
 - Low battery*..... 17
 - Menu*..... 18
- Clinical use** 19
 - Best practice for obtaining FHR*..... 19
 - Causes of unclear or poor signal* 19
 - Obtaining FHR*..... 20
 - Comparing FHR with MHR*..... 22
- Alarms** 24
 - Overview*..... 24
 - Physiological alarms* 26
 - Technical alarms* 28
 - Alarm history* 31

Reprocessing	34
<i>Power off</i>	34
<i>Cleaning and disinfection</i>	35
<i>Inspection</i>	37
<i>Maintenance</i>	37
<i>Waste handling</i>	37
Connectivity	38
Troubleshooting	39
Regulatory	41
<i>EU</i>	41
<i>WEEE statement</i>	42
<i>Warranty</i>	42
Specifications	43
Electromagnetic conformity	45
<i>Electromagnetic emissions tests</i>	45
<i>Electromagnetic immunity tests</i>	46

CLINICAL INDICATIONS

Device description

Moyo is a fetal heart rate (FHR) monitor using Doppler ultrasound technology.

Indication for use

Moyo is indicated to be used during labor to monitor the fetal heart rate.

Intended use

Moyo is intended to monitor the fetal heart rate.

Intended users

Moyo is intended to be used by healthcare professionals specialized in labor and delivery.

Intended use environment

Obstetric healthcare facilities.

Clinical benefit

Positive effect on obstetric management by monitoring the fetal heart rate.

Clinical outcome

Desired outcome of fetal heart rate monitoring is recognition of non-reassuring fetal status to enable appropriate obstetric interventions.

Known side effects

Continuous fetal heart rate monitoring may increase the rate of cesarean sections and operative vaginal deliveries.

NOTE

This side effect is not specific to Moyo.

Contraindications

None known.

Limitations of the device

Not suitable for use on pregnancies with multiple fetuses.

Moyo can be used to measure fetal heart rate either intermittently or continuously.

IMPORTANT INFORMATION

Before using Moyo, read these instructions completely to become familiar with the operation and maintenance of Moyo. Observe all warnings, cautions and instructions in this User Guide.

WARNING

– states a condition, hazard, or unsafe practice that can result in serious personal injury or death.

CAUTION

– states a condition, hazard, or unsafe practice that can result in minor personal injury or damage to the product.

NOTE

– states important information about the product or its operation.

 **WARNINGS**

- Read this User Guide and become familiar with the operation of the device prior to use.
- This product is intended for use by authorized healthcare professionals only.
- Clean and disinfect Moyo and relevant accessories as described in the Reprocessing section before first clinical use.
- Use of Moyo does not replace clinical supervision.
- Do not use Moyo to monitor FHR during a multiple birth. The device is intended for monitoring a single fetus only. In case of monitoring of multiple fetuses, alternative devices or methods should be used.
- Ensure that the ultrasound transducer is not mistakenly picking up maternal heart rate when making obstetric decisions. This is especially relevant when you see a low HR from the ultrasound transducer.
- Turn off Moyo and clean and disinfect between patient uses.
- Additional assessments should be done to confirm the need for C-section.

 **CAUTIONS**


- Do not submerge the Moyo main device or ultrasound transducer in liquid.
- Due to the magnet in the loudspeakers, metal objects like staples can be pulled in and damage them.
- Do not change or modify this equipment without authorization of the manufacturer.

 **NOTES**


- Charge Moyo fully before the first clinical use.
- In the case of device problems, continue following your local procedure with an alternative monitoring device or method.
- Should any serious malfunction, undesirable incident with, or deterioration in the functionality or performance of the device occur, contact Laerdal promptly. The competent authority where the incident took place and/or the device was used should also be notified.

OVERVIEW


Device Front

Back


Accessories

Power supply

Adapters


2x single use transabdominal belts


Charging cable

Transabdominal belt

A transabdominal belt can be used to secure the ultrasound transducer to the mother's abdomen.

Neck strap


To support maternal mobility during fetal heart rate monitoring, Moyo includes an eyelet for secure neck strap attachment.

CHARGING AND SETUP

Charging

Moyo has an internal rechargeable battery, which should be charged between patient uses or if the battery is low and/or empty. Use provided charger, model PSAI05R-050QL6, to charge the battery.

During charging, Moyo shows the battery status on the display.

CAUTION

Do not connect the charger to Moyo while in use on a patient. Moyo does not provide FHR whilst connected to a charger.


NOTES

- Charging time can be up to 4 hours when the battery is empty.
- The screen will dim down to preserve battery.

Set the device language, date and time

When turning on Moyo for the first time, the language, date and time need to be set to reflect your location. Press and hold the **power button** to turn on Moyo.

NOTE

Language, date and time can later be changed by accessing the menu. See **Menu section** for more information.

USER INTERFACE

Moyo has three main views: Basic, Detailed and History. Use the *navigation buttons* to change between the views.

Basic


Detailed

History


Basic view

NOTES


- At least 1 minute of good signal quality is required to first calculate the baseline number.
- Press and hold the right **navigation button** to revert back to the Detailed view.
- When an alarm is triggered, Moyo will automatically revert to the Detailed view.

Detailed view

History view

Moyo can display up to 120 minutes of FHR history (when available) in 15-minute segments. Use the *navigation buttons* to navigate between the different segments.


Signal quality

For Moyo to provide reliable FHR information, it is important to have a good Doppler ultrasound signal. Moyo indicates the quality of the signal as shown below:


Good signal


Graph:
White trace


Unclear signal

Graph:
Grey trace

Poor signal

Graph:
No trace

NOTE

The baseline number is not available after a period of poor signal.

Adjusting the Doppler volume


The Doppler volume from Moyo can be adjusted or muted.


NOTE

The alarm volume will not be affected by the Doppler volume adjustments.

Low battery

Low battery at **start up** (<20%)

Low battery **during use** (<10%)

NOTE

Charge Moyo regularly to ensure that it is always ready for use.

Menu

You can access the Alarm log, Episode log, System log and device specific settings of the device in the menu.

Press the **volume button** for 2 seconds to access the menu of Moyo.

Use the **navigation buttons** to navigate through the menu. To enter one of the submenus, press the **checkmark button**. Press the **volume button** to go back to previous page or back to the FHR monitoring displays.

Menu overview

ALARM LOG	List of alarms from the ongoing sessions.
EPISODE LOG	List of episodes stored on this device.
SYSTEM LOG	List of technical errors on this device.
SETTINGS	
Date	Adjust the date on the device.
Time	Adjust the time on the device.
Language	Adjust the language of the device.
Default Doppler volume	Set the default Doppler volume of the device.
Default screen mode	Set the default screen mode of the device.
Power saving	Set when the screen of the device should dim down.
ABOUT	Information on device firmware, usage, etc.

CLINICAL USE

Best practice for obtaining FHR

- Palpate to identify fetal position.
- Ensure there is gel between the skin and the ultrasound transducer.
- Place the ultrasound transducer over the fetal back.
- Listen for a clear rhythmic Doppler sound (galloping horse).
- Check that the displayed FHR matches the rhythmic Doppler sound.
- Confirm that the ultrasound transducer is not picking up the maternal heart rate (MHR).

Causes of unclear or poor signal



- Dislocated or poor placement of the ultrasound transducer:
- Too little gel.
- Too loose abdominal transducer belt.
- Movement by the mother or fetus.
- Measuring on multiples.
- Unintentionally measuring mother's aorta or other large vessels.

CAUTIONS

- *If the device displays unreliable FHR, ensure that there is enough gel and/or reposition the ultrasound transducer. If FHR is still not reliable or found, continue monitoring according to your medical procedures using alternative devices or methods to obtain FHR.*
- *It may be difficult to obtain a reliable FHR in cases of maternal obesity, polyhydramnios or oligohydramnios.*


Obtaining FHR

1. Palpate the mother's abdomen to identify the fetal position.
2. Apply ultrasound gel on the ultrasound transducer.

3. Press and hold the **power button** to turn on Moyo.

4. Place the ultrasound transducer over the fetal back and listen for a clear rhythmic Doppler sound (galloping horse).
5. Look for a clear FHR number on Moyo. Check that the displayed FHR matches the rhythmic Doppler sound coming from Moyo.
6. Confirm that the ultrasound transducer is not picking up the maternal heart rate.

CAUTION



Only use water-based ultrasound gel.

NOTES

- Listen for the startup sound when turning Moyo on to ensure that the audio system is functioning.
- If the power-on self-test of Moyo passes, the status light lights up steady green. If it fails the status light lights up steady red. Refer to **Troubleshooting section**.


Comparing FHR with MHR

The ultrasound transducer may pick up MHR through signal sources such as the aorta or other large vessels. If you are not sure that the detected HR is from the fetus, instruct the mother to hold the electrode pads to obtain the MHR for comparison.


NOTE





- If having difficulty obtaining a reliable MHR using Moyo (often due to too dry skin), add a drop of water to the fingertips before touching the electrode pads.
- The electrode pads could heat up slightly during operation.

1. **Measuring MHR**
Ensure that the mother is holding the maternal electrodes.

2. **Analyzing heart rates**
Moyo is comparing MHR with HR from ultrasound transducer.

3. **Different results that can appear:**

Heart rate from fetus
Transducer position is good, keep monitoring.

Heart rate from mother
Reposition transducer to achieve a better reading of FHR.

Not able to conclude
Try again or use alternative means.

ALARMS

Overview

When an alarm is triggered, it will be visible at the top of the display accompanied by the alarm sound.

Alarm-bar
Indicate which alarm was triggered and whether the alarm has been acknowledged.

Status light
Changes color when an alarm is triggered.

Legend:

- Alarm present
- Alarm build-up

For detailed information, go to page 31.

Acknowledge alarms
All alarms can be acknowledged (alarm sound muted) by pressing the *checkmark button*. For some alarms, tips will be provided when the alarm is acknowledged.

Moyo is designed to monitor FHR and provide visual and auditory alarms. Moyo provides the following alarms:

ALARM	ALARM CONDITION CATEGORY	MOYO PRIORITY	IEC 60601-1-8 PRIORITY
Critical technical error	Technical	1	Medium
Empty battery	Technical	2	Medium
Poor signal	Technical	3	Low
Maternal HR detected	Technical	4	Low
Bradycardic baseline	Physiological	5	Medium
Tachycardic baseline	Physiological	6	Medium
Low battery	Technical	7	Low

See [Alarm specification section](#) for detailed alarm description including alarm trigger and reset conditions.

NOTES

- Moyo's alarm system is not active until an FHR baseline value is displayed for the first time. (Note that the critical technical error and battery alarms are always active.)
- The alarms will end when the underlying alarm-condition is resolved.
- When an alarm is triggered, Moyo will automatically revert to the Detailed view.
- When an alarm is triggered, Moyo will automatically set the Doppler sound to the default volume.
- The alarm acknowledgment is a "Timed acknowledgement". If the underlying condition is not resolved within 30 minutes after the acknowledgment, the alarm sound will resume.

Physiological alarms


NOTE

Moyo is only able to trigger a physiological alarm for the conditions specified in this user guide. Other non-reassuring FHR indications can be identified by observing the FHR numbers and graph carefully.

Bradycardic baseline

The Bradycardic baseline alarm is triggered if the FHR baseline has been less than 100 bpm for at least 3 minutes or less than 110 bpm for at least 10 minutes.

When the alarm is acknowledged, Moyo will remind you to compare the MHR with the heart rate from the ultrasound transducer to ensure that the heart rate being picked up is from the fetus and not the mother:


Once the alarm is acknowledged, Moyo will remind you to compare the maternal heart rates.

Press **checkmark button** to reopen the tips screen.

Tachycardic baseline

The Tachycardic Baseline alarm is triggered if the FHR baseline has been more than 180 bpm for at least 3 minutes or more than 160 bpm for at least 10 minutes.


Technical alarms

Poor signal

For Moyo to be able to provide reliable FHR information, it needs to get a good signal from the ultrasound transducer. The Poor signal alarm is triggered if the signal picked up by the ultrasound transducer:

- is poor for more than 2 mins.
- is unclear for more than 5 mins.


Once the alarm is acknowledged, Moyo provides some tips that can improve the signal quality.

Press *checkmark button* to reopen the tips screen.


Maternal HR detected

The “Maternal HR detected” alarm is triggered when there is a suspicion that the MHR is being picked up by the ultrasound transducer. Compare the MHR with the heart rate coming from the ultrasound transducer to identify the source of the displayed FHR number.

Once the alarm is acknowledged, Moyo will remind you to compare the MHR with the heart rate from the ultrasound transducer.

Press *checkmark button* to reopen the tips screen.

If Moyo concludes that the heart rate detected is maternal, it reminds you to “Reposition and listen for the FHR sound”.

Critical technical error

Moyo cannot be used due to a technical error. Refer to the [Troubleshooting section](#) to identify the cause of the issue.

Empty battery

Moyo cannot be used due to empty battery.

NOTE


In case of empty battery or critical technical error, Moyo will turn off automatically after 2 minutes.

Alarm history

Moyo stores and displays all alarms that have been triggered during the ongoing episode. These historical alarms can be found in:

1. Alarms history

Alarms are displayed above the graph. When Moyo triggers an alarm, it displays a full line starting at the exact time the alarm was triggered, up until the alarm stops. The lines are colored according to the alarm color:


Additionally, Moyo indicates when it started to notice potential alarm-triggering conditions. This alarm “build-up” is indicated by dotted lines.

2. Alarm log

Moyo also provides an alarm log that can be found in the [Menu](#). The alarm log will only hold the 50 latest entries. A new entry will overwrite the oldest. The content of the alarm log is lost if Moyo is turned off or somehow loses power.

Alarm specification

ALARM	CATEGORY (alarm condition)	MOYO PRIORITY	IEC 60601-1-8 PRIORITY	AUDIO SIGNAL
Critical technical error	Technical	1	Medium	C-C-C (fixed pitch)
Empty battery	Technical	2	Medium	C-C-C (fixed pitch)
Poor signal	Technical	3	Low	E-C ("ding dong")
Maternal HR detected	Technical	4	Low	E-C ("ding dong")
Bradycardic baseline	Physiological	6	Medium	C-E-G (major chord)
Tachycardic baseline	Physiological	7	Medium	C-E-G (major chord)
Low battery	Technical	8	Low	E-C ("ding dong")

NOTE

All alarms, except for the "Low battery" alarm, are mutually exclusive. Only the alarm with the lowest Moyo priority number will trigger and be shown in case multiple alarm conditions exist.

How and when to verify functionality of the alarm system:


When you turn Moyo on, verify that you hear the startup sound, and see that the display turns on as normal without any error message.

ALARM TRIGGER CONDITION	RESET CONDITION	PURPOSE
Technical failure detected (various conditions, e.g. no communication with probe).	N/A	To inform Moyo can no longer be used.
Battery status <5%	N/A	To inform Moyo can no longer be used.
2 min of no FHR or poor quality, or 5 min of unclear quality.	15 s of good signal quality.	To attract attention of users to reposition ultrasound transducer to improve signal quality.
Sudden FHR drop ≥ 30 bpm from baseline lasting ≥ 1 min OR ECG HR is measured with conclusion that Doppler HR is maternal.	15 s of FHR near baseline again OR Measured ECG HR does not match Doppler HR.	To attract attention of users to check maternal heart rate (if confirmed, reposition ultrasound transducer).
Current FHR value <110 bpm AND (The FHR baseline has been <100 bpm for >3 min OR The FHR baseline has been <110 bpm for >10 min)	5 s of baseline FHR within 110-160 bpm	To attract attention of users to observe FHR patterns and consider obstetric actions.
Current FHR value is >160 bpm AND (The FHR baseline has been >180 bpm for >3 min OR The FHR baseline has been >160 bpm for >10 min)	5 s of baseline FHR within 110-160 bpm	To attract attention of users to observe FHR patterns and consider obstetric actions.
Battery status <8%	N/A	To attract attention of users to inform Moyo will soon run out of battery.


REPROCESSING

Power off

Moyo needs to be turned off between each patient use to ensure that FHR history from one patient is not used to make obstetric decisions on the next patient.

1. Press and hold the **power button** for 2 seconds to turn off Moyo

2. Confirm turning Moyo off by pressing the **power button** again.

WARNING

Turn off Moyo and clean and disinfect between patient uses.

NOTE

Moyo will turn off automatically after 5 minutes of inactivity.


Cleaning and disinfection

Clean and disinfect Moyo (main unit, cable and ultrasound transducer) after each patient use to minimize the risk of cross-contamination. Turn off Moyo before starting the cleaning process.


WARNINGS

- Effective disinfection is not possible without first performing a thorough cleaning.
- Care should be taken while handling the product between cleaning and disinfection.

Manual cleaning procedure


1. Clean all soiled surfaces of Moyo main unit, cable and ultrasound transducer using a clean cloth soaked with lukewarm (30 - 40 °C) tap water with mild dishwashing detergent. Ensure all surfaces are wiped a minimum of three times. To remove difficult soil, use a soft bristled brush (e.g. toothbrush) for at least 30 seconds or until visibly free of soil.
2. Rinse under running tap water for 10 seconds.

3. Dry using a clean cloth or by air drying.
4. Visually inspect the device. If any soil is observed, repeat cleaning steps until no visible soil remains.

Disinfecting

1. Disinfect all surfaces of Moyo main unit, cable and ultrasound transducer by wiping with a clean, non-linting cloth with 70% ethanol. Focus on the ultrasound transducer; main unit control panel and maternal heart rate electrodes (front and back).
2. Ensure all surfaces remain visibly wet for a minimum of 12 minutes, re-apply disinfectant as needed.
3. Allow to air dry.

CAUTION

Do not use sodium hypochlorite (bleach) or other chlorine-releasing chemicals to disinfect as this can damage the product.

Inspection


After each use, inspect Moyo's main unit, ultrasound transducer and cable for damage or cracks that could lead to ingress of liquid. If there is any damage, refer to [Troubleshooting section](#).

After you confirm that Moyo has no visual damage, turn on Moyo to verify that you hear the startup sound, and see that the display turns on as normal without any error message. Refer to [Troubleshooting section](#) if there is any error.

CAUTION

Excessive bending of the cable can damage it.

Maintenance

Moyo does not require periodic maintenance or calibration other than noted in the [Inspection section](#). Moyo is not repairable by the user. Refer repair to personnel qualified by Laerdal, or to Laerdal or one of its authorized distributors.

Waste handling

Moyo contains electronic components. Dispose of it at an appropriate recycling facility in accordance with local regulations.

CONNECTIVITY

Moyo allows for secure wireless connectivity using Bluetooth® Low Energy and proprietary communication services. The user of Moyo, or the responsible organization, is not able to, nor intended to use the wireless connectivity functionality of Moyo.

Wireless functionality such as transfer of stored data, troubleshooting commands and firmware upgrade are only available during charging, not during clinical use. These functions are currently limited to non-commercial, proprietary Laerdal applications.

TROUBLESHOOTING

CAUTION

In case of device problems, continue your medical procedure. Do not allow long pauses while troubleshooting.

SYMPTOM	POSSIBLE CAUSE	POSSIBLE SOLUTION
Critical technical error	<p>The error code is the sum of the below individual error codes:</p> <p>0x0002 – Configuration memory checksum</p> <p>0x0004 – Episode storage memory</p> <p>0x0008 – Voltage regulator</p> <p>0x0010 – Real time clock</p> <p>0x0020 – Display communication</p> <p>0x0040 – Ultrasound transducer communication</p> <p>0x0080 – Fuel gauge communication</p> <p>0x0100 – Bluetooth LE communication</p> <p>0x0200 – Audio amplifier communication</p> <p>0x0400 – Speakers</p> <p>0x0800 – Backup voltage regulator</p>	<ul style="list-style-type: none"> • Continue standard procedure without use of the device. • Use alternative means of measuring fetal heart rate, e.g. a Pinard fetoscope. • Restart Moyo, if the problem persists, the device is broken and should not be used. • Report the error by scanning the QR code displayed on the device and email the Laerdal support team.

SYMPTOM	POSSIBLE CAUSE	POSSIBLE SOLUTION
Fetal heart rate is not detected.	<ul style="list-style-type: none"> • Ultrasound transducer is misplaced or not on the mother. • Stimulation/movement/handling of patient temporarily generating too much noise. • Heart rate below 50 bpm. 	<ul style="list-style-type: none"> • If persistent: <ul style="list-style-type: none"> - Palpate to find fetal back and place ultrasound transducer. - Check if enough gel is used. - Check if transabdominal belt is tight. • If the problem continues, use alternative means of measuring fetal heart rate, e.g. a Pinard fetoscope.
Device does not start up or there is power loss during use.	<ul style="list-style-type: none"> • Battery depleted • Power button broken • Water ingress 	<ul style="list-style-type: none"> • Continue standard procedure without use of the device. • After the standard procedure: connect the charger cable to the device, and charge. • If the problem persists, the device is broken and should not be used. • If under warranty, contact Laerdal Global Health.
Nonfunctional or damaged part detected during equipment inspection.		If under warranty, contact Laerdal Global Health.
Not able to measure maternal heart rate.	<ul style="list-style-type: none"> • Skin too dry • Dirty maternal heart rate electrodes 	<ul style="list-style-type: none"> • Add a drop of water to the fingers in contact with the maternal heart rate electrodes. • Clean and rub the maternal heart rate electrodes using ethanol. • Keep fingers still while measuring maternal heart rate.

REGULATORY

SYMBOL GLOSSARY	
	CE-mark with notified body
	CE-mark
	Manufacturer
	Waste Electrical and Electronic Equipment symbol
	Meets IEC type BF applied part leakage current requirements.
IP55	Protected against ingress of dust and jets of water from all directions.
	Follow instructions for use.
	Medical device
	GTIN (Global Trade Item Number) + Serial Number
	Model number
	Temperature limitation
	Atmospheric pressure limitation
	Humidity limitation
	Single patient multiple use

EU

This medical device is in compliance with the general safety and performance requirements of Regulation (EU) 2017/745 on medical devices. This product is in compliance with Council Directive 2014/53/EU on Radio Equipment (RED) and Council Directive 2011/65/EU on restriction of the use of certain hazardous substances (RoHS).

WEEE statement

European Directive 2012/19/EU (WEEE)

This appliance is marked according to the European directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE). By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

Warranty

Refer to the Laerdal Global Warranty for terms and conditions.
For more information visit laerdal.com

SPECIFICATIONS

DIMENSIONS	
Size	100 × 107 × 26 mm (W × H × D, without cable and ultrasound transducer)
Weight	285 g (main unit and ultrasound transducer)
ENVIRONMENTAL	
Operating temperature	5 °C to 40 °C
Short term (24 h) storage/ shipping temperature	-20 °C to 60 °C
Long term storage temperature	15 °C to 25 °C (room temperature)
Atmospheric pressure	700 hPa – 1060 hPa (up to 3000 meters above sea level)
Humidity	15% to 90% relative humidity, non-condensing, but with a water vapor partial pressure ≤50 hPa
Ingress protection	IP55 (both main unit and ultrasound transducer)
MATERIALS	
Main unit	ABS, PC, conductive TPU
Ultrasound transducer	ABS
Transducer cable	TPU
 NOTE	<i>Moyo is not made with natural rubber latex or PVC.</i>
MEASUREMENTS	
Heart rate measurement - ultrasound transducer	Accuracy: short term average ± 2 bpm in the range 50 – 250 bpm. No detectable heart rate, or heart rate outside the range, is displayed as “-?-”
Heart rate measurement - maternal HR electrodes	Accuracy: short term average ± 1 bpm or $\pm 1\%$ bpm in the range 30 – 200 bpm and 0.2 – 5 mV QRS amplitude. No detectable heart rate, or heart rate outside the range, is displayed as “-?-”
Essential performance (EP)	The ability to measure fetal heart rate (FHR) to the specified accuracy when displayed, and with the allowance of occasional spurious FHR numbers (which does not present an unacceptable risk).
User expectation if EP is lost or degraded	The user will perceive unexpected numbers and history tracings generally accompanied by disturbances in the Doppler sound.
Medium priority alarm sound pressure	75 dBA (C-E-G) and 68 dBA (C-C-C) at 1 m

DATA STORAGE	
32 MB internal storage (up to 100 episodes or ca. 600 hours of data).	
Data is stored at 2 Hz and includes heart rates and signal qualities.	
POWER	
Battery	Internal rechargeable battery, Lithium-ion, 3.7 V, 2300 mAh
 NOTES <ul style="list-style-type: none"> Discharged batteries or batteries stored outside room temperature degrade faster. For long-term storage of Moyo, charge the device to at least 60% and recharge every other year. 	
Run time	>10 hours (full charge on new battery) >7 hours (full charge at expected end of battery service life depending on use and if properly maintained)
Charge time	Up to 4 hours (full charge of empty battery)
 NOTE Device cannot be used clinically while charging.	
Power supply (charger)	Input 100-240 V AC, 50-60 Hz, 0.3 A. Output 5 V DC, 1 A. USB connector.
 CAUTION Only charge Moyo using the provided power supply, Pihong PSAI05R-050QL6, or an alternative 5 W USB power supply, 5 V DC, 1 A, that is IEC 60601-1 or IEC 62368-1 certified.	
RADIO	
Type	Bluetooth® Low Energy transmitter
Frequency band	2.400 – 2.4835 GHz
Modulation	Gaussian frequency shift modulation
Maximum radio-frequency power transmitted	1 mW
Effective radiated power	0 dBm
FCC ID/IC ID	Contains FCC ID: X8WBT840F, IC: 4100A-BT840F
ULTRASOUND TRANSDUCER	
Mode	Nine-crystal Pulsed Wave Doppler Ultrasound
Frequency	1 MHz
Thermal indices (TI) and the mechanical index (MI) are always below 0.7.	
LIFETIME PARAMETERS	
Shelf life	2 years
Expected service life	3 years. Tested to 100 reprocessing cycles.

ELECTROMAGNETIC CONFORMITY

Moyo is intended for use in the following environments: health care facilities except for near HF surgical equipment and the RF shielded room for magnetic resonance imaging.

No particular actions are required to maintain safety and performance with regard to electromagnetic disturbances for the expected service life.

List of cables and battery charger with which Moyo is in compliance with the IEC 60601-1-2 EMC standard:

Power supply PSAI05R-050QL6 and USB cable, type A to type C, 1.8 m

CAUTIONS

- *Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.*
- *Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.*
- *Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device. Otherwise, degradation of the performance of this equipment could result.*

Electromagnetic emissions tests

EMISSIONS TEST	STANDARD OR TEST METHOD	COMPLIANCE
RF emissions	CISPR 11	Group 1 Class B
Harmonic emissions	IEC 61000-3-2	Class A
Voltage fluctuations/ flicker emissions	IEC 61000-3-3	Complies

Electromagnetic immunity tests

IMMUNITY TEST	STANDARD OR TEST METHOD	IMMUNITY TEST LEVEL
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Radiated RF EM fields	IEC 61000-4-3	10 V/m 80 MHz – 2.7 GHz 80% AM at 2 Hz in clinical mode and 1 kHz in charging mode
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	380-390 MHz: 27 V/m 430-470 MHz: 28 V/m 704-787 MHz: 9 V/m 800-960 MHz: 28 V/m 1700-1990 MHz: 28 V/m 2400-2470 MHz: 28 V/m 5100-5800 MHz: 9 V/m
Rated power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz
Proximity magnetic fields	IEC 61000-4-39	30 kHz: 8 A/m 134.2 kHz: 65 A/m 13.56 MHz: 7.5 A/m
Electrical fast transients / bursts	IEC 61000-4-4	2 kV 100 kHz repetition frequency
Surges: Line-to-line	IEC 61000-4-5	± 0.5 kV, ± 1 kV
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V; 0.15 MHz – 80 MHz 6 V in ISM and Amateur bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage dips	IEC 61000-4-11	0% of UT for 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° phase angles 0% of UT for 1 cycle and 70% of UT for 25/30 cycles Single phase: at 0
Voltage interruptions	IEC 61000-4-11	0% UT for 250/300 cycles
UT is the a.c. mains voltage prior to application of the test level.		

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© 2026 Laerdal Medical AS. All rights reserved.
Laerdal® is a registered trademark of Laerdal Medical AS.