
Pediatric Pocket Mask


User Guide



Laerdal
helping save lives

www.laerdal.com

CE MD

ENGLISH

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Inspect and verify that this product is complete and without defects before use.

Intended Use

The Laerdal Pediatric Pocket Mask is designed for mouth-to-mask ventilation of a non-breathing child or infant. For use by trained personnel only.

1. Attach one-way valve to mask port by twisting it until seated firmly (figure 1).
2. Open patient's airway.
3. Apply mask over patient's mouth or mouth and nose. If applying over mouth only (as used for children) be sure to pinch nostrils to avoid air leakage (figure 2).
4. Seal as illustrated. Open airway. Blow slowly until chest rises (figure 3).
5. Remove your mouth. Allow patient to exhale.

Notes

- *If patient vomits, remove mask and clear patient's airway. Clear mask thoroughly of foreign material. Blow through valve to verify operation. Reapply mask to patient and continue ventilation.*
- *Should any serious malfunction, undesirable incident with, or deterioration in the functionality or performance of the device occur, contact Laerdal promptly. The competent authority where the incident took place and/or the device was used should also be notified.*

To disassemble, hold mask with one hand, gripping the one-way valve between your thumb and two fingers of the other hand. Separate by twisting valve off.

Patient use

Mask and accessories are single use and should be disposed of after use. Can be used on same patient during a 24-hour period. Laerdal is not responsible for any consequences of re-use.

Training use

Mask can be cleaned and reused only for future training. One-way valve must be disposed of after use.

Cleaning and Disinfection

Wash and scrub in warm soapy water. Rinse in clean water. Submerge for 10 min. in a 0,5% sodium hypochlorite solution. Rinse again and allow to air dry.

Temperature limits

Operating temperatures: -18°C (0,4°F) to 60°C (140°F).
Storage temperatures: -40°C (-40°F) to 70°C (158°F).
Relative humidity: 5% - 95% (non condensing).

Expiratory resistance for patient

Mask & one-way valve: Approximately 4.3 cm H2O
Inspiratory resistance for patient
Mask & one-way valve: Approximately. 3.3 cm H2O

Deadspace 90,5 ml

Performance data are nominal values only.

CE This medical device complies with the general safety and performance requirements of Regulation (EU) 2017/745 for medical devices.

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© 2025 Laerdal Medical AS. All rights reserved