EN

Laerdal BaXstrap Spineboard - User Guide
REF 982500, 982599, 982600 and 982699

Brand Name:	Laerdal BaXstrap Spineboard
Common Name:	Backboard/Spineboard/Longboard
Description:	Full-body length spinal immobilization device

Device description: The Baxstrap spineboard is a rigid, lightweight full-body length spinal immobilization device.

Indication for use: Used for patients with potential spinal injuries to limit movement and likelihood of further injury during transport.

Intended use: A base structure to be used with other adjunct cervical spine and head immobilization devices to facilitate in-line, neutral immobilization and transport of adult and pediatric patients.

Intended users: Healthcare professionals, including paramedics, emergency room staff, and first responders.

Intended use environment: The BaXstrap can be used in any environment, both indoors and outdoors, including water.

Know side effects: Discomfort during transport, potential skin irritation or pressure sores if not monitored properly.

Contraindications: No known

Notes
States important information about the product or its operation.
- *Always follow local protocol for patient movement and immobilization techniques.*
- *Not made with natural rubber latex.*
- *The product will not disturb the images of the spinal column when using X-ray, CT, and MRI.*
- *Should any serious malfunction, undesirable incident with, or deterioration in the functionality or performance of the device occur, contact Laerdal promptly. The competent authority where the incident took place and/or the device was used should also be notified.*

Warnings
States a condition, hazard or unsafe practice that can result in serious personal injury or death.
- *Inspect the Baxstrap prior to initial use and before each patient use.*
- *Pay particular attention to the underside of the head end.*
- *If the product develops punctures, broken pins or cracks then discontinue use immediately.*
- *Ensure that the patient is secured using the strap.*

Cautions
States a condition, hazard or unsafe practice that can result in minor personal injury or damage to the product.
- *Read and thoroughly understand these directions. These directions are not intended as, and do not in any way constitute medical treatment advice, instruction in or teaching of medical practice, or imply or convey any rights, medical direction or authority.*
- *This and other spinal immobilization adjunct devices should only be used by personnel who have received adequate training.*
- *Spine boards should be removed as soon as it is clinically indicated.*
- *Spineboards should only be used when indicated based on local protocol.*
- *The length of the non-adjustable strap can be shortened by wrapping it around the handle of the BaXstrap.*

Directions for Use
1. Immobilize the cervical spine using a proper Cervical Immobilization Device
2. Position the patient on the BaXstrap using appropriate number of personnel.
3. Properly fasten/strap the patient to the BaXstrap to protect against patient movement in any direction (up, down, left or right).
4. Secure the head to the BaXstrap using a proper Head Immobilization Device.
5. Assure in-line neutral immobilization of the head, neck, torso and pelvis.

Cleaning: Non-sterile. Withstands cleaning in soap and water (1:16) and 1:100 household bleach: water.

Technical Specifications	
Operating and storage temperature	-34 °C to +52 °C (-29 °F to +126 °F)
Operating and storage humidity	5-95% RH non-condensing.
Shelf-life	3 years
Expected service-life	5 years 300 cycles of use Note: This is based on an average patient size of 150 kg and the defined reprocessing method.
Maximum load	200 kg
Symbol Glossary	
CE	This medical device complies with the general safety and performance requirements of Regulation (EU) 2017/745 for medical devices.
MD	Medical Device

The product has a 5 year limited warranty. See the enclosed "Laerdal Global Warranty" for terms and conditions. The warranty is also available at www.laerdal.com.

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© 2025 Laerdal Medical AS. All rights reserved.
BaXstrap is a trademark of Laerdal Medical AS.