

Penguin
Newborn Suction
Reusable

© 2025 Laerdal Medical AS.
All rights reserved.

Date of issue: 2025-07

Laerdal Medical AS
Tanke Svilandsgate 30, P.O. Box 377,
4022 Stavanger, Norway
Tel. (+47) 51511700
Fax (+47) 51523557
E-mail: laerdal.norway@laerdal.no

Laerdal
helping save lives

REF Cat. No. 98600060 **EN**
QTY 1
LATEX Not made with natural rubber latex
Discard if damaged
Disinfect after use
MD Medical Device

CE

This medical device complies with the general safety and performance requirements of Regulation (EU) 2017/745 for medical devices.

Note
Should any serious malfunction, undesirable incident with, or deterioration in the functionality or performance of the device occur, contact Laerdal promptly. The competent authority where the incident took place and/or the device was used should also be notified.

Technical specifications:
Suction strength (typical):
100 mm Hg [136 cm H_2O]
Operating temperature:
0 °C (32 °F) to 50 °C (122 °F)
Storage temperature:
-20 °C (-4 °F) to 60 °C (140 °F)
Material: Silicone
Fluid capacity: 75ml

Protected under US Patent 8,696,648 and other international patents; US design patent D634,002 S and other international design registrations. Ownership and all rights reserved.

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EN

This Device is intended for the removal of fluids from nasal and oral cavities of newborns and infants. Also suitable for more general use as an irrigating syringe.

Penguin must be cleaned and high-level disinfected before first time use.

1 To test before use
Squeeze body and block the suction tip to verify there is no air leakage. Body should not expand until tip is unblocked.

2 To Use
a Squeeze and insert
b Release to suction
c Empty contents
d Repeat as required.
Suction mouth first, then nose.

3 To reduce the risk of cross-contamination, disinfect after use
Open head. Wash and brush in soapy water. Flush in clean water. Autoclave (136 °C (277 °F) and 2.0 kg/cm² for 10-20 min); OR boil for 10 min. Dry completely and keep clean.

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