


Thomas Tube Holder - Adult

User Guide

Size 3 - 6

< 8 years

℞ ONLY

REF 600-10000 600-30000

LATEX


6.5 - 21 mm


-34° C to +52° C
(-30° F to +124° F)


CE

MD

Laerdal
helping save lives

Indications for use

Thomas Tube Holder is designed to secure single and double-lumen airway tubes such as Endotracheal Tubes (ETs), Laryngeal Mask Airways (LMAs) and other Supraglottic Airway devices (SGAs), after insertion into the trachea or esophagus to reduce the risk of accidental extubation. It is intended for use on persons of at least 8 years old with a short-term need (less than 24 hours) for airway-tube device management.

Setup Instructions

- Secure the product as shown (Fig. 1 - 3).
- Ensure correct position of head strap to prevent it from sliding (Fig. A).
- Release the tube clamp as shown (Fig. B).
- If necessary, position removable hook tab on strap to accommodate smaller neck size (Fig. C).

Notes

Important information about the product or its operation.

- *The Thomas Tube Holder meets the following standards: EN 1789, ISO 13485, ISO 14971, ISO 10993-1.*
- *Non-sterile product*
- *Should any serious malfunction, undesirable incident with, or deterioration in the functionality or performance of the device occur, contact Laerdal promptly. The competent authority where the incident took place and/or the device was used should also be notified.*

⚠ Warnings

A Warning states a condition, hazard, or unsafe practice that can result in serious personal injury or death.

- *This device should only be used by trained personnel in accordance with the User Guide.*
- *Apply just enough pressure against the tube to hold securely in place. Overtightening of the screw may restrict airflow.*
- *Do not use on children younger than 8 years old.*
- *Designed for single patient use only. Do not re-use. Re-use will lead to increased risk of cross-contamination, degradation of performance and/or device malfunction. Laerdal is not responsible for any consequences of re-use.*
- *Prior to use, inspect all parts to ensure they are free of damage and working properly.*
- *The patient should be constantly monitored while this product is in use.*

⚠ Cautions

A Caution states a condition, hazard, or unsafe practice that can result in minor personal injury or damage to the Thomas Tube Holder.

- *The device should only be attached to the patient after insertion of an advanced airway device.*
- *Do not use the device on a patient if it has been used during training.*
- *Do not disinfect or sterilize the product.*
- *The device is for short term use only. Remove after 24 hours.*
- *Storage temperature: -34 °C to +52 °C (-30 °F to 124 °F).*
- *Operating temperature: -18 °C to +43 °C (-0 °F to 109°F).*

Symbol Glossary

CE	This medical device complies with the general safety and performance requirements of Regulation (EU) 2017/745 for medical devices.
Rx ONLY	Caution: Federal law (US) restricts this device to sale by or on the order of a physician or appropriately licensed practitioner.
	Single Patient Use. Do not re-use.
REF	Reference order and number
UDI	Unique Device Identification
MD	Medical Device
	Not for use on children under 8 years old
	Storage temperature (Maximum and Minimum)
	Size range for supraglottic devices
	Size range for tracheal tubes (Outer Diameter)
LATEX	Not made with natural rubber latex

Materials

Part	Material
Foam	Polyethylene foam - PE
Strap	Polyurethane foam - PU Polyester flannel - PET Polyamide - PA
Mouthpiece	Polypropylene - PP
Other parts	Polyoxymethylene - POM

The materials have not been tested for flammability.

Dispose of according to local protocol.

Thomas Tube Holder is protected by US patent # 8,360,063 and other international patents. Thomas and Laerdal are trademarks owned by Laerdal Medical AS. Ownership and all rights reserved.

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© 2025 Laerdal Medical AS. All rights reserved.