


# SimNewB

## Quick Setup Guide






This Quick Setup Guide describes how to assemble and setup your SimNewB patient simulator. For information about how to use SimNewB, please refer to LLEAP help files and SimNewB Directions for Use.

Visit [www.laerdal.com/downloads](http://www.laerdal.com/downloads) to download the latest Directions for Use and software.



# INCLUDED COMPONENTS


Simulator



BP Cuff



Manuals & Warranty



IO Mandrel Legs



BP Extension Tube



IO Fill Syringe



IO Mandrel Skins



Air/CO<sub>2</sub> Tube



Birth Certificate




Hat and Pants



Eye Set



Level Chart



Blanket



Meconium Aspiration Kit



Blood



Umbilical Cords



CO<sub>2</sub> Regulator



Umbilical Blood



Liquid Soap



Airway Lubricant



Baby Powder



REQUIRED


Instructor PC




or  
SimPad



Link Box




Power Adaptor



Ethernet Cable for PC



Compressor



RECOMMENDED


Patient Monitor



Headset



Ethernet Switch



SpO2 Probe




Ethernet Cable for  
Patient Monitor




1

# CHARGE BATTERIES


Instructor PC




SimPad



Patient Monitor




Link Box




# 2A

# WIRED SETUP




2B

WIRED SETUP WITH PATIENT MONITOR









2c


SIMPAD WIRELESS SETUP




Patient Monitor



SimPad









Compressor



Headset

Link Box



- BP Extension tube 
- Air/CO<sub>2</sub> tube 
- USB cable 
- Simulator cable 
- Ethernet cable 

© 2014 Laerdal Medical AS. All rights reserved.  
Device Manufacturer: Laerdal Medical , P.O. Box 38,  
226 FM 116, Gatesville, Texas 76528 USA

Rev B



[www.laerdal.com](http://www.laerdal.com)



**Laerdal**  
helping save lives