



LLEAP Frequently Asked Questions



What is Laerdal Learning Application- LLEAP?

Laerdal has merged six different software applications into one single application- LLEAP. LLEAP is designed for instructors and simulator operators who use any of Laerdal's PC operated simulator for training and education of personnel and students in health care with the help of simulation training.

LLEAP is based on the SimMan 3G software platform, but has also adapted elements from Laerdal's SimPad interface, like the event registration component and instructor pop up messages. The interface will adapt to show the relevant controls for whichever simulator you are connected, still keeping the same overall layout and making switching between simulators and standardized patients a smooth experience.

Who is eligible to implement LLEAP?

In 2015, customers with a PC operated Laerdal simulator are entitled to one free LLEAP license for each compatible simulator currently owned. Laerdal will provide a quote for any hardware or assistance you may need to transition from your current Instructor Application to LLEAP.

What simulators are compatible with LLEAP?

Currently, there are 11 different Laerdal simulators that can be controlled with LLEAP; in addition, LLEAP can be used with Standardized Patients and task trainers.

- SimMan® 3G
- SimMan® 3G Trauma
- SimMan® 3G Mystic
- SimMan® Essential
- SimMan® Essential Bleeding
- SimMan® Classic
- ALS Simulator
- SimMom®
- SimJunior®
- SimBaby™
- SimNewB®



What are the hardware requirements to operate LLEAP?

- Check your hardware here

What if I do not “upgrade” to LLEAP in 2015?

Laerdal will continue to support the current Simulator Instructor Applications until 2020, but new developments, features and scenarios will only be available for LLEAP.

Can I trial the LLEAP Software?

Yes! LLEAP Lite is the unlicensed, free version of LLEAP. It allows you to:

- Practice operating LLEAP on a virtual manikin
- Bench test pre-programmed scenarios- including SimStore scenarios and patient cases

Adding a LLEAP license provides you full access allowing you to connect to the simulator and Patient Monitor; record sessions with a webcam and save sessions with Session Viewer.

Do I lose any functionality if I upgrade to LLEAP?

No, you do not lose functionality. In fact, you gain features such as:

- Improved simulation home screen with structured overview of all possible configurations and settings
- Automatic updates- when connected to the internet you are prompted about software updates
- Multi language- LLEAP is available in 11 languages
- Faster, more intuitive event registration
- Instructor pop-up messages

