



**Laerdal**<sup>®</sup>  
helping save lives



## Shockingly Realistic

### ShockLink™ Defibrillation Training Tool

ShockLink is a compact solution for standardized CPR-D training that enables users to adapt their own AED or advanced defibrillator for use as a training tool in a realistic and controlled environment. ShockLink allows you to use adhesive training pads on your CPR manikins but also generates multiple shockable heart rhythms that are the foundation for a solid and efficient CPR-D training program.

Laerdal Medical proudly introduces ShockLink™, a new simulation tool for realistic and efficient CPR-D training.



# ShockLink™ Defibrillation Training Tool

## FEATURES

- Adult defibrillation (up to 360 J)
- Child/infant defibrillation (from 10 J)
- Measure defibrillation energy and pacing
- ECG library\*:
  - Sinus Rhythm
  - VT Slow
  - VT Fast
  - SVT
  - VF
  - PEA
  - Asystole
  - Bradycardia

\*When used with Remote Control



## Q CPR®

ShockLink integrates to your Q CPR® platform. When used with Resusci Anne® Q CPR® AED and SimPad® SkillReporter™, heart rhythms and shocks delivered are automatically captured in the screen for future analysis and debriefing.



## SPECIFICATIONS

- Auto power on/off
- Fifty hours of training with one AAA battery
- Maximum defibrillator voltage: +/- 5000V
- Maximum measurable defibrillator voltage: +/- 2000V
- Maximum energy in a shock: 360 J
- Minimum measurable shock energy: 10 J
- Size:
  - 24 cms x 18 cms x 5 cms (l x w x h)
  - 9.5 inch x 7.1 inch x 1.9 inch (l x w x h)
- Weight: 1 lbs.
- Compatible defibrillators:

**PHILIPS**

MRx

FR3

**ZOLL**

Zoll AED Plus®

Zoll R series®

Zoll X series®

\* Zoll® Defibrillator Connector Item Number: 05-10100 required



LIFEPAK® 12

LIFEPAK® 15

LIFEPAK® 20

LIFEPAK® 1000

\* Physio-Control® Quik-Combo™ Connector Item Number: 05-10000 required

## ORDERING INFORMATION

### MAIN PRODUCTS

**185-10050 ShockLink System**  
Includes: ShockLink Unit, Remote Control, Adhesive Training Pads (1 Set), AAA Battery and Directions For Use

### PARTS AND CONSUMABLES

198-80150 Training Pads  
185-21050 Battery Lid  
185-00450 Carrying Pouch  
185-00350 Remote Control



Visit [Laerdal.com/ShockLink](http://Laerdal.com/ShockLink) for a complete listing of parts, accessories, and consumables.

