



Laerdal


helping save lives




Exceptional training makes all the difference

Premature Anne

Premature Anne is a realistically proportioned 25 week preterm manikin developed in collaboration with the American Academy of Pediatrics (AAP). She is designed to facilitate simulation training in teams of healthcare professionals involved in resuscitation and care of preterm infants.



Premature Anne



Airway features

- ET tube insertion
- Right mainstem intubation
- Positive pressure ventilation
- Sellick maneuver
- NG/OG tube insertion
- Cyanosis
- Vocal sounds

Vascular access

- Umbilicus with venous and arterial access
- Internal fluid reservoir

Breathing features

- Chest rise and fall with mechanical ventilation
- Heart sounds
- Lung sounds

Dry pre-ported IV sites

- Right saphenous vein
- Dorsum of left hand
- Left antecubital fossa

Meeting neonatal simulation training requirements.

Anatomically accurate

Realistically proportioned 25 week preterm infant

Realistic airway functionality

Designed to train airway management including the placement of an ET tube

Realistic training experience

Premature Anne facilitates training of resuscitation skills for preterm infants born at the limits of viability

Supports delivery of AAP neonatal resuscitation program

The features of Premature Anne are aligned with the objectives of the NRP course curriculum. The Premature Anne pack consists of the manikin and the SimPad PLUS preloaded with 8 NRP scenarios* specifically for preterm infants.

Lightweight and portable

Premature Anne is designed to be tetherless and highly mobile for use in multiple clinical settings and in transport scenarios

Ease of use - Intuitive SimPad PLUS user interface operates the baby and easily connects to the simulated patient monitor

Catalogue Number	Description
295-55101	AAP Premature Anne Pack
295-00050	Premature Anne Manikin Only
208-30001	SimPad PLUS Premature Anne (US)
208-30003	SimPad PLUS Premature Anne (UK)
208-30005	SimPad PLUS Premature Anne (JA)
208-30015	SimPad PLUS Premature Anne (AUS)
208-30033	SimPad PLUS Premature Anne (IE)
295-50150	Premature Anne Software License
204-30001PP	SimPad PLUS 1 year Protection Plan
204-30002PP	SimPad PLUS 2 year Protection Plan

* Please refer to laerdal.com website for details of the individual scenarios.

