



Laerdal[®]
helping save lives



More Efficient CPR Training

Resusci Anne[®] QCPR

The Resusci Anne[®] QCPR manikin and technology enables healthcare providers to obtain competence to deliver quality CPR in accordance with the latest Guidelines. QCPR Technology is utilized to help both instructors and trainees more effectively monitor and review CPR performance, increasing the efficiency and value of training time.



Resusci Anne® QCPR

In addition to providing high-quality CPR training, the Resusci Anne® QCPR manikin allows for those skills and others to be integrated into team training. The available add-on options and accessories make Resusci Anne QCPR more than just a high-quality CPR training manikin. The manikin can be upgraded to train teams on defibrillator use, airway management, IV administration, and rescue/trauma situations.

Key Training Features

CPR

- Quality CPR feedback (with feedback device)
- Mouth-to-mouth ventilations
- Ventilations with face shield, pocket mask, and BVM

Airway

- Realistic chest rise/fall
- Sellick maneuver
- Positive pressure ventilation
- Airway occlusion
- Supraglottic airway (with Airway Head)

Blood Pressure / pulses

- Manual Carotid pulses

Defibrillation

- Train using a live defibrillator (with ShockLink™)
- Link Technology to detect pad placement (with AED Skin)

General

- Interchangeable pupils (with Airway Head)
- IV initiation and fluid administration and IV drug administration (with IV Arm)
- Rescue/Trauma Simulation (with Trauma Limbs or Rescue Limbs)
- Latex free
- Wired and wireless operation

QCPR Feedback

- Compression rate and depth
- Complete release
- Correct hand position
- Appropriate ventilation volume
- Amount of compressions and ventilations
- Frequency and length of interruptions (not available with SkillGuide)

Operation



206-30001 SimPad® PLUS with SkillReporter™
206-50050 SkillReporter Software License



170-30050
SkillGuide™



202-56050
SkillReporter (PC)

Optional Add-ons



312025 /
150-10150
IV Arm



090050
Bleeding
Control Legs



312050
First Aid/
Trauma Limbs



312000
Rescue Limbs

Starting Configurations

Resusci Anne® QCPR

171-00150	Torso
171-01250	Full Body
172-00150	w/ Airway Head – Torso
172-01250	w/ Airway Head – Full Body
173-00150	AED - Torso
173-01250	AED - Full Body
174-00150	AED w/ Airway Head - Torso
174-01250	AED w/ Airway Head - Full Body

Kits

170-50750	Battery Upgrade Kit
173-50750	AED Skin Upgrade Kit
171-20050	Periodic Maintenance Kit
316000	Airway Trainer Update Kit

Warranty

One-year limited warranty. Refer to Laerdal Global Warranty for terms and conditions.

Information is subject to change. To learn more, visit Laerdal.com