



Laerdal[®]
helping save lives

Realistic and Affordable Infant CPR Training



Baby Anne[®]

The Baby Anne manikin has been developed to provide effective infant CPR training without compromising realism or quality.

Offering all the essential features necessary for learning infant CPR, this low-cost, lifelike infant makes hands-on practice possible for every student.

The latest enhancements to Baby Anne include modernized clothing and improved airway management.



Baby Anne®

Benefits

- Affordability allows increased hands-on practice
- Lightweight and portable
- Durable construction allows unequalled long-term use

Features

- Oral and nasal passages allow realistic nose pinch required for mouth-to-nose ventilation
- Head tilt/chin lift and jaw thrust allow students to practice correctly all maneuvers necessary when resuscitating a real patient
- Realistic airway management, with airway closure during flexion and extension of the neck, and opening only in neutral “sniffing” position
- Realistic chest compliance means students can experience the proper technique required for chest compressions on infants
- Foreign-body airway obstruction feature allows the release of a foreign-body obstruction to be practiced through back blows and chest-thrust techniques
- Economical disposable airways for quick and easy clean up
- Removable and reusable faces for convenient and affordable maintenance



Ordering Information

Baby Anne

130-01050	Baby Anne
130-03050	Baby Anne Dark Skin
131-01050	Baby Anne 4-Pack
131-03050	Baby Anne 4-Pack Dark Skin

Spare Parts and Consumables

130-10050	Baby Anne Face Connector 5-pack
130-10150	Baby Anne Airways 24-pack
130-10250	Baby Anne Jaw
130-10350	Baby Anne Dark Skin Jaw
130-10450	Baby Anne Faces with Face Connectors 6-pack
130-10550	Baby Anne Dark Faces with Face Connectors 6-pack

130-10650	Baby Anne Rib Plate 10-pack
130-10750	Baby Anne Pants with Rib Plate
130-10850	Baby Anne Foreign Objects 10-pack
130-50450	Baby Anne Softpack
131-50450	Baby Anne 4-pack Carry Bag
15120103	Manikin Face Shields (6 rolls of 36)
152400	Manikin Wipes (pkg. 50)
152401	Manikin Wipes (pkg. 1200)

Need spare parts and consumables for the previous version of Baby Anne?

All spare parts and consumables listed above can be used on the previous version of Baby Anne, but some of the parts need to be used with others to work on the manikin.



<i>If you are ordering...</i>	<i>Then you will also need...</i>
Airways 24-pack (130-10150)	Face Connector 5-Pack (130-10050) or Faces 6-Pack (130-10450 light/130-10550 dark) (Faces 6-Pack comes with 6 Face Connectors included)
Face Connector 5-Pack (130-10050)	Airways 24-pack (130-10150)
Faces 6-Pack (130-10450 light/130-10550 dark)	Airways 24-pack (130-10150)