



# TRAININGSMODELL FÜR ABDOMINALE UNTERSUCHUNGEN

60000

Zum Üben von Palpation, Perkussion und Auskultation während einer abdominalen Untersuchung.

Ideal für die Vorbereitung und Beurteilung der OSCE-Prüfung.



*Nie war Lernen realistischer*

[limbsandthings.com](http://limbsandthings.com)

## FÄHIGKEITEN:

- Kenntnisse der abdominalen Regionen und der Anatomie
- Unterscheidung von normalen und anormalen Pathologien
- Identifizierung von Aszites, Überblähungen des Darmes und Darmverschluss
- Abtasten der Nieren

## PRODUKTEIGENSCHAFTEN:

- Rumpf mit Knochenstrukturen, Abdomen, Becken und unterem Thorax
- Austauschbare Organe und Krankheitsbilder in verschiedenen Formen und Größen
- Atembewegungen von Leber und Milz
- Realistische Geräusche durch integrierten MP3-Player
- „Shifting Dullness“ und Fluktuationswellen-Techniken können geübt werden
- Simulation normaler und aneurysmatischer Aorta-Pulswellen
- Latexfrei



- **60005** Abdominal Examination Base Unit
- **60006** Abdominal Examination Skin
- **60007** Abdominal Fluid Bag
- **60008** Abdominal Examination Organ Locator
- **60009** Abdominal Wall
- **60010** Abdominal Examination Livers (Set of 3)
- **60011** Abdominal Examination Spleens (Set of 2)

- **60012** Abdominal Examination Kidneys (x2)
- **60013** Abdominal Examination Bladder
- **60014** Abdominal Examination Pathologies (Set of 6)
- **60015** Abdominal Examination Aortas (Set of 2)
- **60016** Abdominal Examination - Distension Set
- **60053** International 5V DC Power Supply
- AA Batteries (x4)
- Carry case