



**Laerdal**<sup>®</sup>  
helping save lives

# Bystander CPR Saves Lives



## Little Junior QCPR

Providing feedback based on just watching students perform can be challenging for instructors. The next generation of the Little Junior manikin includes QCPR measurement and feedback technology that takes the guesswork out of instruction—all from the palm of your hand. Little Junior QCPR works with an Instructor app that allows you to monitor and coach up to six students at once, and a Learner app that lets students monitor their own performance.

Little Junior QCPR helps instructors improve CPR training quality and efficiency in their courses, while adding a fun and engaging element for students that promotes learning and motivation.<sup>(1)</sup>

**All existing Little Junior manikins can be quickly and easily upgraded to QCPR with the Little Junior QCPR Upgrade Kit.**



# Little Junior QCPR

## Feedback Options

Little Junior QCPR's feedback technology provides:

- Real-time feedback on compressions and ventilations
- A performance summary with an overall score and improvement suggestion
- Details on compression release, depth and rate, ventilation volume, and number of compressions/ventilations/cycles


Training based on the use of QCPR technology guided by an instructor has been shown to be superior to instructor-based feedback training for chest compression technical skill acquisition.<sup>(2)</sup>



### QCPR Instructor App

*For instructors to monitor multiple students*

- Free download on App Store
- Connects to 1-6 manikin(s)
- Shows live and summative feedback of CPR performance
- Provides suggestions for improvement, post-session
- Results are saved in the app
- Race mode



### QCPR Learner App

*For students to monitor their own performance*

- Connects to 1 manikin
- Shows live and summative feedback of CPR performance
- Provides suggestions for improvement, post-session
- Results are saved in the app
- Free download on App store and Google Play



### SkillGuide


- Provides a tethered, "plug and play" option
- Can be used standalone or as a supplement to either of the apps

#### ORDERING INFORMATION

128-01050 Little Junior QCPR  
 128-03050 Little Junior QCPR Dark  
 129-01050 Little Junior QCPR 4-Pack  
 129-03050 Little Junior QCPR 4-Pack Dark  
 128-60750 Little Junior QCPR Upgrade Kit

#### Accessories and Consumables

123-30050 SkillGuide with extension cable  
 123-50550 Extension Cable for SkillGuide  
 183210 Junior Airways (pkg. 25)  
 123211 Junior Airways (pkg. 100)  
 15120103 Manikin Face Shields (6 rolls of 36 ea)  
 152400 Manikin Wipes (pkg. 50)  
 152401 Manikin Wipes (pkg. 1200)  
 183010 Junior Faces (pkg. 6)  
 183015 Junior Faces Dark (pkg. 6)



*All existing Little Junior manikins can be quickly and easily upgraded to QCPR with the Little Junior QCPR upgrade kit.*

Visit [Laerdal.com/LittleJuniorQCPR](http://Laerdal.com/LittleJuniorQCPR) to learn more about Little Junior QCPR and how you can upgrade your existing Little Junior manikins.

