



Practice Monitoring Skills to Improve Patient Outcomes

VitalsBridge[™]

The VitalsBridge monitor interface makes it possible to use your own clinical patient monitor with any Laerdal PC-operated simulator or standardized patient. This highly realistic training experience allows healthcare professionals to more accurately evaluate a patient's vital signs and prepares them to provide the best care with confidence.

VitalsBridge™

Specifications:

VitalsBridge provides the opportunity to train on any scenario using a Laerdal LLEAP-based simulator or standardized patient and is compatible with most major patient monitors.



Features	Model 100	Model 200	Model 300
Temperature	✓ 1 Measurement	✓ 1 Measurement	✓ 2 Measurements
SpO2 Pulse Oximetry	✓	✓	✓
Non Invasive Blood Pressure	✓	✓	✓
ECG (5 lead)	☒	✓	✓
Invasive Blood Pressure	☒	☒	✓ 4 Measurements
Capnography (End Tidal CO2)	☒	☒	✓
Communications : Serial/Bluetooth Capable	✓	✓	✓
WiFi/Ethernet Capable	☒	✓	✓
LLEAP Software Interface Compatible	✓	✓	✓
Battery, Portability & Backup	✓	✓	✓

Ordering Information

- 405-00150 VitalsBridge Model 100 Accessories and Cable Kit
Includes: VitalsBridge Model 100 Unit, Software Plug-in, Cable Kit
- 405-00250 VitalsBridge Model 200 Accessories and Cable Kit
Includes: VitalsBridge Model 200 Unit, Software Plug-in, Cable Kit
- 405-00350 VitalsBridge Model 300 Accessories and Cable Kit
Includes: VitalsBridge Model 300 Unit, Software Plug-in, Cable Kit
- 405-86050 VitalsBridge Implementation Service- On-site Installation and Introduction
- 405-00050EXWI Extended Warranty



View, interact with, and control a vital signs monitor used in actual clinical practice with your new or existing Laerdal simulator operated with LLEAP. VitalsBridge uses the simulated digital vital sign waveforms and values from LLEAP and transforms them to signals that are acceptable for a clinical vital signs monitor.

Visit Laerdal.com/VitalsBridge for a complete list of specifications, parts, and Services.

