



Laerdal[®]
helping save lives




Efficient Teaching for Effective Learning

Enhance Your Little Anne[®] Manikins with Objective Feedback

Your existing Little Anne[®] manikins can now be upgraded to include objective QCPR measurement and feedback technology that takes the guesswork out of instruction—all from the palm of your hand.

Little Anne QCPR helps instructors improve CPR training quality and efficiency in their courses, while adding a fun and engaging element for students that promotes learning and motivation⁽¹⁾.

Upgrading your manikins to QCPR is fast and easy with the Little Anne QCPR Upgrade Kit.



Little Anne® QCPR Upgrade Kit

Enhance your Little Anne manikin with QCPR technology that provides:

- Real-time feedback on compressions and ventilations
- A summary with an overall score and improvement suggestion
- Details on compression release, depth and rate, ventilation volume, number of compressions/ventilations/cycles

When used with your choice of feedback option(s), Little Anne QCPR meets the AHA's requirement for feedback devices

Feedback Options:



QCPR Instructor App

For instructors to monitor multiple students

- Free download on App Store
- Connects to 1-6 manikin(s)
- Shows live and summative feedback of CPR performance
- Provides suggestion for improvement, post-session
- Results are saved in the app
- CPR Race feature reinforces quality CPR skills in a fun and engaging way



QCPR Learner App

For students to monitor their own performance

- Free download on App Store and Google Play
- Connects to 1 manikin
- Shows live and summative feedback of CPR performance
- Provides suggestion for improvement, post-session
- Results are saved in the app



SkillGuide™

- Provides a tethered, "plug and play" option
- Can be used standalone or as a supplement to either of the apps

Little Annes manufactured after February 1999 are compatible with the Upgrade Kit. Visit Laerdal.com for more information about how to identify whether your manikins are compatible.

ORDERING INFORMATION

123-60750 Little Anne QCPR Upgrade Kit

ACCESSORIES AND CONSUMABLES

- 123-30050 SkillGuide with extension cable
- 020300 Little Anne Airways (pkg. 24)
- 020301 Little Anne Airways (pkg. 96)
- 15120103 Manikin Face Shields (6 rolls of 36 ea)
- 152400 Manikin Wipes (pkg. 50)
- 152401 Manikin Wipes (pkg. 1200)
- 310210 Manikin Faces (pkg. 6)
- 310215 Manikin Faces, Dark Skin (pkg. 6)



Visit Laerdal.com/LittleAnneQCPR to learn more about Little Anne QCPR and how to upgrade your Little Anne manikins.