



**Laerdal**<sup>®</sup>  
helping save lives



## Bystander CPR Saves Lives


### Little Junior QCPR

#### **Enhance Your Little Junior Manikins with QCPR Feedback**

Your existing Little Junior manikins can now be upgraded to include objective QCPR measurement and feedback technology that takes the guesswork out of instruction – all from the palm of your hand.

Little Junior QCPR helps instructors improve CPR training quality and efficiency in the courses, while adding a fun and engagement element for students.

Upgrading your manikins to QCPR is quick and easy with the Little Junior QCPR Upgrade Kit.



# Little Junior QCPR

## Feedback Options

Little Junior QCPR's feedback technology provides:

- Real-time feedback on compressions and ventilations
- A performance summary with an overall score and improvement suggestion
- Details on compression release, depth and rate, ventilation volume, and number of compressions/ventilations/cycles


Training based on the use of QCPR technology guided by an instructor has been shown to be superior to instructor-based feedback training for chest compression technical skill acquisition.\*



### QCPR Instructor App

*For instructors to monitor multiple students*

- Free download on App Store
- Connects to 1-6 manikin(s)
- Shows live and summative feedback of CPR performance
- Provides suggestions for improvement, post-session
- Results are saved in the app
- Race mode



### QCPR Learner App

*For students to monitor their own performance*

- Connects to 1 manikin
- Shows live and summative feedback of CPR performance
- Provides suggestions for improvement, post-session
- Results are saved in the app
- Free download on App store and Google Play



### SkillGuide

- Provides a tethered, "plug and play" option
- Can be used standalone or as a supplement to either of the apps

#### ORDERING INFORMATION

128-60750 Little Junior QCPR Upgrade Kit

#### Accessories and Consumables

123-30050 SkillGuide with Extension Cable  
183210 Junior Airways (pkg. 25)  
183211 Junior Airways (pkg. 100)  
15120103 Manikin face shields (6 rolls of 36 ea)  
152400 Manikin Wipes (pkg. 50)  
152401 Manikin Wipes (pkg. 1200)  
183010 Junior Faces (pkg. 6)  
183015 Junior Faces Dark Skin (pkg. 6)



Visit [Laerdal.com/LittleJuniorQCPR](http://Laerdal.com/LittleJuniorQCPR) to learn more about Little Junior QCPR and how to upgrade your Little Junior manikins.