




**Laerdal**<sup>®</sup>  
helping save lives



---

Helping you train for  
the highest level of care

## Nursing Anne Simulator

From basic assessment to advanced practice, Nursing Anne Simulator allows students and nursing professionals the ability to practice core skills, enhance clinical knowledge, and experience highly realistic patient encounters that will prepare them to provide the highest level of care. This engaging and robust platform enables nurses to practice a wide range of competencies both as a skills trainer and as a highly realistic female simulator.

# Nursing Anne Simulator

## AIRWAY AND BREATHING

- Realistic airway with uvula, epiglottis, vocal cords and esophagus
- Various oxygen delivery methods
- Bag-mask ventilation
- Oral/nasal intubation
- Right mainstem intubation
- Trach care and tracheal suctioning (simulated)
- Spontaneous breathing synchronized with selected breath rate (0-60 bpm)
- Anterior and posterior lung sound auscultation sites

## PATIENT HANDLING

- Full articulation of arms and legs for realistic patient handling procedures
- Head can be flexed into chin to chest position and remain flexed until repositioned
- Sits unassisted and bends at waist to tripod position

## PATIENT ASSESSMENT AND CARE

- Palpable anatomy for assessment and site location including clavicle, sternum, spine, anterior superior iliac spine, pubic symphysis, greater trochanter, and scapula
- Wig for hair care procedures
- Ear canal for practice of irrigation and cleaning
- Oral care and hygiene

- Accurate anatomical landmarks for insertion of NG Tube to correct measurement
- Lavage and gavage

## MEDICATION ADMINISTRATION

- Bilateral pre-ported IV arms with capability for intravenous bolus or push through intravenous infusion
- Bilateral deltoid, ventral gluteal, dorsal gluteal and thigh injection sites
- G tube for feeding or medication administration

## URINARY AND BOWEL SKILLS/FEATURES

- Female genitalia with realistic anatomy
- Urinary catheterization with ability for proper fluid return
- Enema administration and simulated rectal suppositories

## CIRCULATORY SKILLS/FEATURES

- Bilateral carotid, brachial, radial, femoral, and pedal pulses
- Pulse palpation is detected and logged
- Bilateral measurement of non-invasive blood pressure (auscultated or palpated)
- Heart sounds synchronized with ECG

## OTHER FEATURES

- Midclavicular central line catheter for site care and infusion procedure
- Pre-recorded vocal sounds can be played at a controllable volume
- Blinking eyes with adjustable blink rate
- Ability to open, close or partially close eyes for consciousness cue
- Interchangeable pupils (normal, dilated, constricted)
- CPR capable
- ECG monitoring capabilities when used with a simulated patient monitor
- Defibrillation training with BLS ECG rhythms when used with ShockLink™
- Four independently controlled auscultation areas for bowel sounds
- Fetal heart tone sounds
- 4 complimentary patient cases available on Laerdal.com

## PATIENT MONITOR

- Touchscreen simulated patient monitor provides concise clinical feedback for physiological parameters
- Simulated parameters for SimPad and SimPad PLUS operation: ECG, etCO<sub>2</sub>, HR, SpO<sub>2</sub>, BP, RR, and Temperature
- LLEAP operation provides additional simulated parameters including CO<sub>2</sub>, CVP, ICP, anesthetic agent, TOF, cardiac output, and more. In addition, the LLEAP monitoring options can display x-ray, custom images, and custom video

## ORDERING INFORMATION

320-05050 Nursing Anne Simulator, Medium

320-04050 Nursing Anne Simulator, Dark

Full body simulator with articulating arms and legs. Includes brunette wig, brown pupil set, ostomy set, blood pressure cuff, lubricant spray, simulated blood, adult gown and quick setup guide.

### Licenses

204-50150 LLEAP for SimPad PLUS Software License

400-01050 LLEAP PC Software License

### Operation and Patient Monitoring Options

204-30101 SimPad PLUS (handheld remote)

400-10201 Laptop Instructor PC-Patient Monitor

400-09201 Tablet Instructor PC-Patient Monitor

400-09501 Rugged Instructor PC-Patient Monitor

400-29301 All in One Panel PC Instructor PC-Patient Monitor

### Accessories

212-29650 Headset and Mic with USB (LLEAP Operation)

200-30950 SimPad Headset

185-10050 ShockLink System

(adapter cables may be required)

375-51001-M Multi-Venous Arm, Adult Female, Left

320-21150 Wig; Blonde

213-20050 Blue Pupil Inserts

320-27050-M Mastectomy Module Medium

320-27050-D Mastectomy Module Dark

320-22550 Subcutaneous Injection Pads, 50

320-23150 Silicone Dressing, IV - 50 Pack

320-23250 Silicone Dressing, 4x4 - 50 Pack

320-23350 Silicone Dressing, 8x10 - 10 Pack

320-24050-M Wound Assessment and Care Kit Medium

320-24050-D Wound Assessment and Care Kit Dark

## Operation

Nursing Anne Simulator can be operated by a SimPad/SimPad® PLUS remote or a LLEAP instructor computer and can utilize a simulated patient monitor.



## SimStore™ Scenarios

Easily integrate Nursing Anne Simulator into your training program by utilizing any of the 80+ validated, standardized scenarios from industry experts such as the National League for Nursing specifically designed to address core nursing competencies including:

- Fundamentals of Nursing
- Pharmacology
- Health Assessment
- Medical-Surgical
- Mental Health

*SimStore scenarios include support materials that contain patient charts and documentation as well as information about preparation for the scenario, learning objectives, student information, debriefing points, and more*

## Services

Laerdal offers comprehensive Educational Services to support your curriculum, and Technical Services to meet your product servicing needs. Whether you are setting up a one-room simulation site or a multi-bed simulation center, Laerdal can help.

Visit [Laerdal.com/NursingAnneSimulator](http://Laerdal.com/NursingAnneSimulator) for a complete list of specifications, parts, accessories and consumables

