



Laerdal[®]

helping save lives



Nursing Anne Simulator

Helping you train for the highest level of care

From basic assessment to advanced practice, the Nursing Anne Simulator platform allows learners and nursing professionals the ability to practice core skills, enhance clinical knowledge, and experience highly realistic patient encounters that will prepare them to provide the highest level of care. This engaging and robust platform enables nurses to practice a wide range of competencies both as a skills trainer and as a highly realistic patient simulator.

Nursing Anne Simulator

Airway and Breathing

- Realistic airway with uvula, epiglottis, vocal cords and esophagus
- Various oxygen delivery methods
- Bag-mask ventilation
- Oral/nasal intubation
- Right mainstem intubation
- Trach care and tracheal suctioning (simulated)
- Spontaneous breathing synchronized with selected breath rate (0-60 bpm)
- Anterior and posterior lung sound auscultation sites

Patient Handling

- Full articulation of arms and legs for realistic patient handling procedures
- Head can be flexed into chin to chest position and remain flexed until repositioned
- Sits unassisted and bends at waist to tripod position

Patient Assessment and Care

- Palpable anatomy for assessment and site location including clavicle, sternum, spine, anterior superior iliac spine, pubic symphysis, greater trochanter, and scapula
- Wig for hair care procedures
- Ear canal for practice of irrigation and cleaning
- Oral care and hygiene
- African Descent and Geriatric modules available for added realism and diversity in patient population.

- Accurate anatomical landmarks for insertion of NG Tube to correct measurement
- Lavage and gavage

Medication Administration

- Bilateral pre-ported IV arms with capability for intravenous bolus or push through intravenous infusion
- Bilateral deltoid, ventral gluteal, dorsal gluteal and thigh injection sites
- NG tube for feeding or medication administration

Urinary and Bowel Skills/Features

- Genitalia with realistic anatomy
- Urinary catheterization with ability for proper fluid return
- Enema administration and simulated rectal suppositories

Circulatory Skills/Features

- Bilateral carotid, brachial, radial, femoral, and pedal pulses
- Pulse palpation is detected and logged
- Bilateral measurement of non-invasive blood pressure (auscultated or palpated)
- Heart sounds synchronized with ECG

Other Features

- Midclavicular central line catheter for site care and infusion procedure
- Pre-recorded vocal sounds can be played at a controllable volume
- Blinking eyes with adjustable blink rate
- Ability to open, close or partially close eyes for consciousness cue
- Interchangeable pupils (normal, dilated, constricted)
- CPR capable
- ECG monitoring capabilities when used with a simulated patient monitor
- Defibrillation training with BLS ECG rhythms when used with ShockLink™
- Four independently controlled auscultation areas for bowel sounds
- Fetal heart tone sounds
- 4 complimentary patient cases available on Laerdal.com

Patient Monitor

- Touchscreen simulated patient monitor provides concise clinical feedback for physiological parameters
- Simulated parameters for SimPad and SimPad PLUS operation: ECG, etCO₂, HR, SpO₂, BP, RR, and Temperature
- LLEAP operation provides additional simulated parameters including CO₂, CVP, ICP, anesthetic agent, TOF, cardiac output, and more. In addition, the LLEAP monitoring options can display x-ray, custom images, and custom video

ORDERING INFORMATION

320-05050 Nursing Anne Simulator, Medium

320-04050 Nursing Anne Simulator, Dark

Full body simulator with articulating arms and legs. Includes brunette wig, brown pupil set, ostomy set, blood pressure cuff, lubricant spray, simulated blood, adult gown and quick setup guide.

Operation and Patient Monitoring Options

204-30101 SimPad PLUS (handheld remote)

400-10201 Laptop Instructor PC-Patient Monitor

400-09201 Tablet Instructor PC-Patient Monitor

400-09501 Rugged Instructor PC-Patient Monitor

400-29301 All in One Panel PC Instructor PC-Patient Monitor

Licenses

204-50150 LLEAP for SimPad PLUS Software License

400-01050 LLEAP PC Software License

Accessories

212-29650 Headset and Mic with USB (LLEAP Operation)

200-30950 SimPad Headset

185-10050 ShockLink System (adapter cables may be required)

375-51001-M Multi-Venous Arm, Adult Female, Left

320-21150 Wig; Blonde

320-21050 Wig; Long Dark

213-20050 Pupil Inserts; Blue

320-27050-M Mastectomy Module; Medium

320-27050-D Mastectomy Module; Dark

320-22550 Subcutaneous Injection Pads - 50 pack

320-28050-M Male Genitalia Kit; Medium

320-24050-M Wound Assessment and Care Kit; Medium

320-24050-D Wound Assessment and Care Kit; Dark

320-25550-M Staged Wound Feet Set; Medium

320-29050-D African Descent Front Skin Dark

320-21250 African Descent Module Wig; Dark Curly

Geriatric Essentials Kit available for pre-order. Contact us for more details.

Operation

Nursing Anne Simulator can be operated by a SimPad/SimPad® PLUS remote or a LLEAP instructor computer and can utilize a simulated patient monitor.



Optimize Your Investment

Our team of educators and service technicians have years of industry related experience that can help you successfully implement and sustain your simulation program. We're here to help you get the most out of your simulation training with support when and where you need it. From comprehensive long-term packages to a full range of individual solutions you may choose from:

- ValuePlus Technical Services Packages
- Educational Pathway Education Services Packages
- On-site technical services
- eLearning courses
- On-site educational courses
- Blended learning workshops
- Help Desk phone and email support
- Self-guided how-to video library
- Virtual courses

Visit Laerdal.com/NursingAnneSimulator for a complete list of specifications, parts, accessories and consumables

