

Airway Management



Laerdal Airway Management Solutions

www.laerdal.com



Laerdal[®]
helping save lives

The BAG™ II Disposable Resuscitator



Performance plus value is in The BAG II ... the bright yellow disposable resuscitator by Laerdal. Temperature stability, versatility, space-saving configuration, and ergonomic body make The BAG II the right choice.

The BAG II Resuscitator's excellent performance and ergonomic design offer all the benefits of Laerdal's reusable resuscitator with the convenience of a disposable, making it an important supplement to our complete airway product offering for EMS, Hospital and the Military.

- Available in Adult, Child, and Infant sizes – allowing clinicians a closer match to the volume needs of the patient without hyperventilation
- Adjustable hook and loop handle design along with dual 360° swivels helps reduce hand fatigue
- Latex-free for latex-sensitive patients and staff



HEALTHCARE




EMS




MILITARY

sable Resuscitator



The BAG II Disposable Resuscitator Features



Complies with ISO 10651-4 2002 and ASTM: F920-93 product standards

Laerdal® Silicone Resuscitator

Laerdal Sili



The latest from the leader in reusable resuscitators for over 40 years, the Laerdal Silicone Resuscitator (LSR) has exceptional bag re-expansion and superior feel for lung compliance. The improved design enhances its convenience and maintains the outstanding quality and economy that healthcare professionals rely on. All the excellence and integrity of the original!

Preferred by healthcare providers for excellent performance, economy, and reliability, Laerdal Silicone Resuscitators have proven to be the perennial gold standard by which other resuscitators are measured.

- Collapsible for compact storage to save space
- Material allows for long product life and has been tested to withstand high-level disinfection and steam sterilization
- Patient Valve allows for spontaneous breathing patients to inhale oxygen with minimal resistance
- Latex-free for latex-sensitive patients and staff



HEALTHCARE



EMS



MILITARY

cone Resuscitator



Laerdal Silicone Resuscitator Features


Private Labeling Available!
Call for Details

Silicone Bag Design Features:

- Quick re-expansion even in extreme temperatures
- Transparent silicone and polysulfone parts helps to identify visual contaminants

Private Label Features:

- Personalized labeling on neck coupling helps prevent loss or theft



Patient Valve Features:

- 35 cm H₂O Pressure Patient Valve is standard on all Paediatric and Preterm models
- Integrated Intake/Reservoir Valve allows for use of ambient air in case of oxygen depletion

Mask Features:

- Transparent to allow for visualization of patient's mouth and nose
- Unique multi-function Mask Cover provides additional support or to secure the Laerdal Silicone Mask to patient's face for two-handed BVM operation (optional headstrap available)

Laerdal® Pocket Mask®

La



The 3M Filtrete hydrophobic filter and low-resistance one-way valve combination acts to protect the rescuer from patient secretions without compromising lifesaving ventilations to the patient. Others may look like the Laerdal Pocket Mask, but there's only one original!

The Laerdal Pocket Mask offers superior protection and durability and is specifically designed to be used by the professional rescuer and first responder.

- Latex-free and powder-free versions of nitrile gloves and antimicrobial wipes* help maximize user protection
- Utilize in CPR instruction courses to teach mouth-to-barrier rescue breathing
- Distance from patient allows rescuer to see chest rise and fall



HEALTHCARE



EMS



MILITARY

Laerdal Pocket Mask




Laerdal Pocket Mask Features

Case Features:

- Hard or Soft case packaging options
- Private labeling for your organization available
- Military green and camouflage

Private Labeling Available!
Call for Details



Both the One-Way Valve and the Disposable Filter can be used with all previous versions of Laerdal Pocket Mask devices.

One-Way Valve Features:

- Oxygen inlet version facilitates delivery of supplemental oxygen to breathing and non-breathing patients
- Intelligently engineered one-way valve and disposable 3M Filtrete hydrophobic filter protects patient and rescuer
- Pre-inflated cuff allows for ease of application and makes effective seal

* Wipes and gloves are expiration dated for user protection.

Laerdal® Paediatric Pocket Mask®



Paediatric Pocket Mask Features

Designed specifically for infants and children, the Laerdal Paediatric Pocket Mask has the same one-way valve as our adult Pocket Mask!

Benefit Features:

- Utilize in CPR instruction courses to teach mouth-to-barrier rescue breathing
- Distance from patient allows rescuer to see chest rise and fall

One-Way Valve Features:

- Oxygen inlet version facilitates delivery of supplemental oxygen to breathing and non-breathing patients
- Intelligently engineered one-way valve helps protect patient and rescuer
- Flanged cuff allows for ease of application and makes effective seal



Glove and Wipe Features:

- Latex-free and powder-free versions of nitrile gloves and antimicrobial wipes help maximize user protection*

**Wipes and gloves are expiration dated for user protection.*

Laerdal® Face Shield



Laerdal Face Shield Features

This single-use CPR barrier is ideal for carrying every day, for inclusion with CPR training courses, and in first aid kits.


Private Labeling Available!
Call for Details

Filter Features:

- 3M Filtrete hydrophobic filter helps reduce risk of blood or saliva contact between rescuer and patient
- Bi-directional filter

Design Features:

- Oversized shield design allows for an unrestricted patient exhalation route without compromising seal quality or effective ventilation achievement



Laerdal® Face Shield Key Ring Features



Design Features:

- Available in multiple colors
- Key Ring clips on belt loop or handbag, and fits in pocket

Imprinted graphic Directions for Use

Thomas™ Tube Holder

Tho



All Laerdal Thomas Tube Holders share features that set them apart from other airway tube holders. For instance, the Laerdal Thomas Tube Holder is the only tube holder on the market that has a quick-set screw for a secure hold after tube placement. Tape is not necessary.

Laerdal Thomas Tube Holder is a unique (and better) tube holder over others on the market, no matter what size or color you choose.

- *Eliminates the need for cumbersome taping or other fixatives*
- *Slide stick facilitates fast, efficient application of the head strap around the patient's neck without movement of the head*
- *Large opening allows for application of a wide range of single and double-lumen airway tubes*
- *Ventilations may continue during application, adjustment and removal of the Thomas Tube Holder*



HEALTHCARE



EMS



MILITARY

Thomas Tube Holder



Thomas Tube Holder Features

Wide Padded Velcro® Head Strap Features:

- Provides secure and comfortable fit, even when wet
- Color-coded for easy identification between adult and pediatric
- Will not stretch

Suctioning Features:

- Larger access opening to allow for ease of suctioning and visual inspection of mouth




Military

Adult

Pediatric

Slide Stick Features:

- Facilitates fast, efficient application of the head strap around the patient's neck without movement of the head



Security Features:

- Set Screw mechanism quickly secures ET tubes and allows for fast readjustment

Enlarged Access Port Features:

- Accommodates multiple size ET tubes, as well as Combitubes®, King Tubes, LMAs, and LMA Supreme™
- BVM does not need to be removed to attach the Thomas Tube Holder
- Allows for in-use suctioning of the oropharynx

Airway Solutions

Laerdal Pocket Mask

- 82001133** Laerdal Pocket Mask
w/Gloves & Wipe in Yellow Hard Case
- 82004033** Laerdal Pocket Mask
w/Gloves & Wipe in Blue Soft Pack
- 82004133** Laerdal Pocket Mask
w/Gloves & Wipe in Black Soft Pack
- 829951** Laerdal Pocket Mask w/Gloves & Wipe
in Olive Green Hard Case
- 82004233** Laerdal Pocket Mask w/Gloves & Wipe
in Camo Soft Pack

Accessories

- 820410** One-Way Valve (each)
- 820610** Filter (each)
- 820700** Update Pack (1 Filter, 1 Valve)
- 830110** Head Strap, Adjustable (pkg. 10)

Laerdal Pocket Mask Oxygen Version

- 83001133** Pocket Mask w/O₂ Inlet and Head Strap
w/Gloves & Wipe in Yellow Hard Case
- 83004033** Pocket Mask w/O₂ Inlet and Head Strap
w/Gloves & Wipe in Blue Soft Pack
- 83004133** Pocket Mask w/O₂ Inlet and Head Strap
w/Gloves & Wipe in Black Soft Pack
- 83001633** Pocket Mask w/O₂ Inlet in Poly Bag

Laerdal Paediatric Pocket Mask

- 820054** Laerdal Paediatric Pocket Mask
w/Gloves & Wipes in Re-Sealable
Poly Bag

Laerdal Face Shield

- 46000001** Laerdal Face Shield CPR Barrier (pkg. 50)

Laerdal Face Shield Key Ring

- 460009** Laerdal Face Shield CPR Barrier Key Ring
in Multiple Colors – 8 Blue, 8 Black,
9 Red (pkg. 25)
- 460017** Laerdal Face Shield CPR Barrier Key Ring,
Camouflage (pkg. 25)

For a complete listing of airway accessories, visit www.laerdal.com
or contact your Authorized Laerdal Distributor

The BAG II Disposable Resuscitator

- 845011** The BAG II Disposable Resuscitator;
w/Adult #5 Mask
- 845021** The BAG II Disposable Resuscitator;
w/Child Mask
- 845031** The BAG II Disposable Resuscitator;
w/Infant Mask
- 845041** The BAG II Disposable Resuscitator;
w/Adult #4 Mask

Accessories

- 845050** Disposable Mask #3 w/Inflation Port (pkg 20)
- 845060** Disposable Mask #4 w/Inflation Port (pkg 20)
- 845070** Disposable Mask #5 w/Inflation Port (pkg 20)
- 845080** Disposable Expiratory Diverter (pkg 10)
- 845040** Disposable PEEP Valve, 5-20cm H₂O (pkg 10)*

**also compatible with the LSR*

Laerdal Silicone Resuscitator

- 87005133** Silicone Resuscitator; Adult, complete in carton
- 86005133** Silicone Resuscitator; Child, complete in carton
- 85005133** Silicone Resuscitator; Preterm, complete in carton


Accessories

- 851500** Infant Silicone Mask No. 0/0
- 851600** Infant Silicone Mask No. 0/1
- 851700** Infant Silicone Mask No.2
- 870220** Adult Silicone Mask No.4/5+ w/Mask Cover
- 860220** Child Silicone Mask No.3/4 w/Mask Cover
- 875200** Multi-Function Mask Cover
for Adult Mask 4-5+
- 865200** Multi-Function Mask Cover
for Child Mask 3-4
- 851350** Patient Valve w/35cm H₂O Pressure Relief
Valve & Lock Clip
- 560200** Patient Valve
- 851250** Patient Valve w/35cm H₂O Pressure Relief Valve
- 850500** Expiration Diverter (OD 30 mm)
- 850900** Manometer Connector

Thomas Tube Holder

- 600-10000** Thomas Tube Holder; Adult
- 600-20000** Thomas Tube Holder; Pediatric
- 600-30000** Thomas Tube Holder; Adult (Green)

For additional information or to order, contact your Authorized Laerdal Distributor!



Laerdal[®]
helping save lives