



Laerdal[®]

helping save lives



How to Use SimPad[®] eLearning Resources

Welcome to SimPad eLearning!

The SimPad eLearning content is designed to help you achieve success quickly.

SimPad® eLearning Resources

Content Organization

The eLearning content may be used in conjunction with Laerdal onsite instruction or as a stand-alone resource. Using a blend of software demonstrations and downloadable job aids, you will learn how to get the SimPad® system started and facilitate simulation sessions using Automatic and Manual modes.

Online topics are organized in two sections.

Getting Started -

This section provides an overview of the SimPad® system starting with components, assembly, and start-up-procedures. The section also reviews topics such as how to register SimPad on SimStore™, synchronize content with SimStore and connect to patient monitor.

Making Things Happen -

This section provides comprehensive tips for using SimPad in both Automatic and Manual modes. This section also contains information about vocal sounds, the patient monitor, and viewing the debriefing log.



How to Use the eLearning Resources



Software Simulations: Narrated tutorials show how to navigate the software and perform tasks. Most tutorials are less than three minutes. View the tutorials with your SimPad® in hand - watch a topic then practice. You may view the tutorials on PC, and MAC operating systems.



Job Aids: Downloadable illustrations provide step-by-step directions for basic and advanced tasks. Bring the job aid directly into the lab when you need just-in-time support.

877-LAERDAL (523-7325)
laerdal.com



Laerdal®
helping save lives