




Changing Your Simulation Perspective

SimView™

Debriefing learners after the simulation is the value-added benefit of any simulation-based learning event. SimView is an integrated debriefing solution that captures and records audio, video, relevant patient monitor data and event log info in one file giving you and your learners every opportunity to evaluate each simulation experience completely.





Capture

Complete integration for one comprehensive solution

- Incorporate synchronized audio, video, patient monitor and event log into one debrief file
- Record training event with simulators, task trainers or standardized patients
- Digital HD (up to 4 streams), analog, and web camera configurations
- Capture comments of other observers during the simulation

Share

Share debriefing data to improve learning outcomes

- Stream video to multiple locations
- Call up debrief files from other computers on the network
- Share files with faculty or students
- Export videos to embed in presentations or store for research

Debrief

Raise the value of debriefing

- Review entire scenario or easily jump to key events within the scenario to reinforce learning objectives
- Annotate the event log for comments, events, or medications, correct/incorrect interventions
- Review debriefing input incorporated from preprogrammed SimStore scenarios
- Upload event log information to student transcripts in SimManager

Implementation, training and support

Whether you are setting up a one-room simulation site or a multi-bed simulation center, the Laerdal Services Team can consult with you to understand your simulation debriefing needs and help configure the optimal SimView solution that works for your current needs and future plans

ORDERING INFORMATION

SimView includes: Desktop PC, LCD Monitor, Keyboard, Mouse, Audio/Video Adapter Cables, 5 Session Viewer Licenses, Directions for Use

- 216-00441 SimView- Digital/Analog
- 216-00442 SimView- Digital/Analog (AVS Upgrade)
- 216-00341 SimView- Digital
- 216-00342 SimView- Digital (AVS Upgrade)

Optional Accessories:

- 216-05033 Digital Fixed Network Camera
- 216-05133 Digital Pan Tilt Zoom (PTZ) Network Camera
- 216-05233 SimView General Area Microphone

Services:

- 216-10050 SimView Elearning
- 216-83050 Basic Setup
- 216-83250 Advanced Setup
- 216-83150 Custom Setup

SimView PC Specifications	
Manufacturer	Hewlett-Packard
Model	HP Z420 WorkStation
Processor/CPU	Intel Xeon E5-1620 3.6 10M 1600-4C
Hard Drive 1	265GB SATA 1st SSD
Hard Drive 2	1 Terabyte 720 RPM SATA 2nd HDD
RAM	8 GB DDR3-1600
Graphic	GeForce GTX 660 3GB PCI Express
Screen	Planar PLL2210W LCD - TET - 22" 1920 x 1080
Mouse	HP USB Optical Scroll
Keyboard	HP USB Standard Keyboard
Operating Software	Windows Server 2012 Standard - License and Media - 2 processors - OEM - DVD - 64bit
DVD	16x SuperMulti DVDRW SATA 1st ODD
Video Capture Driver	Hikvision DS40xx - 64bit/H.264 Compression
Warranty	HP 3/3/3 Warranty Global + 4h 9-5 onsite/DMR add-on