



Laerdal[®]
helping save lives



Bystander CPR Saves Lives

Little Family Pack

The Little Family Pack is a convenient package of age-specific CPR training manikins with a wheeled carry case for ease of transport and storage.

Now, Little Anne and Little Junior come with QCPR feedback technology to help instructors improve CPR training quality, efficiency, and learner engagement, while the updated Baby Anne features modernized clothing and improved airway management for enhanced realism.



Little Family Pack

The Little Family Pack of manikins includes:

- Little Anne QCPR and Little Junior QCPR, each equipped with Bluetooth technology to deliver CPR feedback via free apps on smartphones and tablets
- Baby Anne featuring updated clothing and improved airway management, with airway closure when the neck is hyperextended




Feedback Options for Little Anne QCPR and Little Junior QCPR*



QCPR Instructor App

For instructors to monitor multiple students

- Free download on App Store
- Connects to 1-6 manikin(s)
- Shows live and summative feedback of CPR performance
- Provides suggestions for improvement, post-session
- Results are saved in the app
- Race mode



QCPR Learner App

For students to monitor their own performance

- Connects to 1 manikin
- Shows live and summative feedback of CPR performance
- Provides suggestions for improvement, post-session
- Results are saved in the app
- Free download on App store and Google Play



SkillGuide

- Provides a tethered, "plug and play" option
- Can be used standalone or as a supplement to either of the apps

*At this time, Baby Anne is not QCPR enabled.

ORDERING INFORMATION


126-01050 Little Family Pack
126-03050 Little Family Pack Dark Skin

Accessories

123-30050 SkillGuide with Extension Cable
123-60750 Little Anne QCPR Upgrade Kit
128-60750 Little Junior QCPR Upgrade Kit

Already own a Little Anne** or Little Junior?

Upgrade them to QCPR quickly and easily with QCPR Upgrade Kits.



Visit Laerdal.com/us/LittleFamilyPack for more information about the Little Family Pack, accessories, how to upgrade your existing manikins, and more.

