



**Laerdal**<sup>®</sup>

helping save lives



## Resusci Anne<sup>®</sup> Advanced SkillTrainer

---

Mobile ALS skills trainer for advanced healthcare teams

Resusci Anne Advanced SkillTrainer is designed to deliver high-quality team and skills training in a wide range of ALS scenarios. The combination of features and mobility makes Resusci Anne Advanced SkillTrainer an ideal solution for in-situ training.

## Features:

### Airway

- Airway occlusion, monitor chin lift / jaw thrust
- Ventilation with BVM (Bag Valve Mask)
- Supraglottic airway insertion
- Oropharyngeal and nasopharyngeal airways insertion and suctioning
- Selick maneuver
- Realistic chest rise and fall

### ECG, defibrillation and circulation

- Instructor controlled ECG rhythm
- ECG monitoring (3 lead)
- Defibrillation w/paddles compatible
- Correct defibrillation pad placement confirmation
- Auto conversion of ECG w/defibrillation
- External pacing
- Automatic carotid pulse
- Carotid pulse synchronized with ECG
- Carotid pulse check sensors
- Carotid pulse strength dependent on

blood pressure

- Adjustable BP levels (systolic / diastolic)
- Adjustable pulse strength
- Adjustable heart rate

### Debriefing

- Debriefing through recorded events
- Detailed CPR evaluation

### Circulatory Skills and IV / Drug

#### Administration

- IV insertion
- Fluid infusion
- Realistic flashback

#### CPR

- Anatomical landmarks
- Head tilt/chin lift sensor
- Ventilation measurement and feedback
- Chest compression (Guidelines 2015 depth)

- Compression measurement and feedback
- Abdominal thrust
- Compression counter
- Complete release feedback

### Simulation Administration (SimPad® PLUS)

- Event logging
- Remote control
- Create and edit scenarios
- Pre-programed scenarios
- Simulated patient monitor
- Pre-set patient state level
- Wireless operation (WiFi)

### Other

- Independent stomach
- Interchangeable pupils
- Gender: Adult female full-body manikin
- Optional trauma limbs and module set

## Order Information:

Cat. No.	Product Description
151-22000	RA Adv. SkillTrainer AED Link, IV-arm left, incl. ShockLink®
151-23000	RA Adv. SkillTrainer AED Link, IV-arm right, incl. ShockLink®
151-27000	RA Adv. SkillTrainer Paddle, IV-arm left
151-28000	RA Adv. SkillTrainer Paddle, IV-arm right
204-30101	SimPad® PLUS
204-50150	LLEAP for SimPad® PLUS

\* Cat. no. 151-27000 and 151-28000 designed for defibrillators with paddles.

Note: SimPad® PLUS not included in the manikin configurations - order separately.



**Laerdal®**  
helping save lives

For more information, visit [www.laerdal.com](http://www.laerdal.com)