



Laerdal®
helping save lives



Obstetric Solution **SimMom and MamaBirthie**

Combining high and low fidelity simulation
to help improve patient outcome

SimMom is an advanced full body birthing simulator with accurate anatomy and functionality to facilitate multi-professional obstetric training of birth management, with both manual and automatic delivery modes. MamaBirthie is a simple, easy to use obstetric training device which can be paired with a standardized patient or used as a tabletop skills trainer.

Together, SimMom and MamaBirthie offer a solution for delivering efficient and effective obstetric simulations which meet learning objectives around care, communication, and skills.

SimMom Features

Deliveries and drills

- Normal delivery
- Breech presentation
- Assisted deliveries: forceps and vacuum
- Shoulder dystocia
- Cord prolapse
- Eclampsia and preeclampsia
- Maternal collapse
- Postpartum hemorrhage (PPH)
- Sepsis
- Uterine inversion
- Ruptured uterus

Pelvis anatomy

- Amniotic sac
- Dilating cervix
- Ischial spines
- Anus for post-delivery assessment
- Atonic uterus modules (for PPH, uterine inversion and retained placenta)
- Fluids (e.g., blood, stained amniotic fluid and urine)
- Urine catheterization/instillation
- Placenta: complete and incomplete

Airway

- Difficult airway: laryngospasm, tongue edema and pharyngeal swelling
- Right lung, left lung and both lung blockage
- Head tilt/chin lift, jaw thrust
- Suctioning techniques (oral and nasopharyngeal)
- Bag-Valve-Mask ventilation
- Oropharyngeal and nasopharyngeal intubation
- Combitube, LMA and other airway device placement
- Endotracheal intubation (ET)
- Retrograde intubation
- Nasal and oral fiberoptic intubation
- Trans-tracheal jet ventilation
- Right mainstem intubation
- Surgical and needle cricothyrotomy
- Chest tube insertion
- Vocal chords
- Supports Sellick Maneuver

Breathing

- Simulated spontaneous breathing
- Variable respiratory rates
- Bilateral and unilateral chest rise and fall
- Normal and abnormal breath sounds: 4 anterior auscultation sites, bilateral midaxillary sites

Circulation

- Adjustable BP levels: systolic/diastolic
- BP measured manually by auscultation of Korotkoff sounds

- Bilateral carotid, right side brachial & radial pulses, synchronized with ECG
- Pulse strength variable with BP
- Pulse palpation is detected and logged
- Supports NIBP

Vascular

- Pre-ported IV access in both arms
- Subcutaneous and intramuscular injection sites
- Volume infusion
- Realistic flashback

Cardiac

- CPR: real-time compression and ventilation performance
- CPR compressions generate palpable pulses, blood pressure wave form, and ECG artifacts
- Detection and logging of a series of compressions
- Extensive ECG library
- Heart sounds synchronized with ECG
- ECG rhythm monitoring
- 12 lead ECG display
- Defibrillation, cardioversion and pacing

Movement

- Seizure / convulsions
- Able to position on all fours: realistic rotation of the shoulder and hip joints, legs bend at the knees, arms bend at the elbow
- Other positions: supine, semi-recumbent, left lateral, legs in stirrups, McRoberts position

Other

- Sounds: bowel, heart, fetal heart, lungs
- Interchangeable pupils: normal, blown and constricted
- Patient Voice: pre-recorded and custom sounds, instructor can simulate patient's voice

Birth Baby

- Realistic weight and size of full-term neonate
- Palpable limbs and fontanelles
- Anatomical landmarks to identify presentation of baby
- Natural head flex
- Open mouth

Optional (with optional accessories)

- Automatic delivery
- Non-gravid abdomen skin
- Ultrasound training (LSUS)

MamaBirthing Features

Deliveries and drills

- Normal delivery
- Breech presentation
- Assisted deliveries: forceps and vacuum
- Shoulder dystocia
- Vaginal examination
- Abdominal examination

MamaBirthing Anatomy

- Realistic birth canal and pelvic floor to facilitate natural rotation of head as it descends
- Pelvis with palpable landmarks
- Cervix inserts of 4, 6 and 8cm dilation and effacement

BabyBirthing Anatomy

- Anatomical landmarks to identify presentation of baby
- Realistic neck articulation
- Palpable fontanelles
- Simulation of chignon effect

Other

Table clamp for skills training

Waist strap for wearing on a standardized patient

Optimize Your Investment

We're here to help you get the most out of your simulation training with support when and where you need it, from comprehensive Educational Services to support your curriculum to Technical Services to meet your product servicing needs.

Operation

SimMom can be operated by a LLEAP instructor computer and can utilize a simulated patient monitor. Incorporate SimCapture for a complete debriefing solution.

ORDERING INFORMATION:

SimMom 377-03050

Includes: SimMom, Birthing Baby with placenta, 4 Uteri Modules (Cervix, Amniotic Bag, PPH Module and Inverted Uterus), Consumables and Blood Pressure Cuff. MamaBirthing, BabyBirthing with placenta, Cervix Cups, Accessories and Carry Bag.

Visit Laerdal.com for all accessories and consumables.



Developed in partnership with Limbs & Things

Laerdal has given SimMom an easy-to-use, full-body simulation system, whilst Limbs & Things brings a level of realism, both in terms of anatomy and obstetric experience to the product.

