



Laerdal[®]
helping save lives



VitalsBridge[™]



Practice Monitoring Skills to Improve Patient Outcomes

VitalsBridge[™]

The VitalsBridge monitor interface makes it possible to use your own clinical patient monitor with any Laerdal PC-operated simulator or standardized patient. This highly realistic training experience allows healthcare professionals to more accurately evaluate a patient's vital signs and prepares them to provide the best care with confidence.

Specifications:

VitalsBridge provides the opportunity to train on any scenario using a Laerdal LLEAP-based simulator or standardized patient and is compatible with most major patient monitors.



Features	VitalsBridge ⁴			VitalsBridge FM
	Basic	Standard	Plus	
Temperature	✓ 1 Measurement	✓ 2 Measurements	✓ 2 Measurements	—
Heart Rate	✓	✓	✓	✓
Respiration Rate (ECG impedance, capnography if applicable)	✓	✓	✓	✓
SpO2 Pulse Oximetry	✓	✓	✓	✓
Non-Invasive Blood Pressure	✓	✓	✓	✓
ECG (3, 4, or 5 lead)	✓	✓	✓	✓
Capnography (end tidal CO2) Mainstream / Sidestream	—	✓	✓	—
Arterial Blood Pressure	—	✓	✓	—
Central Venous Pressure	—	✓	✓	—
Pulmonary Artery Pressure/Pulmonary Capillary Wedge Pressure	—	✓	✓	—
Intracranial Pressure	—	—	✓	—
Uterine Activity (Toco)	—	—	—	✓
Fetal Heart Rate (Doppler Ultrasound)	—	—	—	✓
Bluetooth LE/USB Serial Capable	✓	✓	✓	✓
VB4 connector SW, Compatible with iOS, Android, and Windows	✓	✓	✓	✓
Wifi/Ethernet Capable	✓	✓	✓	✓
LLEAP Software Interface Compatible	✓	✓	✓	✓
Battery, Portability & Backup	✓	✓	✓	✓

Ordering Information

405-01400	VitalsBridge ⁴ Basic Accessories and Cable Kit Includes: VitalsBridge Basic Unit, Software Plug-in, Cable Kit
405-01500	VitalsBridge ⁴ Standard Accessories and Cable Kit Includes: VitalsBridge Standard Unit, Software Plug-in, Cable Kit
405-01550	VitalsBridge ⁴ Plus Accessories and Cable Kit Includes: VitalsBridge Plus Unit, Software Plug-in, Cable Kit
405-00500	VitalsBridge FM Accessories and Cable Kit Includes: VitalsBridge Fetal Monitor Unit, Software Plug-in, Cable Kit
405-86050	VitalsBridge Implementation Service - On-site Installation and Introduction
405-00050EXWI	Extended Warranty Purchased in 1-year increments
405-EDVT025	2-hour Virtual Training Delivered by VitalsBridge Support Experts



View, interact with, and control a vital signs monitor used in actual clinical practice with your new or existing Laerdal simulator operated with LLEAP. VitalsBridge uses the simulated digital vital sign waveforms and values from LLEAP and transforms them to signals that are acceptable for a clinical vital signs monitor.

**For more information,
visit Laerdal.com/VitalsBridge**



Laerdal®
helping save lives



VitalsBridge™