



Laerdal
helping save lives



SimCapture Mobile Camera

Run your simulations anywhere, anytime

SimCapture Mobile Camera is a video streaming app which, when configured with a SimCapture account, acts as a live, portable camera for SimCapture Cloud sessions.

Capture your training on the move - transport cases, on site, and in situ simulations - with your mobile device, while benefiting from full access to SimCapture functionalities from anywhere. Enhance your experience with offline recording functionality, allowing true mobility in areas with limited or no internet coverage.

Available on both Android and iOS.

Never miss a training opportunity - capture everything with ease

Extend the reach of your SimCapture platform by turning any mobile phone or tablet into a capture device and record your on-the-move sessions. SimCapture includes additional, complementary recording modalities to meet your program needs - contact us to get more information.



Truly mobile

Capture simulations from anywhere, in situ, transport cases, and more, regardless of network coverage.



Up to six mobile devices

Pair up to six mobile devices and capture every angle in one or multiple locations. Additionally, capture up to two patient monitors.

A maximum of four feeds per session is enabled (offline recording is single feed).



Live-stream to SimCapture cloud

Record to the cloud without storage concerns. Annotate videos in the SimCapture interface during simulations, and access them immediately after the session for efficient debriefing. The offline recording feature lets you record and annotate without internet, then uploads to the cloud once the connection is restored.



Fully compatible with LLEAP

Include your patient monitor and events in the simulation data for a complete simulation session that supports efficient debriefing and assessment. LLEAP treats the Mobile Camera the same as a stationary SimCapture camera. With the livestream mode, sessions can be started and stopped from within LLEAP.



Synchronized mobile devices

Recording starts from the SimCapture interface, allowing all mobile devices to start and stop simultaneously. This simplifies the beginning of sessions and ensures high video quality. This also offers higher security against confidentiality breaches.




Ultra quick setup

Log into SimCapture, then enter four letters into your mobile device, and your camera is paired. Simple and easy to use!

Leverage SimCapture's full features with mobile freedom

- Secure, mobile-friendly, and cloud-based
 - Use your mobile device to stream sessions
 - Uncapped data storage is included at no additional charge with the subscription purchase
 - Offline recording and annotations features for sessions with limited/no internet coverage
 - Simulator data capture and visualization
 - Dedicated learner interface with assignments, schedule, and portfolio results
 - Single Sign-On module
 - Center sign-in - directing and tracking users
 - Video annotation and session self-reflection
 - Debrief from anywhere using just a browser
 - Learner, faculty, and facility usage reports
 - Customizable scenarios, roles and permissions
 - Resource tracking and reporting
 - Inventory management features
 - Scheduling and reservation request process
 - Checklist builder and custom assessments
 - Courses and curriculum tracking
 - Learner-driven features and checklists for peer-to-peer and self-directed training
 - Robust suite of assessment and course performance reports
 - Integration with Learning Management Systems
 - Best-in-class training and support
- * Features above are dependent on the subscription tier*



For more information, visit Laerdal.com/SimCapture

©2025 Laerdal Medical. All rights reserved. 2503-1610204763



Laerdal
helping save lives