



**Laerdal**<sup>®</sup>  
helping save lives



## Newborn Anne

Clinical skills training for neonatal healthcare professionals

Newborn Anne represents a full term (40 week), 50th percentile newborn female. With anatomical accuracy and product features designed to focus on the critical resuscitation skills required in the first ten minutes of a newborn's life, Newborn Anne meets the key components of the Neonatal Resuscitation Program<sup>®</sup> and most other neonatal clinical training curricula. Lightweight and portable design allows the manikin to be highly mobile for use in multiple clinical settings.



# Features

## Airway

- Anatomically accurate, realistic airway
- Positioning the newborn to simulate opening the airway via head tilt, chin lift or jaw thrust
- Positive pressure ventilation
- ET tube insertion
- LMA insertion
- Orogastric tube insertion
- Right mainstem intubation
- Sellick maneuver
- Stomach distension (when ET is misplaced)
- Suctioning (Fluids should not be introduced into the airway)
- Meconium module for suction removal

## Breathing

- Bilateral and unilateral (with mainstem intubation) chest rise and fall with mechanical ventilation
- Pneumothorax - Needle thoracentesis left mid axillary

## Cardiac & Circulation

- Anatomical landmarks
- Manual chest compression
- Manual umbilical pulse

## Vascular Access

- Umbilical Vein/ Artery access via patent umbilical cord
- Simulated blood flashback upon cannulation
- IO access in left and right lower leg

## Other Features

- Articulation in all four limbs


## ORDERING INFORMATION

### 220-25050 Newborn Anne manikin

Full-body neonatal manikin, meconium module set, umbilical cord, clamp, airway lubricant, simulated blood, liquid soap, baby powder, IV bag connector tube, pulse bulb tube, baby pants, carry case, directions for use and warranty.

## ACCESSORIES (OPTIONAL)

- 200-03050 Eye kit (interchangeable pupils with normal, dilated & constricted pupils)
- 220-00350 Right & left IO Mandrel and Leg skin kit
- 220-00150 Umbilical cord (4)
- 240-00250 Umbi simulated blood
- 220-00250 IV Bag & Tubing set



## Optimize Your Investment

Our team of educators and service technicians have years of industry related experience that can help you successfully implement and sustain your simulation program. We're here to help you get the most out of your simulation training with support when and where you need it. From comprehensive long-term packages to a full range of individual solutions.

## The tools you need

Be more efficient and achieve lasting results with Laerdal's collection of apps and services. We offer everything you need to plan, operate and scale your simulation with solutions that can be custom-made to fit your training needs.

Visit [Laerdal.com](https://www.laerdal.com) for more information and a complete overview of parts, accessories and consumables

